# CREATININE

Cat. No.	Pack Name	Packaging (Content)
BLT00022	CREA 500	R1: 4 × 100 mL, R2: 1 × 100 mL, R3 zero calibrator: 1 × 10 mL, instruction for use

## INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of creatinine in human serum, plasma and urine on various automatic systems. In combination with other parameters it is intended for screening, monitoring of renal function, diagnosis of renal diseases. For professional use in clinical laboratory only.

## CLINICAL SIGNIFICANCE

Creatinine is a waste product excreted by the kidneys mainly by glomerular filtration. The concentration of creatinine in plasma of a healthy individual is fairly constant, independent from water intake, exercise and rate of urine production. Therefore, increased plasma creatinine values always indicate decreased excretion, i.e. impaired kidney function. Creatinine clearance is a good indicator for the glomerular filtration rate (GFR) which allows better detection of kidney diseases and monitoring of renal function. For this purpose, creatinine is measured simultaneously in serum and urine collected over a defined time period. The serum creatinine levels do not start to rise until renal function has decreased by at least 50 %.

## PRINCIPLE

Creatinine reacts with alkaline picrate to produce an orange-red complex (Jaffe reaction). This is a non-specific reaction and is given by many other substances. Specificity of the assay has been improved by the introduction of a kinetic method, however, the cephalosporin antibiotics are still major interferences<sup>1,2,3,4,5,6</sup>.

The rate of complex formation or change of absorbance, measured at 490–510 nm, is proportional to the creatinine concentration in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

<b>R1</b>	
Sodium hydroxide	240 mmol/L
<b>R2</b>	
Picric acid	26 mmol/L
<b>R3</b>	
Zero calibrator	

## COMPOSITION OF REACTION MIXTURE

Sodium hydroxide	183 mmol/L
Picric acid	5 mmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 490–510 nm may be used, general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034  
 XL MULTICAL 10×3, Cat. No. XSYS0122  
 ERBA NORM 4×5, Cat. No. BLT00080  
 ERBA NORM 10×5, Cat. No. XSYS0123  
 ERBA PATH 4×5, Cat. No. BLT00081  
 ERBA PATH 10×5, Cat. No. XSYS0124

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

### Two reagents method

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

### Monoreagent method

Stability of working reagent: 4 weeks at 2–25 °C in dark

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K<sub>2</sub>-EDTA plasma.

Urine: Collect urine without using additives. If urine must be collected with a preservative for other analytes, only hydrochloric acid (14 to 47 mmol/L urine, e.g. 5 mL 10 % HCl or 5 mL 30 % HCl per liter urine) or boric acid (81 mmol/L, e.g. 5 g per liter urine) may be used. Dilute urine samples using distilled water in 1 + 19 ratio and multiply results by 20.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

<b>Stability in serum / plasma<sup>7</sup>:</b>	7 days at	15–25 °C
	7 days at	2–8 °C
	3 months at	-20 °C

<b>Stability in urine<sup>8</sup>:</b>	2 days at	15–25 °C
	6 days at	2–8 °C
	6 months at	-20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL is recommended. 2 point calibration (blank and calibrator); distilled water or reagent R3 (zero calibrator) are recommended as blank.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS.

## ASSAY PROCEDURE

Wavelength: 492 (490–510) nm

Cuvette: 1 cm

### Two reagents method

	Reagent blank	Calibrator	Sample
Reagent 1	0.80 mL	0.80 mL	0.80 mL
Sample	–	–	0.05 mL
Calibrator	–	0.05 mL	–
Distilled water (R3*)	0.05 mL	–	–

Mix and after 1–5 min. incubation (at 37 °C) add:

Reagent 2	0.20 mL	0.20 mL	0.20 mL
-----------	---------	---------	---------

Mix and after exactly 1 minute incubation read the initial absorbance for blank A<sub>bl</sub>, sample A<sub>sam</sub> and calibrator A<sub>cal</sub>. Exactly after 2 minutes read the final absorbance of blank A<sub>bl</sub>, sample A<sub>sam</sub> and calibrator (standard) A<sub>cal</sub>. Calculate resulting absorbance as the difference between the final and initial absorbance (ΔA/min). \* see Note.

## Monoreagent method

	Reagent blank	Calibrator (Standard)	Sample
Working reagent	1.00 mL	1.00 mL	1.00 mL
Sample	–	–	0.05 mL
Calibrator (Standard)	–	0.05 mL	–
Distilled water (R3*)	0.05 mL	–	–

Mix and after exactly 1 minute incubation read the initial absorbance for blank A<sub>bl</sub>, sample A<sub>sam</sub> and calibrator A<sub>cal</sub>. Exactly after 2 minutes read the final absorbance of blank A<sub>bl</sub>, sample A<sub>sam</sub> and calibrator (standard) A<sub>cal</sub>. Calculate resulting absorbance as the difference between the final and initial absorbance (ΔA/min). \* see Note.

## CALCULATION

$$\text{Creatinine (mg/dL)} = \frac{\Delta A_{\text{sam}}/\text{min} - \Delta A_{\text{bl}}/\text{min}}{\Delta A_{\text{cal}}/\text{min} - \Delta A_{\text{bl}}/\text{min}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

## Note

The determination of serum/plasma creatinine passed on Jaffe method is not fully specific. Alkaline picrate reacts also with some others substances present in serum matrix. Therefore, to correct the matrix effects, it is recommended to perform two-point calibration with XL MULTICAL calibrator and zero calibrator (reagent R3). Zero calibrator is used as a blank instead of water or saline.

## UNIT CONVERSION

mg/dL × 88.4 = μmol/L

## EXPECTED VALUES

Serum <sup>10</sup>			
Cord	0.6–1.2 mg/dL	18–60 y:	Male 0.9–1.3 mg/dL
Newborn (1–4 d)	0.3–1.0 mg/dL	Female	0.6–1.1 mg/dL
Infant	0.2–0.4 mg/dL	60–90 y:	Male 0.8–1.3 mg/dL
Child	0.3–0.7 mg/dL	Female	0.6–1.2 mg/dL
Adolescent	0.5–1.0 mg/dL	>90 y:	Male 1.0–1.7 mg/dL
		Female	0.6–1.3 mg/dL

## Urine<sup>11</sup>

Infant	8–20 mg/kg/day	Adult:	Male 14–26 mg/kg/day
Child	8–22 mg/kg/day	Female	11–20 mg/kg/day
Adolescent	8–30 mg/kg/day		

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

### Limit of quantification:

Serum/plasma 0.020 mg/dL  
 Urine 0.38 mg/dL  
 Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

### Linearity:

Serum/plasma 30.0 mg/dL  
 Urine 600 mg/dL  
 Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

### Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.34	0.015	1.10	Sample 1	56.2	0.73	1.31
Sample 2	3.34	0.042	1.24	Sample 2	125.2	1.70	1.36

Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.40	0.042	3.02	Sample 1	61.8	1.89	3.05
Sample 2	3.56	0.129	3.62	Sample 2	134.3	4.07	3.03

## Accuracy

Two different validated control materials for serum and urine were used. Determined bias is -1.9 % at the target value 1.38 mg/dL, -7.5 % at the target value 3.99 mg/dL for serum, -12.2 % at the target value 66.9 mg/dL and -14.2 % at the target value 141.4 mg/dL for urine.

## Comparison

A comparison between XL-640 automatic system CREATININE (y) and a commercially available test (x) using 149 samples (serum) gave following results:

Linear regression:  
 $y = 0.955x + 0.131 \text{ mg/dL} \quad r = 0.996$   
 Passing-Bablok<sup>12</sup>:  
 $y = 0.973x + 0.111 \text{ mg/dL} \quad r = 0.992$

## Interferences

Criterion: Recovery within ±10 % of initial value of creatinine concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 10 g/L, bilirubin up to 4 mg/dL, triglycerides up to 850 mg/dL.  
 Cefoxitin, Cephalothin and Cyanokit (Hydroxocobalamin) may cause false-positive results in serum or plasma<sup>13,14,15,16</sup>.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximal allowable absorbance of the reagent blank measured at 505 nm against the distilled water is 0.15.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of creatinine. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

R1  
 UFI: KWFQ-3J78-FE72-MFJX

### Warning

**Hazard statements:**  
 H315 Causes skin irritation.  
 H319 Causes serious eye irritation.

### Precautionary statements:

P280 Wear protective gloves / protective clothing / eye protection.  
 P302 + P352 IF ON SKIN: Wash with plenty of water.  
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### R2, R3

Reagents are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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
- Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Ham erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 10: 391-400, 1886.
- Fabiny DL, Ertinghausen G. Automated reaction-rate method for determination of serum creatinine with the CentrifChem Clin Chem. 17: 696-700, 1971.
- Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 32: 81-85, 1971.
- Myers, G. L., Greg Miller, W., Coresh, J., Fleming, J., Greenberg, N. et al.: Recommendations for Improving Serum Creatinine Measurement, Clin. Chem. 52: 5-18, 2006.
- Fridecký B. Program zlepšovani kvality měření sérového kreatininu, Klin. Biochem. Metab., 14 (35), No.3: 173-176, 2006.
- Bowers, L. D., Wong, E. T. Kinetic serum creatinine assays. II. A critical evaluation and review, Clin. Chem. 26: 555, 1980.
- Guder W, Fonseca-Wollheim W, Ehret W, et al. Die Qualität Diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.
- Newman DJ, Price CP. Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 1204-1270, 1999.
- Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 366-374, 1998.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis CA. Ashwood. ER. Bruns. DE; 5th edition. WB Saunders Comp.. 2012.
- Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. Rifai N. Horvath AR. Wittwer CT; 8th edition. Elsevier. 2019.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
- Kroll MH, Elin RJ. Mechanism of cefoxitin and cephalothin interference with the Jaffé method for creatinine, Clin Chem 29, 12, 2044-2048, 1983.
- Wang P, Commentary on Impact of Cyanokit Administration on Point-of-Care Glucose and Laboratory Test Accuracy, Clin Chem 68, 4, 509-510, 2022.
- Green AJE, Halloran SP, Mould GP, Barbour HM, Pritchard JL, Halloran MJ, Labib M, Interference by newer cephalosporins in current methods for measuring creatinine, Clin Chem 36, 12, 2139-2140, 1990.

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## LITERATÚRA

1. Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Ham erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 10: 391-400, 1886.
2. Fabiny DL, Ertinghausen G. Automated reaction-rate method for determination of serum creatinine with the CentrifChem Clin Chem. 17: 696-700, 1971.
3. Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 32: 81-85, 1971.
4. Myers, G. L., Greg Miller, W., Coresh, J., Fleming, J., Greenberg, N. et al.: Recommendations for Improving Serum Creatinine Measurement, Clin. Chem. 52: 5-18, 2006.
5. Fridecký B. Program zlepšování kvality měření sérového kreatininu, Klin. Biochem. Metab., 14 (35), No.3: 173-176, 2006.
6. Bowers, L. D., Wong, E. T. Kinetic serum creatinine assays. II. A critical evaluation and review, Clin. Chem. 26: 555, 1980.
7. Guder W, Fonseca-Wollheim W, Ehret W, et al. Die Qualität Diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.
8. Newman DJ, Price CP. Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 1204-1270, 1999.
9. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 366-374, 1998.
10. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis CA. Ashwood. ER. Bruns. DE; 5th edition. WB Saunders Comp.. 2012.
11. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. Rifai N. Horvath AR. Wittwer CT; 8th edition. Elsevier. 2019.
12. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
13. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
14. Kroll MH, Elin RJ. Mechanism of cefoxitin and cephalothin interference with the Jaffé method for creatinine, Clin Chem 29, 12, 2044-2048, 1983.
15. Wang P, Commentary on Impact of Cyanokit Administration on Point-of-Care Glucose and Laboratory Test Accuracy, Clin Chem 68, 4, 509-510, 2022.
16. Green AJE, Halloran SP, Mould GP, Barbour HM, Pritchard JL, Halloran MJ, Labib M, Interference by newer cephalosporins in current methods for measuring creatinine, Clin Chem 36, 12, 2139-2140, 1990.

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