GAMMAGLUTAMYLTRANSFERASE

Cat. No.	Pack Name	Packaging (Content)
BLT00023	GGT 100	R1: 4 x 20 mL, R2: 1 x 20 mL instruction for use
BLT00024	GGT 250	R1: 4 x 50 mL, R2: 1 x 50 mL instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of gamma glutamyl transferase in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases and damage, detection of jaundice, cholangitis, cholecystitis. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Although gamma glutamyl transferase (GGT) is present in a variety of tissues, the serum enzyme appears to be primarily from the hepatobiliary system. Consequently, GGT is elevated in all forms of liver disease or damage. It is clinically useful in detecting obstructive jaundice, cholangitis and cholecystitis. Elevated levels are also observed with drug use (alcohol, sedatives, anticonvulsants and tranquilizers).

PRINCIPLE

The products are based on the principle of photometric determination according to Persijn & van der Slik¹ and standardized to methods recommended by IFCC² and Szasz¹. The performance claims and data presented here are independent from the standardization. GGT present in the sample catalyzes the transfer of the glutamyl group from the substrate L-γ-glutamyl-3-carboxy-4-nitroanilide to glycylglycine forming L-γ-glutamylglycylglycine and 5-amino-2-nitrobenzoate.

The rate of formation of 5-amino-2-nitrobenzoate is proportional to the activity of GGT present in the sample and can be measured kinetically at 400–420 nm.

REAGENT DESCRIPTION AND COMPOSITION

R1	R2
Tris buffer, pH 8.25	L-γ-glutamyl-3-carboxy-4-nitroanilide 20 mmol/L
Glycylglycine	125 mmol/L

COMPOSITION OF REACTION MIXTURE

Tris buffer	81 mmol/L
Glycylglycine	81 mmol/L
L-γ-glutamyl-3-carboxy-4-nitroanilide	3.7 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 400–420 nm may be used, general laboratory equipment.

- XL MULTICAL 4x3, Cat. No. XSYS0034
- XL MULTICAL 10x3, Cat. No. XSYS0122
- ERBA NORM 4x5, Cat. No. BLT00080
- ERBA NORM 10x5, Cat. No. XSYS0123
- ERBA PATH 4x5, Cat. No. BLT00081
- ERBA PATH 10x5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method – substrate start

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

Monoreagent method – sample start

Stability of working reagent: 2 weeks at 20–25 °C in dark
6 weeks at 2–8 °C in dark

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.
Plasma: Li-heparin and K₂-EDTA plasma.
The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.
Centrifuge samples containing precipitates before performing the assay.
See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma^{3,4}: 7 days at 15–25 °C
7 days at 2–8 °C
12 months at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
2 point calibration (blank and calibrator); distilled water is recommended as blank
Calibration frequency: it is recommended to do a calibration
• after reagent lot change
• as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC² or Szasz¹ recommended methods.

ASSAY PROCEDURE

Wavelength: 405 (400–420) nm
Cuvette: 1 cm

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.100 mL
Calibrator	–	0.100 mL	–
Distilled water	0.100 mL	–	–

Mix and after 1 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.100 mL
Calibrator	–	0.100 mL	–
Distilled water	0.100 mL	–	–

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

CALCULATION

$$1. \text{GGT (U/L)} = \frac{\Delta A_{\text{sam}}/\text{min.}}{\Delta A_{\text{cal}}/\text{min.}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \text{Using factor: GGT (U/L)} = f \times \Delta A/\text{min} \quad f = \text{factor}$$

Substrate Start:

standardized against	Szasz	IFCC
factors at 405 nm and 37 °C	1421	1606

Sample Start:

standardized against	Szasz	IFCC
factors at 405 nm and 37 °C	1477	1669

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic	Reaction direction	Increasing
Wavelength (nm)	405	Normal Low U/L	0
Sample Volume (μL)	50/100	Normal High U/L	55
Working Reagent Volume (μL)	500/1000	Linearity Low U/L	1.95
Lag time (sec.)	60	Linearity High U/L	1114
Kinetic interval (sec.)	60	Blank with	Water
No. of readings	3	Absorbance limit (max.)	1.5
Kinetic factor	1477	Units	U/L
Reaction temperature (°C)	37		

UNIT CONVERSION

U/L × 0.0167 = μkat/L

EXPECTED VALUES⁵

At 37 °C

Females: <38 U/L

Males: <55 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 1.95 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 1114 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	42.2	0.51	1.22
Sample 2	189.5	1.57	0.83

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	41.9	1.25	2.98
Sample 2	186.8	3.27	1.75

Accuracy

Two different validated control materials were used. Determined bias is -1.8 % at the target value 130.2 U/L and -0.7 % at the target value 190.5 U/L.

Comparison

A comparison between XL-640 automatic system GGT (y) and a commercially available test (x) using 54 samples gave following results:

Linear regression: y = 0.993x - 3.35 U/L r = 0.998

Passing-Bablok⁶: y = 0.993x - 1.71 U/L r = 0.991

Interferences

Criterion: Recovery within ±10 % of initial value of GGT activity in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 4.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁷.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 400–420 nm against the distilled water is 1.5.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of GGT. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
- Persijn JP, van der Slik W. A new Method for the Determination of γ -Glutamyltransferase. J Clin Chem Clin Biochem 1976;4: 421.
- Schumann G, Bonora R, Ceriottiet F et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 6. Reference Procedure for the Measurement of Catalytic Activity Concentrations of gamma-glutamyltransferase. Clin Chem Lab Med 2002; 40 (7): 734-738.
- Szasz G. Methods of Enzymatic Analysis. 2nd English ed. New York, Academic Press, Inc. 1974; 717.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company 1995; 286.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.

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LITERATÚRA

1. Persijn JP, van der Slik W. A new Method for the Determination of γ -Glutamyltransferase. J Clin Chem Clin Biochem 1976;4: 421.
2. Schumann G, Bonora R, Ceriottiet F et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 6. Reference Procedure for the Measurement of Catalytic Activity Concentrations of gamma-glutamyltransferase. Clin Chem Lab Med 2002; 40 (7): 734-738.
3. Szasz G. Methods of Enzymatic Analysis. 2nd English ed. New York, Academic Press, Inc. 1974; 717.
4. Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company 1995; 286.
5. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
6. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
7. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.

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