# HDL CHOLESTEROL DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00030	HDL C DIRECT 240	R1: 3 × 60 mL, R2: 3 × 20 mL, instruction for use

## INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of HDL cholesterol in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of coronary heart disease. For professional use in clinical laboratory only.

## CLINICAL SIGNIFICANCE


High-density lipoproteins (HDL) compose one of the major classes of plasma lipoproteins. They are synthesized in liver as complexes of apolipoprotein and phospholipid and are capable of picking up cholesterol and carrying it from arteries to the liver, where the cholesterol is converted to bile acids and excreted into the intestine.

An inverse relationship between HDL Cholesterol (HDL C) levels in serum and the incidence/prevalence of coronary heart disease (CHD) has been demonstrated in a number of epidemiological studies. The importance of HDL C as a risk factor for CHD is now recognized.

Accurate measurement of HDL C is of vital importance when assessing patient's risk for CHD.

## PRINCIPLE

Antihuman  $\beta$ -lipoproteins antibody in the R1 binds to lipoproteins (LDL, VLDL and chylomicrons) other than HDL. The antigen-antibody complexes formed blocked enzyme reaction when R2 is added. Cholesterol esterase (CHE) and cholesterol oxidase (CHO) in R2 react only with HDL cholesterol. Hydrogen peroxide ( $H_2O_2$ ) produced by the enzyme reactions with HDL cholesterol yields a blue colour complex upon oxidative condensation of F-DAOS [N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-3,5-dimethoxy-4-fluoroaniline, sodium salt] and 4-aminoantipyrine (4-AAP) in the presence of peroxidase (POD)<sup>1,2,3</sup>.

The absorbance of the blue dye at 600 nm is proportional to the HDL cholesterol concentration in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

R1		
Good's buffer (pH 7.0)	30	mmol/L
4-aminoantipyrine	0.9	mmol/L
Peroxidase	2.4	kU/L
Ascorbate oxidase	2.7	kU/L
Antihuman $\beta$ -lipoproteins antibody		

R2		
Good's buffer (pH 7.0)	30	mmol/L
Cholesterol esterase	4	kU/L
Cholesterol oxidase	2.0	kU/L
F-DAOS	0.8	mmol/L

## COMPOSITION OF REACTION MIXTURE

Good's buffer (pH 7.0)	29.8	mmol/L
4-aminoantipyrine	0.67	mmol/L
F-DAOS	0.2	mmol/L
Peroxidase	1.8	kU/L
Ascorbate oxidase	2.0	kU/L
Cholesterol esterase	1.0	kU/L
Cholesterol oxidase	5.0	kU/L
Antihuman $\beta$ -lipoproteins antibody		

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of  $37 \pm 0.5$  °C that is capable of reading absorbance at 600/700 nm may be used, general laboratory equipment.

HDL C CAL, Cat. No. BLT00072  
 ERBA NORM 4x5, Cat. No. BLT00080  
 ERBA NORM 10x5, Cat. No. XSYS0123  
 ERBA PATH 4x5, Cat. No. BLT00081  
 ERBA PATH 10x5, Cat. No. XSYS0124

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma.

Fasting and non-fasting samples can be used<sup>4,5</sup>.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma <sup>6</sup> :	2 days at	20–25 °C
	7 days at	4–8 °C
	3 months at	-20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with HDL C CAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, HDL/LDL calibrator and controls ERBA NORM and ERBA PATH have been standardized against the US CDC Reference Method UCMH.

## ASSAY PROCEDURE

Wavelength: 600/700 nm

Cuvette: 1 cm

	Reagent blank	Calibrator	Sample
Reagent 1	0.900 mL	0.900 mL	0.900 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and after 5 min. incubation read the initial absorbance for blank  $A_{bl}$ , sample  $A_{sam}$  and calibrator  $A_{cal}$ .

	0.300 mL	0.300 mL	0.300 mL
Reagent 2			

Mix and after 5 min. incubation read the final absorbance for blank  $A_{bl}$ , sample  $A_{sam}$  and calibrator  $A_{cal}$ . Calculate resulting absorbance like the difference between the final and initial absorbance  $A = (A_{FINAL} - A_{INITIAL})$ .

## CALCULATION

$$\text{HDL cholesterol (mg/dL)} = \frac{A_{sam} - A_{bl}}{A_{cal} - A_{bl}} \times C_{cal}$$

$C_{cal}$  = calibrator concentration

## UNIT CONVERSION

mg/dL × 0.026 = mmol/L

## EXPECTED VALUES<sup>15</sup>

In serum:	Male	Female
5–9 y	38–75	36–73 mg/dL
10–14 y	37–74	37–70 mg/dL
15–19 y	30–63	35–74 mg/dL
20–24 y	30–63	33–79 mg/dL
25–29 y	31–63	37–83 mg/dL
30–34 y	28–63	36–77 mg/dL
35–39 y	29–62	34–82 mg/dL
40–44 y	27–67	34–88 mg/dL
45–49 y	30–64	34–87 mg/dL
50–54 y	28–63	37–92 mg/dL
55–59 y	28–71	37–91 mg/dL
60–64 y	30–74	38–92 mg/dL
65–69 y	30–75	35–96 mg/dL
>69 y	31–75	33–92 mg/dL

## Adult Treatment Panel III classification<sup>16</sup>:

Low <40 mg/dL  
 High >59 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.48 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

**Linearity:** 208 mg/dL

Linearity is the highest measured activity with recovery within  $\pm 10$  % from theoretical value.

## Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	41.7	0.29	0.70
Sample 2	79.9	0.47	0.59

Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	41.1	1.05	2.55
Sample 2	79.6	1.84	2.31

## Accuracy

Two different validated control materials were used. Determined bias is -3.7 % at the target value 55.4 mg/dL and -6.9 % at the target value 83.8 mg/dL.

## Comparison

A comparison between XL-640 automatic system HDL CHOLESTEROL DIRECT (y) and a commercially available test (x) using 130 samples gave following results:

Linear regression:  
 $y = 1.047x - 4.550$  mg/dL  $r = 0.982$

Passing-Bablok<sup>17</sup>:  
 $y = 1.047x - 4.654$  mg/dL  $r = 0.985$

## Interferences

**Criterion:** Recovery within  $\pm 10$  % of initial value of cholesterol concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

N-acetylcysteine, Metamizole and Acetaminofen (Paracetamol) including its metabolite N-acetyl-p-benzoquinone imine may cause false-negative results<sup>18,19</sup>.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 600 nm against the distilled water is 0.3.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of HDL cholesterol. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Health protection

Used animal-source materials are accompanied by an official state veterinarian health certificate. The animal components were collected from healthy animals under the care of a registered veterinarian and were apparently free from infectious and contagious diseases and injurious parasites. However, this material should be handled as though capable of transmitting infectious diseases.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### R1

UF1: PAWG-4WA6-3J5S-PM0Q

### Warning

Contains: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one

### Hazard statement:

H317 May cause an allergic skin reaction.

### Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

### R2

EUH208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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


1. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Wiebe DA, Warnick GR. Measurement of high-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 127-144, 1997.
3. Nauck M, Maerz W, Wieland H. New immunoseparationbased homogenous assay for HDL-cholesterol compared with three homogenous and two heterogeneous methods for HDL-cholesterol, Clin Chem 44: 1443-1451, 1998.
4. Sidhu D, Naugler C. Fasting time and lipid levels in a community-based population: A cross sectional study. Arch. Intern. Med. Dec 10, 172(22): 1707-1710, 2012.
5. Ontario Community Laboratory Guideline for Adult Lipid Testing (CLP017) 2013.
6. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
7. Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
8. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
9. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
10. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
11. National Institutes of Health publication No. 93-3095, September, 1993.
12. Williams P, et al., High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
13. Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
14. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
15. Anonymous, Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), 285, 2486-2497, JAMA 2001.
16. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
17. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
18. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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## LITERATÚRA

1. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Wiebe DA, Warnick GR. Measurement of high-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 127-144, 1997.
3. Nauck M, Maerz W, Wieland H. New immunoseparationbased homogenous assay for HDL-cholesterol compared with three homogenous and two heterogeneous methods for HDL-cholesterol. Clin Chem 44: 1443-1451, 1998.
4. Sidhu D, Naugler C. Fasting time and lipid levels in a community-based population: A cross sectional study. Arch. Intern. Med. Dec 10, 172(22): 1707-1710, 2012.
5. Ontario Community Laboratory Guideline for Adult Lipid Testing (CLP017) 2013.
6. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
7. Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
8. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
9. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
10. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
11. National Institutes of Health publication No. 93-3095, September, 1993.
12. Williams P, et al., High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
13. Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
14. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
15. Anonymous, Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), 285, 2486-2497, JAMA 2001.
16. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
17. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
18. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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