# LACTATE DEHYDROGENASE-L

Cat. No.	Pack Name	Packaging (Content)
BLT00038	LDH-L 100	R1: 4 × 20 mL, R2: 1 × 20 mL, instruction for use
BLT00039	LDH-L 250	R1: 4 × 50 mL, R2: 1 × 50 mL, instruction for use

## INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of lactate dehydrogenase in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases, myocardial infarction, renal damage, anemias. For professional use in clinical laboratory only.

## CLINICAL SIGNIFICANCE

Lactate dehydrogenase (LDH) is an enzyme, consisting of five different isoenzymes which catalyze the interconversion of L-lactate and pyruvate. LDH is present in the cytoplasm of all human tissues with higher concentrations in liver, heart and skeletal muscle, and lower values in erythrocytes, pancreas, kidney and stomach. Elevated serum levels of LDH have been observed in a variety of disease states. The highest levels are seen in patients with megaloblastic anemia, disseminated carcinoma and shock. Moderate increases occur in muscular disorders, nephrotic syndrome and cirrhosis. Mild increases in LDH activity have been reported in cases of myocardial or pulmonary infarction, leukemia, hemolytic anemia and non-viral hepatitis.

## PRINCIPLE

This assay follows the recommendations of the IFCC<sup>1,2</sup>. LDH catalyzes the conversion of L-lactate to pyruvate, NAD<sup>+</sup> is reduced to NADH in the process.

The rate of formation of NADH is proportional to the activity of LDH present in the sample and can be measured kinetically at 340 nm.

## REAGENT DESCRIPTION AND COMPOSITION

R1	R2
N-Methyl-D-glucamine, pH 9.40	NAD <sup>+</sup>
L-lactate	

COMPOSITION OF REACTION MIXTURE	50 mmol/L
N-Methyl-D-glucamine	320 mmol/L
L-lactate	49 mmol/L
NAD <sup>+</sup>	10 mmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 340 nm may be used, general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034  
 XL MULTICAL 10×3, Cat. No. XSYS0122  
 ERBA NORM 4×5, Cat. No. BLT00080  
 ERBA NORM 10×5, Cat. No. XSYS0123  
 ERBA PATH 4×5, Cat. No. BLT00081  
 ERBA PATH 10×5, Cat. No. XSYS0124

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light – particularly R2, and without any contamination.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable. Serum.

Plasma: Li-heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay. See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma <sup>2</sup> :	7 days at	15–25 °C
	4 days at	2–8 °C
	6 weeks at	-20 °C

In connection with certain diseases (e.g. hepatopathy, diseases of skeletal muscle, malignant tumors), the LDH-4 and LDH-5 isoenzyme portions are increased and unstable in cooled and frozen samples; this may lead to an incorrect LDH value in samples collected from patients suffering from such diseases. Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL is recommended. 2 point calibration (blank and calibrator); distilled water is recommended as blank. Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method<sup>2</sup>.

## ASSAY PROCEDURE

Wavelength: 340 nm  
 Cuvette: 1 cm

## Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.020 mL
Calibrator	–	0.020 mL	–
Distilled water	0.020 mL	–	–

Mix and after 1–5 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

## CALCULATION

$$1. \text{LDH (U/L)} = \frac{\Delta A_{\text{sam}}/\text{min.}}{\Delta A_{\text{cal}}/\text{min.}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \text{Using factor (f):} \quad \text{LDH (U/L)} = f \times \Delta A/\text{min} \quad f = 10080 \text{ (at 340 nm)}$$

## UNIT CONVERSION

U/L × 0.0167 = μkat/L

## EXPECTED VALUES<sup>3</sup>

At 37 °C

Adults 125–250 U/L

Children (24 mo–12 y) 180–360 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 14.8 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV < 20 % (n = 30).

**Linearity:** 1516 U/L

Linearity is the highest measured activity with recovery within ± 10 % from theoretical value.

## Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	179.3	1.76	0.98
Sample 2	350.8	5.07	1.44

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	179.4	6.50	3.74
Sample 2	348.1	10.37	2.98

## Accuracy

Two different validated control materials were used. Determined bias is -2.1 % at the target value 260.6 U/L and -2.7 % at the target value 442.7 U/L.

## Comparison

A comparison between XL-640 automatic system LDH-L (y) and a commercially available test (x) using 48 samples gave following results:

Linear regression:  
 $y = 1.036x - 6.26 \text{ U/L} \quad r = 0.988$

Passing-Bablok<sup>3</sup>:  
 $y = 1.019x - 2.44 \text{ U/L} \quad r = 0.967$

## Interferences

Criterion: Recovery within ± 10 % of initial value of LDH activity in the sample without interfering substance. Following substances do not interfere: haemoglobin up to 1 g/l, bilirubin up to 40 mg/dl, triglycerides up to 850 mg/dl. Significant hemolysis may increase LDH concentration because of high levels of LDH in the erythrocytes.

Drugs: No interference was found at therapeutic concentrations using common drug panels<sup>10</sup>.

## Limitations:

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of LDH. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents of the kit are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Van der Heiden C, Bais R, Gerhardt W, et al. Approved recommendation on IFCC methods for the measurement of catalytic concentration of enzymes. Part 8. IFCC method for lactate dehydrogenase. Eur J Clin Chem Clin Biochem 1994;32:639-655.
2. Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 3. Reference Procedures for the Measurement of Catalytic Concentrations of Lactate Dehydrogenase. Clin Chem Lab Med 2002;40(6):643-648.
3. Thomas L, ed. Labor und Diagnose, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft 1992.
4. Moss DW, Henderson AR, Kachmar JF. Enzymes. In: Tietz NW, ed. Fundamentals of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders 1987;346-421.
5. Zimmerman HJ, Henry JB In: Henry JB, ed. Clinical Diagnosis and Management by Laboratory Methods. 17th ed. Philadelphia, PA: WB Saunders 1984;251-282.
6. Moss DW, Henderson AR. Clinical enzymology In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company;1999.617-721.
7. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
8. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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