
# LDL CHOLESTEROL DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00042	LDL C DIRECT 240	R1: 3 × 60 mL, R2: 3 × 20 mL, instruction for use

## INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of LDL cholesterol in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of coronary heart disease and lipoprotein metabolism disorders. For professional use in clinical laboratory only.


## CLINICAL SIGNIFICANCE

Low Density Lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on triglyceride-rich Very Low Density Lipoproteins (VLDLs). Specific LDL receptors exist to facilitate the elimination of LDL from plasma by liver parenchymal cells. It has been shown that most of the cholesterol stored in atherosclerotic plaques originates from LDL. For this reason the LDL Cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL Cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

## PRINCIPLE

When a sample is mixed with Reagent 1, the protective reagent binds to LDL and protects LDL from enzyme reactions. Cholesterol esterase (CHE) and cholesterol oxidase (CHO) reacts with non-LDL lipoproteins (chylomicrons, VLDL, HDL). Hydrogen peroxide produced by the enzyme reactions with non-LDL cholesterol is decomposed by catalase in Reagent 1. When Reagent 2 is added, the protective reagent is removed from LDL and catalase is inactivated by sodium azide. In this second process, CHE and CHO react only with LDL cholesterol. Hydrogen peroxide produced by the enzyme reactions with LDL cholesterol yields a colour complex upon oxidative condensation with H-DAOS [N-(2-hydroxy-3-sulfo-propyl)-3,5-dimethoxyaniline] and 4-aminoantipyrine in the presence of peroxidase (POD)<sup>1,2</sup>.

The absorbance of the blue dye at 600 nm is proportional to the LDL cholesterol concentration in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

R1	
Good's buffer (pH 6.8)	25 mmol/L
H-DAOS	0.64 mmol/L
Cholesterol esterase	5 KU/L
Cholesterol oxidase	5 KU/L
Catalase	1 KU/L

R2	
Good's buffer (pH 7.0)	25 mmol/L
4-aminoantipyrine	3.4 mmol/L
Peroxidase	20 KU/L
Sodium azide	<0.1 %

COMPOSITION OF REACTION MIXTURE	
Good's buffer	25 mmol/L
4-aminoantipyrine	0.84 mmol/L
H-DAOS	0.47 mmol/L
Cholesterol esterase	3.7 KU/L
Cholesterol oxidase	3.7 KU/L
Catalase	0.7 KU/L
Peroxidase	5 KU/L
Sodium azide	<0.025 %

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 600/700 nm may be used, general laboratory equipment.  
 LDL C CAL, Cat. No. BLT00073  
 ERBA NORM 4x5, Cat. No. BLT00080  
 ERBA NORM 10x5, Cat. No. XSYS0123  
 ERBA PATH 4x5, Cat. No. BLT00081  
 ERBA PATH 10x5, Cat. No. XSYS0124

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. After opening, reagents are stable for 60 days at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.  
 For specimen collection and preparation only use suitable tubes or collection containers.  
 Only the specimens listed below were tested and found acceptable.

Serum.  
 Plasma: Li-heparin plasma.  
 Fasting and non-fasting samples can be used.  
 The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.  
 Centrifuge samples containing precipitates before performing the assay.  
 See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma <sup>2</sup> :	
	1 days at 20–25 °C
	7 days at 4–8 °C
	3 months at -20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with LDL C CAL is recommended.  
 2 point calibration (blank and calibrator); distilled water is recommended as blank  
 Calibration frequency: it is recommended to do a calibration  
 • after reagent lot change  
 • as required by internal quality control procedures

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.  
 The control intervals and limits should be adapted according to each individual laboratory's require-

ments. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, LDL C CAL calibrator and controls ERBA NORM and ERBA PATH have been standardized against the US CDC Reference Method BQ.

## ASSAY PROCEDURE

Wavelength: 600/700 nm  
 Cuvette: 1 cm

	Reagent blank	Calibrator	Sample
Reagent 1	0.900 mL	0.900 mL	0.900 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and after 5 min. incubation read the initial absorbance for blank  $A_{bl}$ , sample  $A_{sam}$  and calibrator  $A_{cal}$ .  
 Then add:

Reagent 2	0.300 mL	0.300 mL	0.300 mL
-----------	----------	----------	----------

Mix and after 5 min. incubation read the final absorbance for blank  $A_{bl}$ , sample  $A_{sam}$  and calibrator  $A_{cal}$ . Calculate resulting absorbance like the difference between the final and initial absorbance  $A = (A_{FINAL} - A_{INITIAL})$ .

## CALCULATION

$$LDL \text{ cholesterol (mg/dL)} = \frac{A_{sam} - A_{bl}}{A_{cal} - A_{bl}} \times C_{cal} \quad C_{cal} = \text{calibrator concentration}$$

## UNIT CONVERSION

mg/dL × 0.025 = mmol/L

## EXPECTED VALUES<sup>10</sup>

In serum:	Male	Female
5–9 y	63–129	68–140 mg/dL
10–14 y	64–133	68–136 mg/dL
15–19 y	62–130	59–137 mg/dL
20–24 y	66–147	57–159 mg/dL
25–29 y	70–165	71–164 mg/dL
30–34 y	78–185	70–156 mg/dL
35–39 y	81–189	75–172 mg/dL
40–44 y	87–186	74–174 mg/dL
45–49 y	97–202	79–186 mg/dL
50–54 y	89–197	88–201 mg/dL
55–59 y	88–203	89–210 mg/dL
60–64 y	83–210	100–224 mg/dL
65–69 y	98–210	92–221 mg/dL
>69 y	88–186	96–206 mg/dL

## Coronary heart disease risk:

Optimal	<100 mg/dL
Near/above optimal	100–129 mg/dL
Borderline high	130–159 mg/dL
High	160–189 mg/dL
Very high	>189 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.65 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

**Linearity:** 413 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

## Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	92.7	0.52	0.57
Sample 2	163.2	0.85	0.52

Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	94.3	1.86	1.98
Sample 2	165.1	3.45	2.09

## Accuracy

Two different validated control materials were used. Determined bias is 2.9 % at the target value 128.8 mg/dL and -1.7 % at the target value 188.5 mg/dL.

## Comparison

A comparison between XL-640 automatic system LDL CHOLESTEROL DIRECT (y) and a commercially available test (x) using 140 samples gave following results:

Linear regression:  
 $y = 0.947x + 10.719 \text{ mg/dL} \quad r = 0.976$

Passing-Bablok<sup>11</sup>:  
 $y = 0.962x + 8.385 \text{ mg/dL} \quad r = 0.974$

## Interferences

Criterion: Recovery within ±10 % of initial value of LDL cholesterol concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 36 mg/dL, triglycerides up to 850 mg/dL.

N-acetylcysteine, Metamizole and Acetaminofen (Paracetamol) including its metabolite N-acetyl-p-benzoquinone imine may cause false-negative results<sup>12,13</sup>.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 600 nm against the distilled water is 0.3.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of LDL cholesterol. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

**R1**  
 EUH208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.  
 EUH210 Safety data sheet available on request.

**R2**  
 EUH208 Contains Alcohols, C11-15-secondary, ethoxylated. May produce an allergic reaction.  
 EUH210 Safety data sheet available on request.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl Lab 31: 325-330, 1985.
3. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
4. Bachorik PS. Measurement of low-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press, 145-160, 1997.
5. Rifai N, Warnick GR, McNamara JR, et al. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 38: 150-160, 1992.
6. Rifai N. and Warnick GR, Dominiczak MHEd, Hand book of Lipoprotein Testing, AACC Press, Washington, DC, USA, 1997.
7. Friewald, W. T., Levy, R. I. and Frederickson, D. S.: Estimation of the concentration of low density lipoprotein cholesterol in plasma without use of ultracentrifuge. Clin. Chem. 18, 449–502, 1972.
8. The Expert Panel. Report of the National Cholesterol Education Program Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults. Arch. Intern. Med., 148, 36–69, 1988.
9. The Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults. Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II). JAMA 269, 3015–3023, 1993.
10. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
11. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
12. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
13. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100-104, 2016.

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## LITERATÚRA

1. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Arztl Lab 31: 325-330, 1985.
3. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
4. Bachorik PS. Measurement of low-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press, 145-160, 1997.
5. Rifai N, Warnick GR, McNamara JR, et al. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 38: 150-160, 1992.
6. Rifai N. and Warnick GR, Dominiczak MHEd, Hand book of Lipoprotein Testing, AACC Press, Washington, DC, USA, 1997.
7. Friedwald, W. T., Levy, R. I. and Frederickson, D. S.: Estimation of the concentration of low density lipoprotein cholesterol in plasma without use of ultracentrifuge. Clin. Chem. 18, 449–502, 1972.
8. The Expert Panel. Report of the National Cholesterol Education Program Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults. Arch. Intern. Med., 148, 36–69, 1988.
9. The Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults. Summary of the Second Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II). JAMA 269, 3015–3023, 1993.
10. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
11. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
12. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
13. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100-104, 2016.

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