MAGNESIUM

Cat. No.	Pack Name	Packaging (Content)
BLT00049	MG 250	R1: 1 × 250 mL, R2 standard: 1 × 5 mL, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of magnesium in human serum, plasma and urine on various automatic systems. Intended for screening, monitoring and diagnosis of impairment of neuromuscular functions, gastrointestinal and cardiac symptoms, treatment monitoring. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE


Magnesium is an essential nutrient which is involved in many biochemical functions. It has a structural role in nucleic acids and ribosomal particles, required as an activator for many enzymes and has a role in energy producing oxidative phosphorylation.

The normal body contains between 21 to 28 g magnesium, more than 50 % of which is complexed with calcium and phosphate in bones. Only approximately 1 % of the total magnesium is found in the extracellular fluid; hence, it tends to enter and leave cells under the same conditions as potassium. Approximately 35 % of plasma magnesium is protein-bound, mainly to albumin, and therefore changes in albumin concentration may affect magnesium.

Hypomagnesemia results in the impairment of neuromuscular function and may develop in severe prolonged diarrhea, malabsorption syndromes, hyperaldosteronism, and diuretic therapy. Hypermagnesemia is seen in renal glomerular failure and diabetic coma.

PRINCIPLE

Magnesium reacts with Xylidyl Blue to form a purple colored complex in alkaline solution. Interference by calcium is prevented by the use of glycoethylerdiaminetetraacetic acid (GEDTA)^{2,3}.

The resulting absorbance at 500–550 nm is proportional to the magnesium concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2 standard	
Ethanolamine (pH 11.0)	1.0 µmol/L	Magnesium (II)	See bottle label
Xylidyl Blue	110 µmol/L		
GEDTA	60 µmol/L		

COMPOSITION OF REACTION MIXTURE

Ethanolamine (pH 11.0)	0.99 µmol/L
Xylidyl Blue	109 µmol/L
GEDTA	59 µmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 500–550 nm may be used, general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagent is stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin. Chelating anticoagulants such as EDTA, fluoride and oxalate must be avoided.

Urine: Acidify urine with some drops of conc. HCl to pH 3–4 to prevent precipitation of magnesium ammonium phosphate. Collect urine samples in metal-free container². Dilute urine samples using redistilled water in 1 + 4 ration and multiply results by 5.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma ^a :	7 days at	15–25 °C
	7 days at	2–8 °C
	1 year at	-20 °C

Stability in urine ^a :	3 days at	15–25 °C
	3 days at	2–8 °C
	1 year at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL, R2 standard and controls ERBA NORM and ERBA PATH have been standardized against the SRM 956 reference material.

ASSAY PROCEDURE

Wavelength: 520 (500–550) nm

Cuvette: 1 cm

	Reagent blank	Calibrator (Standard)	Sample
Reagent 1	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and incubate 5 min. at 37 °C. Measure absorbance of the sample A_{sam} and calibrator (standard) A_{cal} against reagent blank. The coloration is stable during one hour.

CALCULATION

$$\text{Magnesium (mg/dL)} = \frac{A_{\text{sam}}}{A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator (standard) concentration}$$

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point	Normal High (mg/dL)	2.6
Wavelength 1 (nm)	505	Linearity Low (mg/dL)	0.11
Sample Volume (µL)	5/10	Linearity High (mg/dL)	6.8
Reagent Volume (µL)	500/1000	Concentration of Standard	See bottle label
Incubation time (min.)	5	Blank with	Reagent
Reaction temperature (°C)	37	Absorbance limit (max.)	0.7
Reaction direction	Increasing	Units	mg/dL
Normal Low (mg/dL)	1.6		

UNIT CONVERSION

mg/dL × 0.4114 = mmol/L

EXPECTED VALUES⁹

Serum or plasma:

Newborn, 2–4 days	1.5–2.2 mg/dL
5 months–6 years	1.7–2.3 mg/dL
6–12 years	1.7–2.1 mg/dL
>12 years	1.6–2.6 mg/dL

Urine: 12–291 mg / 24 h

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

Serum / plasma	0.11 mg/dL
Urine	0.45 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

Linearity:

Serum / plasma	6.80 mg/dL
Urine	34.0 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.96	0.042	2.14	Sample 1	2.00	0.068	3.39
Sample 2	3.61	0.033	0.93	Sample 2	3.61	0.065	1.80

Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	2.92	0.050	1.71	Sample 1	2.05	0.045	2.17
Sample 2	5.21	0.040	0.77	Sample 2	4.21	0.141	3.36

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is 0.5 % at the target value 1.90 mg/dL, -5.1 % at the target value 5.16 mg/dL for serum, 1.7 % at the target value 6.32 mg/dL and -1.3 % at the target value 12.0 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system MAGNESIUM (y) and a commercially available test (x) using 150 samples (serum) gave following results:

Linear regression:

$$y = 0.917x + 0.222 \text{ mg/dL} \quad r = 0.973$$

Passing-Bablok¹⁰:

$$y = 0.929x + 0.198 \text{ mg/dL} \quad r = 0.965$$

Interferences

Criterion: Recovery within ±10 % of initial value of magnesium concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 9 g/L, bilirubin up to 24 mg/dL, triglycerides up to 850 mg/dL.

Drugs:

Serum: No interference was found at therapeutic concentrations using common drug panels¹¹.

Urine: No interference was found at therapeutic concentrations using common drug panels except Acetylcysteine, Gentamicin and Methyldopa¹¹.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 505 nm against the distilled water is 0.7.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of magnesium. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFF: MXWU-1WQR-TJ5W-RXPE

Danger

Contains: 2-aminoethanol

Hazard Statements:

H315 Causes skin irritation.

H318 Causes serious eye damage.

Precautionary Statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water and soap.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

R2

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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7. Sitzmann FC. Normalwerte. München: Hans Marseille Verlag GmbH, 166, 1986.