# TOTAL PROTEIN

Cat. No.	Pack Name	Packaging (Content)
BLT00054	TP 250	R1: 5 × 50 mL, R2 standard: 1 × 3 mL instruction for use
BLT00055	TP 500	R1: 2 × 250 mL, R2 standard: 1 × 3 mL instruction for use

## INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of total protein in human serum and plasma on various automatic systems. Intended for monitoring of gross changes in protein levels caused by various disease states. For professional use in clinical laboratories only.

## CLINICAL SIGNIFICANCE


Plasma proteins derive primarily from synthesis in the liver, plasma cells, lymph nodes, spleen, and bone marrow. In disease states both the total plasma protein level and the ratio of the individual fractions may be dramatically altered from their normal values.

Hypoproteinaemia may be caused by such conditions as nephrotic syndrome, extensive bleeding, sprue (deficient protein absorption), severe burns, salt retention syndromes, and Kwashiorkor (acute protein starvation).

Hyperproteinaemia may be observed in cases of severe dehydration and disease states such as multiple myeloma. Changes in the proportions of the plasma proteins may occur in one or several of the protein fractions and often without alterations in the quantity of the total protein. The A/G ratio has commonly been used as an index of the distribution between the albumin and globulin fractions. This ratio can be significantly altered in such conditions as cirrhosis of the liver, glomerulonephritis, nephrotic syndrome, acute hepatitis, lupus erythematosus, and in some acute and chronic infections.

## PRINCIPLE

Biuret method<sup>1,2,3,4,5</sup>. The peptide bonds of protein react with copper (II) ions in alkaline solution to form a blue violet ion complex, the so called biuret reaction, each copper ion complexing with 5 or 6 peptide bonds. Tartrate is added as a stabiliser whilst iodide is used to prevent auto-reduction of the alkaline copper complex.

Absorbance of the blue violet complex measured at 520–560 nm is proportional to the protein concentration in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

R1	R2 standard	
Copper sulphate	12 mmol/L	See bottle label
Potassium sodium tartrate	31.9 mmol/L	
Potassium iodide	30.1 mmol/L	
Sodium hydroxide	600 mmol/L	

## COMPOSITION OF REACTION MIXTURE

Copper sulphate	11.8 mmol/L
Potassium sodium tartrate	31.3 mmol/L
Potassium iodide	29.5 mmol/L
Sodium hydroxide	588 mmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 520–560 nm may be used, general laboratory equipment.

- XL MULTICAL 4 × 3, Cat. No. XSYS0034
- XL MULTICAL 10 × 3, Cat. No. XSYS0122
- ERBA NORM 4 × 5, Cat. No. BLT00080
- ERBA NORM 10 × 5, Cat. No. XSYS0123
- ERBA PATH 4 × 5, Cat. No. BLT00081
- ERBA PATH 10 × 5, Cat. No. XSYS0124

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagent is stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.  
Plasma: Li-heparin and K<sub>2</sub>-EDTA plasma.  
The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.  
Centrifuge samples containing precipitates before performing the assay.  
See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma:	
	6 days at 20–25 °C
	4 weeks at 4–8 °C
	1 year at -20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.  
2 point calibration (blank and calibrator); distilled water is recommended as blank.  
Calibration frequency: it is recommended to do a calibration  
• after reagent lot change  
• as required by internal quality control procedures

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, calibrator XL MULTICAL, R2 standard and controls ERBA NORM and ERBA PATH have been standardized against the SRM 927 reference material.

## ASSAY PROCEDURE

Wavelength: 520–560 nm  
Cuvette: 1 cm

	Reagent blank	Calibrator (Standard)	Sample
Reagent 1	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.020 mL
Calibrator	–	0.020 mL	–
Distilled water	0.020 mL	–	–

Mix and incubate for 10 minutes (in case of automatic procedure incubate for 5 minutes) incubation in the dark. Absorbance of the sample  $A_{sam}$  and the calibrator (standard)  $A_{cal}$  against reagent blank is read in interval 30 minutes.

## CALCULATION

$$\text{Total protein} = \frac{A_{sam}}{A_{cal}} \times C_{cal} \quad C_{cal} = \text{calibrator (standard) concentration}$$

## ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point	Linearity Low (mg/dL)	0.26
Wavelength 1 (nm)	546	Linearity High (mg/dL)	14.0
Sample Volume (μL)	10/20	Concentration of Standard	See bottle label
Reagent Volume (μL)	500/1000	Blank with	Reagent
Incubation time (min.)	10	Absorbance limit (max.)	0.4
Reaction temperature (°C)	37	Units	g/dL
Reaction direction	Increasing		
Normal Low (g/dL)	6.4		
Normal High (g/dL)	8.3		

## UNIT CONVERSION

$$\text{g/dL} \times 10 = \text{g/L}$$

## EXPECTED VALUES\*

### In Serum:

Adults	6.4–8.3 g/dL
Premature	3.6–6.0 g/dL
Newborn	4.6–7.0 g/dL
1 week	4.4–7.6 g/dL
7–12 months	5.1–7.3 g/dL
1–2 years	5.6–7.5 g/dL
>2 years	6.0–8.0 g/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 0.26 g/dL  
Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV < 20 % (n = 30).

**Linearity:** 14.0 g/dL  
Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

**Precision:**  
Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	6.06	0.056	0.93
Sample 2	9.39	0.041	0.43

Intermediate precision	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	5.99	0.114	1.91
Sample 2	9.35	0.158	1.69

## Accuracy

Two different validated control materials were used. Determined bias is -2.4 % at the target value 7.21 g/dL and 1.7 % at the target value 8.43 g/dL.

## Comparison

A comparison between XL-640 automatic system TOTAL PROTEIN (y) and a commercially available test (x) using 150 samples gave following results:

Linear regression:  
 $y = 0.955x + 0.183 \text{ g/dL} \quad r = 0.987$   
 Passing-Bablok:  
 $y = 0.964x + 0.135 \text{ g/dL} \quad r = 0.978$

## Interferences

Criterion: Recovery within ± 10 % of initial value of total protein concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 3 g/L, bilirubin up to 28 mg/dL, triglycerides up to 320 mg/dL.  
Drugs: No interference was found at therapeutic concentrations using common drug panels<sup>10</sup>.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 546 nm against the distilled water is 0.4.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of total protein. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazard identification in accordance with Regulation (EC) No 1272/2008

R1  
 UFI: ARHK-VJ1C-KE7U-9PC4

## Danger

Contains: sodium hydroxide

## Hazard statement:

- H314 Causes severe skin burns and eye damage.
- H412 Harmful to aquatic life with long lasting effects.

## Precautionary statement:

- P260 Do not breathe vapours.
- P280 Wear protective gloves / protective clothing / eye protection / face protection.
- P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
- P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

## R2

Reagent is not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Weichselbaum TE., An accurate and rapid method for the determination of proteins in small amounts of blood serum and plasma. Am J Clin Pathol,10, 40-49 (1946)
2. Cornall, A. G., Bardawill, C. J., David, M. M., Determination of serum proteins by means of the biuret reaction, J. Biol. Chem. 177, 751, 1949.
3. Dumas, B. T., Bayse, D. D., Carter, R. J., Peters, T., Jr., and Schaffer, R., A candidate Reference Method for determination of total protein in serum. I. Development and validation. Clin. Chem. 27, 1642-1650 (1981).
4. Chromy, V., and Fischer, J., Photometric determination of total protein in lipemic sera. Clin. Chem. 23, 754-756 (1977).
5. Chromý, V., Fischer, J., Vozníček, Assay of total serum protein by the biuret method, reaction rate and lipemia, J. Z. Med. Labor.-Diagn. 21, 333-337, 1980.
6. Brobeck JR, ed. Physiological Basis of Medical Practice, 9th ed. Baltimore, MD: Wilkins and Wilkins 1973;4-7.
7. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
8. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11): 783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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