CREATININE ENZYMATIC

Cat. No.	Pack Name	Packaging (Content)
BLT00065	CREA ENZ 204	R1: 3 × 50 mL, R2: 3 × 18 mL, instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of creatinine in human serum, plasma and urine on various automatic systems. In combination with other parameters it is intended for screening, monitoring of renal function, diagnosis of renal diseases. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Creatinine is a waste product excreted by the kidneys mainly by glomerular filtration. The concentration of creatinine in plasma of a healthy individual is fairly constant, independent from water intake, exercise and rate of urine production. Therefore, increased plasma creatinine values always indicate decreased excretion, i.e. impaired kidney function. Creatinine clearance is a good indicator for the glomerular filtration rate (GFR) which allows better detection of kidney diseases and monitoring of renal function. For this purpose, creatinine is measured simultaneously in serum and urine collected over a defined time period. The serum creatinine levels do not start to rise until renal function has decreased by at least 50 %.

PRINCIPLE

In the first reaction, creatinase and sarcosine oxidase are used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide (H₂O₂) that is eliminated by catalase. Creatinase and 4-aminoantipyrine (4-AAP) are added, and only the creatine generated from creatinine by creatinase is hydrolysed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly formed hydrogen peroxide is measured in a coupled reaction catalysed by peroxidase, with N-ethyl-N-sulphopropyl-m-toluidine (TOPS) as a chromogen^{1,2}.

The absorbance of the produced quinone dye at 530–560 nm is proportional to the creatinine concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
MOPS pH (7.5)	25 mmol/L	MOPS (pH 7.5)	90 mmol/L
TOPS	0.5 mmol/L	Creatinase	30 kU/L
Creatinase	10 kU/L	Peroxidase	10 kU/L
Sarcosine oxidase	5 kU/L	Sodium azide	0.5 g/L
Catalase	3 kU/L		
EDTA	1 mmol/L		

COMPOSITION OF REACTION MIXTURE

MOPS pH (7.5)	41 mmol/L
TOPS	0.37 mmol/L
Creatinase	7.4 kU/L
Sarcosine oxidase	3.7 kU/L
Catalase	2.2 kU/L
Creatinase	7.4 kU/L
Peroxidase	2.5 kU/L
EDTA	0.7 mmol/L
Sodium azide	0.1 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 530–560 nm may be used, general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Reagents are ready to use. After opening, reagent is stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

Urine: Collect urine without using additives. If urine must be collected with a preservative for other analytes, only hydrochloric acid (14 to 47 mmol/L urine, e.g. 5 mL 10 % HCl or 5 mL 30 % HCl per liter urine) or boric acid (81 mmol/L, e.g. 5 g per liter urine) may be used. Dilute urine samples using redistilled water in 1 + 19 ratio and multiply results by 20.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma³:	7 days at	15–25 °C
	7 days at	2–8 °C
Stability in urine³:	3 months	-20 °C
	2 days at	15–25 °C
	6 days at	2–8 °C
	6 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS.

ASSAY PROCEDURE

Wavelength: 546 (530–560) nm

Cuvette: 1 cm

	Reagent blank	Calibrator	Sample
Reagent 1	0.90 mL	0.90 mL	0.90 mL
Sample	–	–	0.02 mL
Calibrator	–	0.02 mL	–
Distilled water	0.02 mL	–	–

Mix and incubate 3–5 min. at 37 °C. Measure absorbance A1 of the sample (A_{sam}), calibrator (A_{cal}) and blank (A_{bl}). Then add:

Reagent 2	0.30 mL	0.30 mL	0.30 mL
-----------	---------	---------	---------

Mix and incubate 5–10 min. at 37 °C. Measure absorbance A2 of the sample (A_{sam}), calibrator (A_{cal}) and blank (A_{bl}).

CALCULATION

$$\text{Creatinine} = \frac{\Delta A_{\text{sam}} - \Delta A_{\text{bl}}}{\Delta A_{\text{cal}} - \Delta A_{\text{bl}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

UNIT CONVERSION

mg/dL × 88.4 = μmol/L

EXPECTED VALUES

Serum⁶:

0–1 y	0.04–0.33 mg/dL
2–5 y	0.04–0.45 mg/dL
6–9 y	0.20–0.52 mg/dL
10 y	0.22–0.59 mg/dL
Adult:	
Male	0.62–1.10 mg/dL
Female	0.45–0.75 mg/dL

Urine⁷:

Infant	8–20 mg/kg/day
Child	8–22 mg/kg/day
Adolescent	8–30 mg/kg/day
Adult:	
Male	14–26 mg/kg/day
Female	11–20 mg/kg/day

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

Serum/plasma	0.024 mg/dL
Urine	0.49 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV < 20 % (n = 30).

Linearity:

Serum/plasma	85.8 mg/dL
Urine	1716 mg/dL

Linearity is the highest measured activity with recovery within ± 10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.08	0.009	0.86	Sample 1	106.5	0.90	0.85
Sample 2	3.36	0.018	0.53	Sample 2	222.1	0.94	0.43

Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.13	0.040	3.58	Sample 1	70.2	1.25	1.78
Sample 2	3.59	0.092	2.57	Sample 2	158.6	1.76	1.11

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is -3.2 % at the target value 1.33 mg/dL, -1.5 % at the target value 3.94 mg/dL for serum, 3.0 % at the target value 67.8 mg/dL and 2.0 % at the target value 143.7 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system CREATININE ENZYMATIC (y) and a commercially available test (x) using 147 samples (serum) gave following results:

Linear regression: y = 1.052x - 0.033 mg/dL r = 0.997

Passing-Bablok⁸: y = 1.043x - 0.010 mg/dL r = 0.996

Interferences

Criterion: Recovery within ± 10 % of initial value of creatinine concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Ca-dobesilate, Levodopa, Ethamsylate (Dicynone), N-acetylcysteine, Metamizole and Acetaminofen (Paracetamol) including its metabolite N-acetyl-p-benzoquinone imine may cause false-negative results in serum or plasma^{9,10,11,12,13}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 546 nm against the distilled water is 0.2.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of creatinine. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
- Börner U. Szaz G. et. Al. A specific fully enzymatic method for creatinine reference values in serum. J. Clin. Chem. Clin. Biochem 17: 679-882. 1979.
- Badiou S. Dupuy AM. Descamps B. Cristolead. JP. Comparison between the enzymatic vitros assay for creatinine determination and three other methods adapted on the Olympus analyzer. Journal of Clinical Laboratory Analysis 17: 235-240. 2003.
- Guder W. Fonseca-Wollheim W. Ehret W. et al. Die Qualität Diagnostischer Proben. 6. Aufl. Heidelberg: BD Diagnostics. 2009.
- Newman DJ. Price CP. Renal function and nitrogen metabolites. In: Burtis CA. Ashwood ER. editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 1204-1270. 1999.
- Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 366-374. 1998.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis CA. Ashwood. ER. Bruns. DE; 5th edition. WB Saunders Comp. 2012.
- Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. Rifai N. Horvath AR. Wittwer CT; 8th edition. Elsevier. 2019.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
- McCudden C, Clark EG, Akbari A, Kong J, Kanji S, Hiremath S, N-Acetylcysteine Interference with Creatinine Measurement: An In Vitro Analysis, Kidney Int Rep 6 (7), 1973 -1976, 2021.
- Dastyh M, Wiewiorka O, Benovska M. Ethamsylate (Dicynone) Interference in Determination of Serum Creatinine, Uric Acid, Triglycerides, and Cholesterol in Assays Involving the Trinder Reaction; In Vivo and In Vitro. Clin Lab 60, 1373-1376, 2014.
- Guinna A, Singha C, Piggotta Z, Palatnick W, Significance of falsely low creatinine values in diagnosing massive acetaminophen ingestion, Clinical Toxicology 60 (52), 65, 2022.
- Steinbach D, Racek J, Rajdl D, Interference of natural metabolites and drugs in enzymatic determination of creatinine and uric acid, Klin. Biochem. Metab., 29 (50), 132-138, 2021.

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LITERATÚRA

1. Börner U. Szaz G. et. Al. A specific fully enzymatic method for creatinine reference values in serum. J. Clin. Chem. Clin. Biochem 17: 679-882. 1979.
2. Badiou S. Dupuy AM. Descamps B. Cristolead. JP. Comparison between the enzymatic vitros assay for creatinine determination and three other methods adapted on the Olympus analyzer. Journal of Clinical Laboratory Analysis 17: 235-240. 2003.
3. Guder W. Fonseca-Wollheim W. Ehret W. et al. Die Qualität Diagnostischer Proben. 6. Aufl. Heidelberg: BD Diagnostics. 2009.
4. Newman DJ. Price CP. Renal function and nitrogen metabolites. In: Burtis CA. Ashwood ER. editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 1204-1270. 1999.
5. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 366-374. 1998.
6. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis CA. Ashwood. ER. Bruns. DE; 5th edition. WB Saunders Comp. 2012.
7. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. Rifai N. Horvath AR. Wittwer CT; 8th edition. Elsevier. 2019.
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
10. McCudden C, Clark EG, Akbari A, Kong J, Kanji S, Hiremath S, N-Acetylcysteine Interference with Creatinine Measurement: An In Vitro Analysis, Kidney Int Rep 6 (7), 1973 -1976, 2021.
11. Dastych M, Wiewiorka O, Benovska M. Ethamsylate (Dicynone) Interference in Determination of Serum Creatinine, Uric Acid, Triglycerides, and Cholesterol in Assays Involving the Trinder Reaction; In Vivo and In Vitro. Clin Lab 60, 1373-1376, 2014.
12. Guinna A, Singha C, Piggotta Z, Palatnick W, Significance of falsely low creatinine values in diagnosing massive acetaminophen ingestion, Clinical Toxicology 60 (52), 65, 2022.
13. Steinbach D, Racek J, Rajdl D, Interference of natural metabolites and drugs in enzymatic determination of creatinine and uric acid, Klin. Biochem. Metab., 29 (50), 132-138, 2021.

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