CHOLINESTERASE

Cat. No.	Pack Name	Packaging (Content)
BLT00066	CHE 120	R1: 2 × 50 mL, R2: 1 × 20 mL, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of cholinesterase (CHE) in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Two related enzymes have the ability to hydrolyze acetylcholine. One is acetylcholinesterase, which is called true CHE, or CHE I. True CHE is found in erythrocytes, in the lungs and spleen, in nerve endings, and in the gray matter of the brain. The other CHE is acylcholine acylhydrolase, it is usually called pseudocholinesterase, benzoyl cholinesterase, or CHE II. Although it is found in the liver, pancreas, heart, white matter of the brain, and serum, its biological role is unknown.

The serum enzyme is the one whose assay is clinically useful. Serum cholinesterase measurements are useful as an indicator of insecticide poisoning. Patients exposed to organophosphorus compounds exhibit decreased serum cholinesterase levels. CHE measurements are also used in the detection of patients with inherited abnormal variants of the enzyme. Some of these variants lead to prolonged apnea in patients receiving the anesthetic succinylcholine¹. In addition, the measurement of CHE activity is used in the determination of the synthetic capacity of the liver. In the absence of genetic causes or known inhibitors, any decrease in activity in serum reflects impaired synthesis of the enzyme by the liver. A 30–50 % decrease in level is observed in acute hepatitis and in chronic hepatitis of long duration. Decreases of 50–70 % occur in advanced cirrhosis and carcinoma with metastases to the liver. Essentially normal levels are seen in chronic hepatitis, mild cirrhosis, and obstructive jaundice².

PRINCIPLE

The method uses butyrylthiocholine as the specific substrate for CHE^{3,4}. CHE catalyses the hydrolysis of butyrylthiocholine substrate forming butyrate and thiocholine. Thiocholine reduces hexacyanoferrate (III) to hexacyanoferrate (II):

The decrease of absorbance at 405 nm is proportional to the activity of CHE in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	R2
Pyrophosphate buffer, pH 7.6 (37 °C)	Butyrylthiocholine
Hexacyanoferrate (III+)	
92 mmol/L	91 mmol/L
2.5 mmol/L	

COMPOSITION OF REACTION MIXTURE

Pyrophosphate buffer	75 mmol/L
Hexacyanoferrate (III+)	2.0 mmol/L
Butyrylthiocholine	15 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 405 nm may be used, general laboratory equipment.

- XL MULTICAL 4×3, Cat. No. XSYS0034
- XL MULTICAL 10×3, Cat. No. XSYS0122
- ERBA NORM 4×5, Cat. No. BLT00080
- ERBA NORM 10×5, Cat. No. XSYS0123
- ERBA PATH 4×5, Cat. No. BLT00081
- ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After the first opening the vials, reagents are stable for 60 days at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

- Serum.
- Plasma: Li-heparin plasma.
- Fresh serum, plasma not haemolyzed and promptly separated from the red blood cells⁵. Do not use sodium fluoride as an anticoagulant because it inhibits CHE.
- The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.
- See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ^{2,5,6} :	6 hours at	15–25 °C
	7 days at	2–8 °C
	1 year at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended. 2 point calibration (blank and calibrator); distilled water is recommended as blank. Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to DGKC method⁷.

ASSAY PROCEDURE

Wavelength: 405 (400–420) nm
Cuvette: 1 cm

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.020 mL
Calibrator	–	0.020 mL	–
Distilled water	0.020 mL	–	–

Mix and after 5 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix and incubate for 90 seconds at 37 °C and read initial absorbance of sample, calibrator and reagent blank. Repeat the reading after exactly 30, 60 and 90 seconds. Calculate the mean 30 seconds absorbance change (ΔA_{30s}).

CALCULATION

- CHE (U/L) = $\frac{\Delta A_{\text{sam}}/30s - \Delta A_{\text{blank}}/30s}{\Delta A_{\text{cal}}/30s - \Delta A_{\text{blank}}/30s} \times C_{\text{cal}}$ C_{cal} = calibrator concentration
- Using factor:
CHE(U/L) = $f \times (\Delta A_{\text{sam}}/30s - \Delta A_{\text{blank}}/30s)$ $f = 131600$ (at 405 nm)

UNIT CONVERSION

U/L × 0.0167 = $\mu\text{kat/L}$

EXPECTED VALUES⁸

At 37 °C

Females: 4000 – 12600 U/L
Males: 5100 – 11700 U/L

In infants up to 6 months of age, cholinesterase activity is 40 % to 50 % higher than in adults. In young women, the enzyme activity is approximately 64 % to 74 % higher than in adult males. The activity decreases in pregnancy. It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 168 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 16200 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	6585	35.4	0.54
Sample 2	9736	49.0	0.50

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	6499	134.1	2.06
Sample 2	9718	150.1	1.54

Accuracy

Two different validated control materials were used. Determined bias is 0.5 % at the target value 8400 U/L and 2.2 % at the target value 10200 U/L.

Comparison

A comparison between XL-640 automatic system CHE (y) and a commercially available test (x) using 46 samples gave following results:

Linear regression: $y = 0.914x - 331$ U/L $r = 0.990$

Passing-Bablok⁹: $y = 0.915x - 367$ U/L $r = 0.985$

Interferences

Criterion: Recovery within ±10 % of initial value of CHE activity in the sample without interfering substance. Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 28 mg/dL, triglycerides up to 850 mg/dL. Drugs: No interference was found at therapeutic concentrations using common drug panels⁹.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 405 nm against the distilled water is 0.55. Reagents must be limp; do not use if turbid.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of CHE. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

The reagent is not classified as dangerous.

R2

UFI: XFJU-7W7J-SJ5M-36KS

Warning

Contains: maleic acid

Hazard statement:
H317 May cause an allergic skin reaction.

Precautionary statement:
P280 Wear protective gloves/protective clothing/eye protection.
P302 + P352 IF ON SKIN: Wash with plenty of water and soap.
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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


1. Kaplan LA, Pesce AJ. Clinical Chemistry, Theory, Analysis and Correlation. Ladig D, Kasper R (ed), St Louis, CV Mosby Co 1984;1108-1109.
2. Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company 1995;132-133..
3. Deutsche Gesellschaft für Klinische Chemie: Proposal of standard methods for determination of enzyme catalytic concentration in serum and plasma at 37 °C II. Cholinesterase (acetylcholine acylhydrolase, E.C.3.1.1.8.). Eur.J.Clin.Chem. Biochem 30, 163 (1992).
4. Schmidt E, Gerhardt W, Henkel E, et al. Proposal of Standard Methods for the Determination of Enzyme Catalytic Concentrations in Serum and Plasma at 37 °C. Eur J Clin Chem Clin Biochem 1992;30:163-170.
5. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
6. Huizenga JR, van der Belt K, Gips CH. The Effect of Storage at Different Temperatures on Cholinesterase Activity in Human Serum. J Clin Chem Clin Biochem 1985;24:283-385.
7. Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996)
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.

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