CHOLESTEROL

Cat. No.	Pack Name	Packaging (Content)
XSYS0009	CHOL 440	R1: 10 × 44 mL, RFID tag instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of cholesterol in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of coronary artery disease, liver function. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Cholesterol is essential structural component of cell membranes and precursor of bile acids and all steroids hormones. This is why cholesterol has enormous significance for organism normal functioning. But there is also well established association between blood cholesterol concentration and coronary heart disease. Measurement of cholesterol serum level is valuable in prevention and monitoring cardiovascular disease. This determination is useful also for evaluation of intestine absorption, liver and gallbladder function.

PRINCIPLE

Enzymatic, colorimetric method with cholesterol esterase (CHE) and cholesterol oxidase (CHO). Cholesterol is determined after enzymatic hydrolysis and oxidation. Cholesterol esters are cleaved by the action of CHE to yield free cholesterol and fatty acids. CHO then catalyzes the oxidation of cholesterol to cholest-4-en-3-one and hydrogen peroxide (H₂O₂). In the presence of peroxidase (POD), the hydrogen peroxide formed effects the oxidative coupling of phenol and 4-amino-phenazone to form a red quinoneimine dye (Trinder's reaction)^{1,2,3,4,5,6}.

Absorbance of Quinoneimine dye measured at 505 nm is proportional to the concentration of cholesterol in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
Good's Buffer (pH 6.7)	50 mmol/L
Phenol	5.0 mmol/L
4-aminoantipyrine	0.3 mmol/L
Cholesterol esterase	≥200 U/L
Cholesterol oxidase	≥50 U/L
Peroxidase	≥3 KU/L
Sodium azide	0.5 g/L

COMPOSITION OF REACTION MIXTURE

Good's Buffer (pH 6.7)	49.5 mmol/L
Phenol	4.95 mmol/L
4-aminoantipyrine	0.3 mmol/L
Cholesterol esterase	≥198 U/L
Cholesterol oxidase	≥50 U/L
Peroxidase	≥3 KU/L
Sodium azide	0.5 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4 × 3, Cat. No. XSYS0034
 XL MULTICAL 10 × 3, Cat. No. XSYS0122
 ERBA NORM 4 × 5, Cat. No. BLT00080
 ERBA NORM 10 × 5, Cat. No. XSYS0123
 ERBA PATH 4 × 5, Cat. No. BLT00081
 ERBA PATH 10 × 5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
 On board stability: min. 60 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.
 For specimen collection and preparation only use suitable tubes or collection containers.
 Only the specimens listed below were tested and found acceptable.

Serum.
 Plasma: Li-heparin and K₂-EDTA plasma.
 Do not use citrate, oxalate or fluoride.
 Fasting and nonfasting samples can be used⁹.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.
 See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ^{9,10} :	7 days at	15–25 °C
	7 days at	4–8 °C
	3 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank
 Calibration frequency: 30 days
 Calibration is needed:
 • after reagent lot change
 • as required by internal quality control procedures
 • calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.
 The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS¹³.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	505/700 nm
Reading time	10 min. after adding of R1
Reaction direction	Increase
Unit	mg/dL (mmol/L)
Reagent volumes	
R1	200 µL
Sample volumes	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

UNIT CONVERSION

mg/dL × 0.026 = mmol/L

EXPECTED VALUES¹⁴

In serum:	Male	Female
0–4 y	114–203	112–200 mg/dL
5–9 y	125–189	131–197 mg/dL
10–14 y	124–204	125–205 mg/dL
15–19 y	118–191	119–208 mg/dL
20–24 y	118–212	121–237 mg/dL
25–29 y	130–234	130–231 mg/dL
30–34 y	142–258	133–227 mg/dL
35–39 y	147–267	139–249 mg/dL
40–44 y	150–260	146–259 mg/dL
45–49 y	163–275	148–268 mg/dL
50–54 y	156–274	163–281 mg/dL
55–59 y	161–280	167–294 mg/dL
60–64 y	163–287	172–300 mg/dL
65–69 y	166–288	167–291 mg/dL
>69 y	144–265	173–280 mg/dL

Coronary heart disease risk:

Adult	
Desirable	<200 mg/dL
Borderline high	200–239 mg/dL
High	>239 mg/dL
Child	
Desirable	<170 mg/dL
Borderline high	170–199 mg/dL
High	>199 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification:

1.4 mg/dL
 Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity:

750 mg/dL
 Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	142	0.81	0.57	Sample 1	143	2.61	1.83
Sample 2	258	1.09	0.42	Sample 2	259	5.56	2.15

Accuracy

Two different validated control materials were used. Determined bias is and -2.7 % at the target value 149 mg/dL and -3.7 % at the target value 197 mg/dL.

Comparison

A comparison between XL-640 automatic system CHOLESTEROL (y) and a commercially available test (x) using 120 samples gave following results:

Linear regression:
 $y = 0.964x + 6.229 \text{ mg/dL}$ $r = 0.997$
 Passing-Bablok¹⁵:
 $y = 0.967x + 5.517 \text{ mg/dL}$ $r = 0.996$

Interferences

Criterion: Recovery within ±10 % of initial value of cholesterol concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 4.5 g/L, bilirubin up to 6 mg/dL, triglycerides up to 850 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels except N-acetylcysteine, Metamizole and Acetaminofen (including metabolite N-acetyl-p-benzoquinone imine)^{16,17}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of cholesterol. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1
 Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
- Artiss JD, Zak B, Measurement of cholesterol concentration. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 99–114, 1997.
- Deeg R, Ziegenhorn J, Kinetic enzymatic method for automated determination of total cholesterol in serum. Clin Chem 29: 1798-1802, 1983.
- Flegg HM, An., Preparation and properties of a cholesterol oxidase from Nocardia sp. and its application to the enzymatic assay of total cholesterol in serum. Clin Chem 19(12): 1350-1356, 1973.
- Allain CC, Poon LS, Chan CS, et al. Enzymatic determination of total serum cholesterol. Clin Chem 20(4): 470-475, 1974.
- Roeschlau P, Bernt E, Gruber WA. Enzymatic determination of total cholesterol in serum. Z Klin Chem Klin Biochem 12: 226, 1974.
- Barham D, Trinder P, An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
- Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press p.176.
- Pisani T, GebSKI CP, Leary ET, et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med Dec 119(12): 1127-1135, 1995.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company: 130-131, 1995.
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
- Rifai N, Bachorik PS, Albers JJ. Lipids, lipoproteins and apolipoproteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia: W.B Saunders Company: 809-61, 1999.
- Recommendation of the Second Joint Task Force of European and other Societies on Coronary Prevention. Prevention of coronary heart disease in clinical practice. Eur Heart J 19: 1434-503, 1998.
- Siekmann L, Hüskes KP, Breuer H. Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279: 145-146, 1976.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
- Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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LITERATÚRA

1. Artiss JD, Zak B, Measurement of cholesterol concentration. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 99–114, 1997.
2. Deeg R, Ziegenhorn J, Kinetic enzymatic method for automated determination of total cholesterol in serum. Clin Chem 29: 1798-1802, 1983.
3. Flegg HM, An., Preparation and properties of a cholesterol oxidase from Nocardia sp. and its application to the enzymatic assay of total cholesterol in serum. Clin Chem 19(12): 1350-1356, 1973.
4. Allain CC, Poon LS, Chan CS, et al. Enzymatic determination of total serum cholesterol. Clin Chem 20(4): 470-475, 1974.
5. Roeschlau P, Bernt E, Gruber WA. Enzymatic determination of total cholesterol in serum. Z Klin Chem Klin Biochem 12: 226, 1974.
6. Barham D, Trinder P, An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
7. Nader R, Dufour DR, Cooper GR. Preanalytical Variation in Lipid, Lipoprotein, and Apolipoprotein Testing. In: Rifai N, Warnick GR, and Dominiczak MH, editors. Handbook of Lipoprotein Testing. 2nd ed. Washington: AACC press p.176.
8. Pisani T, Gebiski CP, Leary ET, et al. Accurate Direct Determination of Low-density Lipoprotein Cholesterol Using an Immunoseparation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med Dec 119(12): 1127-1135, 1995.
9. Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company: 130-131, 1995.
10. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
11. Rifai N, Bachorik PS, Albers JJ. Lipids, lipoproteins and apolipoproteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 809-61, 1999.
12. Recommendation of the Second Joint Task Force of European and other Societies on Coronary Prevention. Prevention of coronary heart disease in clinical practice. Eur Heart J 19: 1434-503, 1998.
13. Siekmann L, Hüskes KP, Breuer H. Determination of cholesterol in serum using mass fragmentography - a reference method in clinical chemistry. Z Anal Chem 279: 145-146, 1976.
14. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
15. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
16. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
17. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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