PHOSPHORUS

Cat. No.	Pack Name	Packaging (Content)
XSYS0015	PHOS 120	R1: 10 × 12 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of phosphorus in human serum, plasma and urine on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of renal diseases and liver diseases (bone metabolism). For professional use in clinical laboratory only.


CLINICAL SIGNIFICANCE

More than 80 % of the phosphorus contained in the body is localized in bone in the form of calcium phosphate as the apatite. The remainder is involved in intermediary carbohydrate metabolism and in physiologically important substances such as phospholipids, nucleic acids and ATP. Phosphorus occurs in blood in the form of inorganic phosphate and in organically bound phosphoric acid. The small amount of extracellular organic phosphorus is found almost exclusively in the form of phospholipids. The ratio of phosphate to calcium in the blood is approximately 6:10. An increase in the level of phosphorus causes a decrease in the calcium level.

The mechanism is influenced by interactions between parathormone and vitamin D. Hypoparathyroidism, vitamin D intoxication and renal failure with decreased glomerular phosphate filtration give rise to hyperphosphatemia. Hypophosphatemia occurs in rickets, hyperparathyroidism and Fanconi's syndrome.

PRINCIPLE

Molybdate UV. Inorganic phosphate forms an ammonium phosphomolybdate complex with ammonium molybdate in the presence of sulfuric acid^{1,2}.

The concentration of phosphomolybdate formed is directly proportional to the inorganic phosphate concentration and is measured as increasing of absorbance at 340 nm.

REAGENT DESCRIPTION AND COMPOSITION

R1

Ammonium molybdate	1.00 mmol/L
Sulphuric acid	336 mmol/L

COMPOSITION OF REACTION MIXTURE

Ammonium molybdate	0.98 mmol/L
Sulphuric acid	329 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.
 Urine: Collect in detergent-free containers. Acidify with hydrochloric acid after collection (pH <3)^{3,4}. Dilute urine 1 + 9 with dist. water before determination and multiply the result by 10. The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay. See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma⁵:	24 hours at	15–25 °C
	4 days at	4–8 °C
	1 year at	-20 °C

Stability in urine^{3,4}:	6 months at	2–8 °C (when acidified)
	24-hour urine:	Store cooled during collection

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended. 2 point calibration (blank and calibrator); distilled water is recommended as blank. Calibration frequency: 14 days
 Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the primary reference material.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	340/700 nm
Reading time	10 min. after adding of R1
Reaction direction	Increase
Unit	mg/dL (mmol/L)
Reagent volumes	
R1	200 µL
Sample volumes	4 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

UNIT CONVERSION

mg/dL × 0.32 = mmol/L

EXPECTED VALUES¹¹

Serum:

Adult: 2.5–4.5 mg/dL

Children: 4.0–7.0 mg/dL

Urine:

Adult: 0.4–1.3 g/24 h

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification:

Serum / plasma	0.062 mg/dL
Urine	0.67 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity:

Serum / plasma	20 mg/dL
Urine	200 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	4.37	0.027	0.61
Sample 2	6.63	0.035	0.53

Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	4.45	0.168	3.77
Sample 2	6.95	0.204	2.94

Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	26.3	0.27	1.01
Sample 2	53.6	0.38	0.72

Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	27.5	0.43	1.56
Sample 2	52.6	0.95	1.80

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is 5.4 % at the target value 3.22 mg/dL, 4.0 % at the target value 5.38 mg/dL for serum, -6.2 % at the target value 22.4 mg/dL and -3.2 % at the target value 44.4 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system PHOSPHORUS (y) and a commercially available test (x) using 116 samples (serum) gave following results:

Linear regression:

$$y = 0.997x + 0.320 \text{ mg/dL} \quad r = 0.988$$

Passing-Bablok¹²:

$$y = 0.990x + 0.301 \text{ mg/dL} \quad r = 0.960$$

Interferences

Criterion: Recovery within ±10 % of initial value of phosphorus concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 2 g/L, bilirubin up to 4 mg/dL, triglycerides up to 570 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels in serum or in urine¹³.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of phosphorus. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UF1: 9YDF-8EDQ-8H7H-KTFP

Danger

Hazard Statements:

- H314 Causes severe skin burns and eye damage.
- H412 Harmful to aquatic life with long lasting effects.

Precautionary Statements:

- P280 Wear protective gloves / protective clothing / eye protection.
- P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
- P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Henry R ed. Clinical Chemistry: Principles and Technics, 2nd ed. New York, NY: Harper & Row: 723, 1974.
2. Daly JA, Ertingshausen G, Direct Method for Determining Inorganic Phosphate in Serum with the „CentrifChem“, Clin Chem 18: 263, 1972.
3. Tietz NW, ed. Clinical Guide to Laboratory Tests, 4th ed. Philadelphia.WB Saunders Co;852-855, 2006.
4. NCCLS GP-16A2, Urineanalysis and Collection, Transportation and Preservation of Urine specimens, 2nd edition 2001.
5. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
6. Külpmann WR, Stummvoll HK, Lehmann P. Elektrolyte, Klinik und Labor. Heidelberg: Verlag Klinisches Labor 1993.
7. Tietz NW, ed. Fundamentals of Clinical Chemistry Philadelphia, PA: WB Saunders Company: 901, 1976.
8. Fiske CH, Subbarow Y. The colorimetric determination of phosphorus. J Biol Chem 66: 375-400, 1925.
9. Taussky HH, Schoor EA. A microcolorimetric method for the determination of inorganic phosphorus. J Biol Chem 202: 675, 1953.
10. Henry R ed. Clinical Chemistry: Principles and Technics, 2nd ed. New York, NY: Harper & Row: 723, 1974.
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
12. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
13. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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LITERATÚRA

1. Henry R ed. Clinical Chemistry: Principles and Technics, 2nd ed. New York, NY: Harper & Row: 723, 1974.
2. Daly JA, Ertingshausen G, Direct Method for Determining Inorganic Phosphate in Serum with the „CentrifChem“, Clin Chem 18: 263, 1972.
3. Tietz NW, ed. Clinical Guide to Laboratory Tests, 4th ed. Philadelphia.WB Saunders Co;852-855, 2006.
4. NCCLS GP-16A2, Urineanalysis and Collection, Transportation and Preservation of Urine specimens, 2nd edition 2001.
5. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
6. Külpmann WR, Stummvoll HK, Lehmann P. Elektrolyte, Klinik und Labor. Heidelberg: Verlag Klinisches Labor 1993.
7. Tietz NW, ed. Fundamentals of Clinical Chemistry Philadelphia, PA: WB Saunders Company: 901, 1976.
8. Fiske CH, Subbarow Y. The colorimetric determination of phosphorus. J Biol Chem 66: 375-400, 1925.
9. Tausky HH, Schoor EA. A microcolorimetric method for the determination of inorganic phosphorus. J Biol Chem 202: 675, 1953.
10. Henry R ed. Clinical Chemistry: Principles and Technics, 2nd ed. New York, NY: Harper & Row: 723, 1974.
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
12. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
13. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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REF

Katalógové číslo

LOT

Číslo šarže

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eFU:
www.erba.com**IVD**Diagnostický zdravotnícky prostriedok *in vitro*

Výrobca

Obmedzenie teploty

CONT

Obsah

QUALITY SYSTEM CERTIFIED
ISO 13485, IVDRErba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com

CC/IFU/063/26/A

Dátum revízie: DD. MM. YYYY