TRIGLYCERIDES

Cat. No.	Pack Name	Packaging (Content)
XSYS0041	TG 440	R1: 10 × 44 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of triglycerides in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of hyperlipidemias. For professional use in clinical laboratories only.


CLINICAL SIGNIFICANCE

Triglycerides are esters of the trihydric alcohol glycerol with 3 long-chain fatty acids. They are partly synthesized in the liver and partly ingested in food.

Measurement of triglycerides is important in the diagnosis and management of hyperlipidemias. These diseases can be genetic or secondary to other disorders including nephrosis, diabetes mellitus and endocrine disturbances. Elevation of triglycerides has been identified as a risk factor for atherosclerotic disease.

PRINCIPLE

Triglycerides are enzymatically hydrolyzed by lipoproteinlipase (LPL) to free acids and glycerol. The glycerol is phosphorylated by adenosine triphosphate (ATP) with glycerol kinase (GK) to produce glycerol-3-phosphate and adenosine diphosphate (ADP). Glycerol-3-phosphate is oxidized to dihydroxy-acetone phosphate (DAP) by glycerol-3-phosphate oxidase (GPO) producing hydrogen peroxide (H₂O₂). In a Trinder type colour reaction¹ catalyzed by peroxidase (POD), the H₂O₂ reacts with 4-aminoantipyrine (4AAP) and 4-chlorophenol to produce a red coloured quinoneimine dye^{2,3,4}.

Absorbance of coloured quinoneimine complex measured at 505 nm is proportional to the concentration of triglycerides in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		
Good's buffer (pH 7.2)	50	mmol/L
4-chlorophenol	4	mmol/L
4-aminoantipyrine	0.5	mmol/L
Mg ²⁺	15	mmol/L
ATP	2	mmol/L
Glycerolkinase	≥0.4	kU/L
Peroxidase	≥2.0	kU/L
Lipoproteinlipase	≥2.0	kU/L
Glycerol-3-phosphate oxidase	≥0.5	kU/L

COMPOSITION OF REACTION MIXTURE

Good's buffer (pH 7.2)	49.5	mmol/L
4-chlorophenol	3.96	mmol/L
4-aminoantipyrine	0.495	mmol/L
Mg ²⁺	14.9	mmol/L
ATP	1.98	mmol/L
Glycerolkinase	≥0.396	kU/L
Peroxidase	≥1.98	kU/L
Lipoproteinlipase	≥1.98	kU/L
Glycerol-3-phosphate oxidase	≥0.495	kU/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 60 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma ⁵ :		
	2 days at	20–25 °C
	7 days at	4–8 °C
	1 year at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	505/700 nm
Reading time	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (mmol/L)

Reagent volumes

R1 200 µL

Sample 2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1: sample does not change.

UNIT CONVERSION

mg/dL × 0.0113 = mmol/L

EXPECTED VALUES⁶

Recommended triglycerides levels for adults:

Normal: <150 mg/dL

High: 150–199 mg/dL

Hypertriglyceridemic: 200–499 mg/dL

Very high: >499 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 1.64 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV < 20% (n = 30).

Linearity: 1250 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	106.4	3.13	2.94	Sample 1	105.2	3.34	3.17
Sample 2	168.7	0.94	0.56	Sample 2	167.2	3.78	2.26

Accuracy

Two different validated control materials were used. Determined bias is and -0.8 % at the target value 91.3 mg/dL and -1.5 % at the target value 156.7 mg/dL.

Comparison

A comparison between XL-640 automatic system TRIGLYCERIDES (y) and a commercially available test (x) using 149 samples gave following results:

Linear regression:

$$y = 0.997x + 16.64 \text{ mg/dL} \quad r = 0.994$$

Passing-Bablok⁷:

$$y = 1.025x + 13.99 \text{ mg/dL} \quad r = 0.993$$

Interferences

Criterion: Recovery within ±10 % of initial value of triglycerides concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 5 g/L, bilirubin up to 8 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels except Ca-dobesilate, N-acetylcysteine, Metamizole and Acetaminofen (including metabolite N-acetyl-p-benzoquinone imine)^{8,9}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.

- High concentration of haemoglobin and bilirubin in sample can interfere with determination of triglycerides. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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


1. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. *Analyst* 1972; 97: 142-5.
2. Guder WG, Zawta B et al. *The Quality of Diagnostic Samples*. 1st ed. Darmstadt: GIT Verlag; 2001; 30-31.
3. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. 750-808.
4. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002; 48: 436-72.
5. Tietz NW, (Ed.), *Textbook of Clinical Chemistry*. Burtis CA and Ashwood ER, Fifth Edition, 2012.
6. *Tietz Textbook of Clinical Chemistry*. 5th ed. Edited by Burtis CA, Ashwood ER, Bruns DE. WB Saunders Company, Philadelphia 2006, 25: 837-907.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11): 783-790.
8. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 38, 376-385, 2001.
9. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, *Clin Biochem* 49, 100 -104, 2016.

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