MAGNESIUM

Cat. No.	Pack Name	Packaging (Content)
XSYS0040	MG 88	R1: 2 x 44 mL, RFID tag, instruction for use, chimney: 1 pc

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of magnesium in human serum, plasma and urine on automatic systems ERBA XL. Intended for screening, monitoring and diagnosis of impairment of neuromuscular functions, gastrointestinal and cardiac symptoms, treatment monitoring. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE


Magnesium is an essential nutrient which is involved in many biochemical functions. It has a structural role in nucleic acids and ribosomal particles, required as an activator for many enzymes and has a role in energy producing oxidative phosphorylation.

The normal body contains between 21 to 28 g magnesium, more than 50 % of which is complexed with calcium and phosphate in bones. Only approximately 1 % of the total magnesium is found in the extracellular fluid; hence, it tends to enter and leave cells under the same conditions as potassium. Approximately 35 % of plasma magnesium is protein-bound, mainly to albumin, and therefore changes in albumin concentration may affect magnesium.

Hypomagnesemia results in the impairment of neuromuscular function and may develop in severe prolonged diarrhea, malabsorption syndromes, hyperaldosteronism, and diuretic therapy. Hypermagnesemia is seen in renal glomerular failure and diabetic coma.

PRINCIPLE

Magnesium reacts with Xylidyl Blue to form a purple colored complex in alkaline solution. Interference by calcium is prevented by the use of glycoethyldiaminetetraacetic acid (GEDTA)^{1,2,3}.

The resulting absorbance at 505 nm is proportional to the magnesium concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
Ethanolamine (pH 11.0)	1.0 µmol/L
Xylidyl Blue	110 µmol/L
GEDTA	60 µmol/L

COMPOSITION OF REACTION MIXTURE

Ethanolamine (pH 11.0)	0.99 µmol/L
Xylidyl Blue	109 µmol/L
GEDTA	59 µmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4x3, Cat. No. XSYS0034
 XL MULTICAL 10x3, Cat. No. XSYS0122
 ERBA NORM 4x5, Cat. No. BLT00080
 ERBA NORM 10x5, Cat. No. XSYS0123
 ERBA PATH 4x5, Cat. No. BLT00081
 ERBA PATH 10x5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated and if chimneys are used.

Absorption of atmospheric CO₂ in open reagent bottle R1 leads to the reduced stability of reagents. To minimize the rate of absorption of CO₂ the use of "chimney" provided with the kit is recommended. Use of chimney has shown improvement of calibration stability, provided the storage and use of reagents is as per instructions. Chimney provided with the kit is to be placed directly into the reagent bottle immediately after opening the cap. Chimney can be re-used again for the next bottle of the same kit. Before using the chimney in a new bottle, it is recommended to wash chimney well with distilled water and dry it thoroughly.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin. Chelating anticoagulants such as EDTA, fluoride and oxalate must be avoided. Urine: Acidify urine with some drops of conc. HCl to pH 3–4 to prevent precipitation of magnesium ammonium phosphate. Collect urine samples in metal-free container². Dilute urine samples using redistilled water in 1 + 4 ration and multiply results by 5.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ^a :	
	7 days at 15–25 °C
	7 days at 2–8 °C
	1 year at -20 °C

Stability in urine ^a :	
	3 days at 15–25 °C
	3 days at 2–8 °C
	1 year at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 12 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the SRM 956 reference material.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	505 nm
Reading time	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (mmol/L)
Reagent volumes	
R1	200 µL
Sample	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

UNIT CONVERSION

mg/dL × 0.4114 = mmol/L

EXPECTED VALUES^a

Serum or plasma:

Newborn, 2–4 days	1.5–2.2 mg/dL
5 months–6 years	1.7–2.3 mg/dL
6–12 years	1.7–2.1 mg/dL
>12 years	1.6–2.6 mg/dL

Urine:

12–291 mg / 24 h

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification:

Serum / plasma	0.11 mg/dL
Urine	0.45 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity:

Serum / plasma	6.80 mg/dL
Urine	34.0 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.96	0.042	2.14	Sample 1	2.00	0.068	3.39
Sample 2	3.61	0.033	0.93	Sample 2	3.61	0.065	1.80

Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	2.92	0.050	1.71	Sample 1	2.05	0.045	2.17
Sample 2	5.21	0.040	0.77	Sample 2	4.21	0.141	3.36

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is 0.5 % at the target value 1.90 mg/dL, -5.1 % at the target value 5.16 mg/dL for serum, 1.7 % at the target value 6.32 mg/dL and -1.3 % at the target value 12.0 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system MAGNESIUM (y) and a commercially available test (x) using 150 samples (serum) gave following results:

Linear regression:

$$y = 0.917x + 0.222 \text{ mg/dL} \quad r = 0.973$$

Passing-Bablok¹⁰:

$$y = 0.929x + 0.198 \text{ mg/dL} \quad r = 0.965$$

Interferences

Criterion: Recovery within ±10 % of initial value of magnesium concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 9 g/L, bilirubin up to 24 mg/dL, triglycerides up to 850 mg/dL.

Drugs:

Serum: No interference was found at therapeutic concentrations using common drug panels¹¹.

Urine: No interference was found at therapeutic concentrations using common drug panels except Acetylcysteine, Gentamicin and Methyldopa¹¹.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.

Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of magnesium. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: MXWU-1WQR-TJ5W-RXPE

Danger

Contains: 2-aminoethanol

Hazard Statements:

H315 Causes skin irritation.

H318 Causes serious eye damage.

Precautionary Statements:

P280 Wear protective gloves / protective clothing / eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water and soap.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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7. Sitzmann FC. Normalwerte. München: Hans Marseille Verlag GmbH, 166, 1986.