# BICARBONATE

Cat. No.	Pack Name	Packaging (Content)
XSYS0100	CO2	R1: 4 × 34 mL, R2 standard: 1 × 3 mL, RFID tag, instruction for use

## INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of bicarbonate in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of disorders associated with disturbances of the metabolic and respiratory systems. For professional use in clinical laboratories only.

## CLINICAL SIGNIFICANCE

Bicarbonate is the second largest fraction of the anions in plasma. Included in this fraction are the bicarbonate ( $\text{HCO}_3^-$ ) and carbonate ( $\text{CO}_3^{2-}$ ) ions, as well as the carbamino compounds. At the physiological pH of blood, the concentration of carbonate is 1/1000 that of bicarbonate. The carbamino compounds are also present in such low quantities that they are generally not mentioned specifically. The bicarbonate content of serum or plasma is a significant indicator of electrolyte dispersion and anion deficit. Together with pH determination, bicarbonate measurements are used in the diagnosis and treatment of numerous potentially serious disorders associated with acid-base imbalance in the respiratory and metabolic systems.

## PRINCIPLE

Bicarbonate reacts with phosphoenolpyruvate (PEP), in the presence of phosphoenolpyruvate carboxylase (PEP-C), to form oxaloacetate and phosphate. The oxaloacetate is then converted to malate by the action of malate dehydrogenase (MDH) and reduced nicotinamide adenine dinucleotide analog (NADH-analog)<sup>1,2</sup>.

The decrease in absorbance at 405 or 415 nm resulting from the oxidation of NADH-analog is proportional to the amount of  $\text{CO}_2$  in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

R1	R2 standard	
Buffer, pH 7.5	Bicarbonate	30 mmol/L
PEP		12.5 mmol/L
PEC-C		>400 U/L
MDH		>4100 U/L
NADH-analog		0.6 mmol/L

## COMPOSITION OF REACTION MIXTURE

Buffer, pH 7.5	
PEP	12.4 mmol/L
PEC-C	>396 U/L
MDH	>4060 U/L
NADH-analog	0.594 mmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Erba XL analysers: XL-200, Cat. No. INS00002  
XL-640, Cat. No. INS00008  
XL-1000, Cat. No. INS00010

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin.

The preferred specimen is from venous blood collected anaerobically in the usual manner for bicarbonate analysis. Bicarbonate content in uncapped tubes decreases approximately 4 mmol/L after one hour<sup>6</sup>. It has been reported that alkalized serum stored in open cups is stable for up to 4 hours<sup>6</sup>.

Separate from erythrocytes and store tightly stoppered.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

<b>Stability in serum / plasma<sup>7</sup>:</b>	1 day at	20–25 °C
	7 days at	4–8 °C
	2 weeks at	-20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with the R2 standard included in the kit is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 14 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

## QUALITY CONTROL

For quality control all control solutions with  $\text{CO}_2$  values determined by this method can be used. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method has been standardized against a primary standard on basis of sodium carbonate.

## ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see [www.erbalachema.com](http://www.erbalachema.com).

## Assay parameters for ERBA XL automatic systems

Assay type	2-Point
Curve type	Linear
Wavelength (prim. / sec.)	405 (415) / 505 nm
Reading time 1	1 min after adding of R1
Reading time 2	9 min after adding of R1
Reaction direction	Decrease
Unit	mmol/L (mEq/L)
Reagent volumes	
R1	200 µL
Sample	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

## UNIT CONVERSION

mmol/L = mEq/L

## EXPECTED VALUES<sup>8</sup>

### Serum, plasma:

Newborn:	13–22 mmol/L (mEq/L)
Infant, child:	20–28 mmol/L (mEq/L)
Adults:	23–29 mmol/L (mEq/L)
Adults >60r:	23–31 mmol/L (mEq/L)

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on [www.erba.com](http://www.erba.com).

**Limit of quantification:** 0.80 mmol/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

**Linearity:** 50 mmol/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

## Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mmol/L)	SD (mmol/L)	CV (%)	Intermediate precision	Mean (mmol/L)	SD (mmol/L)	CV (%)
Sample 1	32.2	0.19	0.60	Sample 1	29.9	0.96	3.23
Sample 2	15.5	0.16	1.02	Sample 2	14.3	0.42	2.91

## Accuracy

Two different validated control materials were used. Determined bias is and 0.5 % at the target value 20.0 mmol/L and 1.7 % at the target value 30.0 mmol/L.

## Comparison

A comparison between XL-640 automatic system BICARBONATE (y) and a commercially available test (x) using 103 samples gave following results:

Linear regression:

$$y = 0.970x - 0.041 \text{ mmol/L} \quad r = 0.995$$

Passing-Bablok<sup>9</sup>:

$$y = 0.968x - 0.002 \text{ mmol/L} \quad r = 0.983$$

## Interferences

Criterion: Recovery within ±10 % of initial value of  $\text{CO}_2$  concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 3.5 g/L, bilirubin up to 35 mg/dL, triglycerides up to 1650 mg/dL.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.

Quality of reagents is monitored on automatic systems ERBA XL by checking of the minimal permissible absorbance value of blank.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of  $\text{CO}_2$ . See paragraph Interferences.

Drugs: No interference was found at therapeutic concentrations using common drug panels<sup>10</sup>.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### R1, R2

Reagents of the kit are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Wilson W, Jesyk P, Rand R, et al. Use of Vickers discrete analyzer for enzymatic determination of the bicarbonate content of serum. Clin Chem, 19 (6), 640, 1973.
2. Norris KA, Atkinson AR, Smith WG, et al. Colorimetric enzymatic determination of serum total carbon dioxide, as applied to the Vickers Multichannel 300 discrete analyzer. Clin Chem, 21, 1093-1101, 1975.
3. Scott MG, Heusel JW, LeGrys VA, et al. Electrolytes and blood gases, in Tietz NW. Textbook of Clinical Chemistry. 3rd ed. Philadelphia, PA: WB Saunders Co, 1065-1066, 1999.
4. Müller-Plathe O. Acid base balance and blood gases. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 318-329, 1998.
5. Bergmeyer HU. Methods of enzymatic analysis, third edition, Vol. VII, 572, 1987.
6. Gambino SR, Schreiber H. The measurement of CO<sub>2</sub> content with the autoanalyzer. A comparison with 3 standard methods and a description of a new method (alkalinization) for preventing loss of CO<sub>2</sub> from open cups. Am J Clin Path, 45, 406, 1966.
7. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag, 18-19, 2001.
8. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11): 783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

## USED SYMBOLS / POUŽITÉ SYMBOLY / POUŽITÉ SYMBOLY / SÍMBOLOS UTILIZADOS / ВИКОРИСТАНІ ПОЗНАЧКИ / SYMBOLES UTILISÉS / SÍMBOLOS USADOS

<b>REF</b>	Catalogue number Katalogové číslo Katalógové číslo Número de catálogo Каталожний номер Numéro de catalogue Número de catálogo	<b>LOT</b>	Lot number Číslo šarže Číslo šarže Número de lote Номер партії Numéro de lot Número de lote		Expiry date Datum expirace Datum expirácie Fecha de caducidad Термін придатності Date d'expiration Data de validade		Consult instructions for use Čtěte návod k použití Čítajte návod k použitiu Consulte las instrucciones de uso Перед використанням уважно вивчіть інструкцію Consulter la notice d'utilisation Veja as instruções de uso
<b>IVD</b>	<i>In vitro</i> diagnostic medical device Diagnostický zdravotnícký prostriedek <i>in vitro</i> Diagnostický zdravotnícký prostriedok <i>in vitro</i> Dispositivo médico para diagnóstico <i>in vitro</i> <i>In vitro</i> діагностика Dispositif médical de diagnostic <i>in vitro</i> Diagnóstico <i>in vitro</i>		Manufacturer Výrobce Výrobca Fabricante Виробник Fabricant Fabricante		Temperature limit Omezení teploty Obmedzenie teploty Limite de temperatura Температура зберігання Limites de température Temperatura de armazenamento	<b>CONT</b>	Content Obsah Obsah Contenido Вміст Contenu Conteúdo

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