CREATINE KINASE MB

Cat. No.	Pack Name	Packaging (Content)
BLT00018	CK MB 100	R1: 4 × 20 mL, R2: 1 × 20 mL, instruction for use

EN

CE 2797 IVD

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of creatine kinase MB in human serum and plasma on various automatic systems. In combination with other parameters it is intended for monitoring and diagnosis of the people with a suspected myocardial infarction. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Creatine kinase (CK) appears as three isoenzymes which are dimers composed of two types of monomer subunits. The isoenzymes comprise all three combinations of monomers, M (for skeletal muscle derived) and B (for brain derived), as represented by the notations MM, MB, and BB. CK-MB is found in a high concentration in the myocardium (between 14 to 42 %) and to a lesser extent in skeletal muscle. Damage to the myocardium, as will occur in acute myocardial infarction (AMI), will result in increased circulating levels of the CK-MB isoform. Typically CK-MB levels become elevated 4 to 6 hours after the onset of chest pain, peak between 12–24 hours and return to baseline within 48 hours. Creatine kinase MB isoform (CK-MB) is less sensitive and less specific than some other biomarkers for the evaluation of myocardial injury¹.

PRINCIPLE

After immunoinhibition with antibodies to the CK-M subunit², the CK-B activity is determined with a standardized method for the determination of CK with activation by NAC as recommended by the German Society for Clinical Chemistry (DGKC)³ and the International Federation of Clinical Chemistry (IFCC)^{4,5} in 1977 and 2002 respectively. This assay meets the recommendations of the IFCC and DGKC. Specific antibodies against CK-M inhibit the complete CKMM activity and the CK-M subunit of CKMB. Only CK-B activity is measured. The remaining CK-B activity, corresponding to half the CK-MB activity, is determined by the total CK method.

As the CK-BB isoenzyme only rarely appears in serum and the catalytic activity of the CK-M and CK-B subunits hardly differ, the catalytic activity of the CK-MB isoenzyme can be calculated from the measured CK-B activity by multiplying the result by 2. The rate of formation of NADPH is therefore proportional to the activity of CK-MB present in the sample and can be measured kinetically at 340 nm.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Imidazole buffer, pH 6.1	125 mmol/L	Imidazole buffer, pH 8.9	125 mmol/L
Glucose	25 mmol/L	ADP	15.2 mmol/L
Magnesium acetate	12.5 mmol/L	D-glucose-6-phosphate-dehydrogenase	>8800 U/L
EDTA	2 mmol/L	Creatine phosphate	250 mmol/L
N-acetyl-L-cysteine	25 mmol/L	AMP	25 mmol/L
NADP ⁺	2.4 mmol/L	Diadenosine pentaphosphate	103 µmol/L
Hexokinase	>6800 U/L		
Anti-CK antibodies (goat) blocking capacity up to 2000 U/L CK-MB			

COMPOSITION OF REACTION MIXTURE

Imidazole buffer	120 mmol/L
Glucose	19.2 mmol/L
Magnesium acetate	9.6 mmol/L
EDTA	1.5 mmol/L
N-acetyl-L-cysteine	19.2 mmol/L
NADP ⁺	1.8 mmol/L
Hexokinase	>5230 U/L
Anti-CK antibodies (goat) blocking capacity up to 2000 U/L CK-MB	
ADP	2.9 mmol/L
D-glucose-6-phosphate-dehydrogenase	>1692 U/L
Creatine phosphate	48 mmol/L
AMP	4.8 mmol/L
Diadenosine pentaphosphate	19.8 µmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 340 nm may be used, general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method – substrate start

Reagents are ready to use. After the first opening the vials, reagents are stable for 30 days at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

Monoreagent method – sample start

Stability of working reagent: 2 days at 15–25 °C in dark
 2 weeks at 2–8 °C in dark

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.^{8,9,10}

Serum.

Plasma: Li-heparin and K₂-, K₃-EDTA plasma.

Li-heparin in the usual concentration does not interfere with the test, but IFCC warns against its use⁴.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum ¹¹ :	8 hours at	20–24 °C
	8 days at	2–8 °C
	4 weeks at	-20 °C

Stability in heparin plasma ¹¹ :	8 hours at	20–24 °C
	5 days at	2–8 °C
	8 days at	-20 °C

Stability in EDTA plasma ¹⁰ :	2 days at	20–25 °C
	7 days at	4–8 °C
	1 year at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank
 Calibration frequency: it is recommended to do a calibration
 • after reagent lot change
 • as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.
 The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method⁴ with addition of antibodies.

ASSAY PROCEDURE

Wavelength: 340 nm

Cuvette: 1 cm

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.040 mL
Calibrator	–	0.040 mL	–
Distilled water	0.040 mL	–	–

Mix and after 3 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix, incubate 3 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.040 mL
Calibrator	–	0.040 mL	–
Distilled water	0.040 mL	–	–

Mix, incubate 3 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

CALCULATION

$$1. \text{CK-MB (U/L)} = \frac{\Delta A_{\text{sam}} / \text{min.}}{\Delta A_{\text{cal}} / \text{min.}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \text{Using factor (f):} \quad \text{CK-MB (U/L)} = f \times \Delta A / \text{min} \quad f = 8254 \text{ (at 340 nm)}$$

Note: for CHEM 7 it is recommended to use factor 8360.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic	Reaction direction	Increasing
Wavelength (nm)	340	Normal Low U/L	0
Sample Volume (µL)	40	Normal High U/L	2.5
Working Reagent Volume (µL)	1000	Linearity Low U/L	5.23
Lag time (sec.)	180	Linearity High U/L	1200
Kinetic interval (sec.)	60	Blank with	Water
No. of readings	3	Absorbance limit (max.)	0.7
Kinetic factor	8254	Units	U/L
Reaction temperature (°C)	37		

UNIT CONVERSION

U/L × 0.0167 = µkat/L

EXPECTED VALUES^{12,13,14}

At 37 °C

<25 U/L

Myocardial infarction: the risk of myocardial infarction is high if following three conditions are fulfilled:

- CK (Men) >190 U/L
CK (Women) >167 U/L
- CK-MB >24 U/L
- CK-MB activity is between 6 and 25 % of total CK activity.

If myocardial infarction is suspected and the conditions are not fulfilled, the infarction may be fresh. In this case the measurements should be repeated after 4 hours with fresh samples. In healthy individuals different values are found depending on race and age^{13,14}. It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

5.23 U/L
 Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

Linearity:

1200 U/L
 Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	37.5	0.75	2.00
Sample 2	81.4	2.76	3.40

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	39.3	1.26	3.20
Sample 2	85.8	2.91	3.39

Comparison

A comparison between XL-640 automatic system CK-MB (y) and a commercially available test (x) using 48 samples gave following results:

Linear regression:

$$y = 1.005x - 0.62 \quad r = 0.998$$

Passing-Bablok¹⁵:

$$y = 1.016x - 0.77 \quad r = 0.925$$

Interferences

Critereon: Recovery within ±10 % of initial value of CK-MB activity in the sample without interfering substance. Following substances do not interfere: haemoglobin up to 4.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.
 Drugs: No interference was found at therapeutic concentrations using common drug panels¹⁶.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 340 nm against the distilled water is 0.7.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of CK-MB. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
 The Summary of Safety and Performance is available at Eudamed see: <https://ec.europa.eu/tools/eudamed>.

Health protection

Used animal-source materials are accompanied by an official state veterinarian health certificate. The animal components were collected from healthy animals under the care of a registered veterinarian and were apparently free from infectious and contagious diseases and injurious parasites. However, this material should be handled as thought capable of transmitting infectious diseases.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1
 UFI: X3UU-AWXN-QJ5K-WPF9

R2
 UFI: N6UU-UWN2-1J53-J11C
 R1, R2

Contains: imidazole
Hazard statement:
 H360D May damage the unborn child.
Precautionary statement:
 P201 Obtain special instructions before use.
 P280 Wear protective gloves / protective clothing / eye protection / face protection.
 P308 + P313 IF exposed or concerned: Get medical advice / attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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


1. Thygesen K, Alpert JS, Jaffe AS, Chaitman BR, Bax JJ, Morrow DA, White HD, The Executive Group on behalf of the Joint ESC/ACC/AHA/WHF Task Force for the Universal Definition of Myocardial Infarction; Fourth universal definition of myocardial infarction, 2018.
2. Würzburg U, Hennrich N, Lang H, et al. Determination of creatine kinase-MB in serum using inhibiting antibodies. *Klin Wschr* 1976;54(8):357-360.
3. Bergmeyer HU, Breuer H, Büttner H, et al. Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Standard-Methode zur Bestimmung der Aktivität der Creatin-Kinase. *J Clin Chem Clin Biochem* 1977;15:249-254.
4. Hørdler M, Elser RC, Gerhardt W, et al. IFCC methods for the measurement of catalytic concentration of enzymes. Provisional recommendation IFCC method for creatine kinase Appendix A. *J Int Fed Clin Chem* 1990;2:26-35.
5. Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine Kinase. *Clin Chem Lab Med* 2002;40(6):635-642.
6. Moss DW, Henderson AR, Kachmar JF. Enzymes. In: Tietz NW, ed. *Fundamentals of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders 1987;346-421.
7. Lott JA, Stang JM. Serum enzymes and isoenzymes in the diagnosis and differential diagnosis of myocardial ischemia and necrosis. *Clin Chem* 1980;26:1241-1250.
8. Castaldo AM, Ercolini P, Forino F, Basevi A, Vrenna L, Castaldo P, D'Ambrosio V, Castaldo A, Plasma Myoglobin in the Early Diagnosis of Acute Myocardial Infarction, *Eur J Clin Chem Clin Biochem* 1994; 32, 349-353.
9. Dominici R, Infusino I, Valente C, Moraschinelli I, Franzini C, Plasma or serum samples; measurements of cardiac troponin T and other analytes compared, *Clin Chem Lab Med* 2004; 42(8), 945-951.
10. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
11. Braun S, Rösenthaller F, Jarausch J, et al. Analyte Stability of CK-MB Activity and cTnT in ICU Patient Serum and Heparin Plasma. Poster presented at Medica 2004, Düsseldorf. (Roche Diagnostics GmbH No. 04587979990).
12. Stein W. Creatine kinase (total activity), creatine kinase isoenzymes and variants. In: Thomas L, ed. *Clinical laboratory diagnostics*. Frankfurt: TH-Books Verlagsgesellschaft;1998.p.71-80
13. Stein W. Strategie der klinisch-chemischen Diagnostik des frischen Myokardinfarkts. *Med Welt* 1985;36:572-577.
14. Myocardial infarction redefined – a consensus document of the Joint European society of Cardiology/America College of Cardiology Committee for the redefinition of myocardial infarction. *Eur Heart J* 2000;21:1502-1513.
15. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
16. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.

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