Erba H7100 Dil Haematology Analyser Diluent

Cat. No.	Product Name	Packaging (Content)
HEM00044	Erba H7100 Dil	20 L, Instruction for use, RFID card

INTENDED USE

Erba H7100 Dil is an *in vitro* diagnostic reagent for automated dilution of human blood for counting and sizing of blood cells on haematology analysers H7100. It is designed to report on the haematological status of patient samples, supporting the monitoring and diagnosis of pathological conditions. For professional use in clinical laboratory only.

PRINCIPLE

The diluent is a solution with certain ion strength and conductivity, which can dilute the blood and form the sheath flow, providing a stable environment for blood cell counting.

DESCRIPTION AND COMPOSITION

Sodium Chloride; Sodium Sulphate; Buffering Agent; Antibacterial agent; Preservative.

PREPARATION

This reagent is ready for use and can be applied straight from the container; no special reagent preparation is necessary. Avoid reagent contaminations.

Before using, leave the reagent at room temperature for at least 24 hours.

For reagent connection to the analyser, refer to the H7100 User manual.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser H7100, Cat. No. INS00087

Reagents: Erba H7100 Dil Ret, Cat. No. HEM00045, HEM00065

Erba H7100 Lyse1, Cat. No. HEM00046, HEM00047, HEM00066

Erba H7100 Lyse2, Cat. No. HEM00048, HEM00049, HEM00050, HEM00067

Erba H7100 Fluro Ret, Cat. No. HEM00051, HEM00052, HEM00068

Erba H7100 Fluro Diff, Cat. No. HEM00053, HEM00054

SPECIFICATION

The product is a clear liquid without sediment, particles or floccules. Background result: WBC $\leq 0.1 \times 10^9/L$, RBC $\leq 0.02 \times 10^{12}/L$, HGB $\leq 1 \text{ g/L}$, PLT $\leq 5 \times 10^9/L$, HCT $\leq 0.5 \%$.

STABILITY AND STORAGE

The product is stable till the expiry date mentioned on the label when stored at (+2 to +30) °C. Discard opened container after 60 days.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: W4NN-8388-XN46-257Y

Warning

Contains: formaldehyde, reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)

Hazard statement:

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Precautionary statement:

P261 Avoid breathing mist/vapours/spray.

P273 Avoid release to the environment.

P280 Wear protective gloves.

P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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