# URINORM XL

Cat. No.	Product Name	Packaging (Content)
REG00060	URINORM XL	URINORM XL N: 1 × 120 mL, URINORM XL P: 1 × 120 mL, instruction for use

## INTENDED USE

URINORM XL control urines are liquid synthetic material with concentrations of analytes and urine sediment within normal (URINORM XL N) and pathological (URINORM XL P) levels. URINORM XL control urines are designed for verification and confirmation of precision and accuracy of PHAN AUTO diagnostic urine strips as well as LAURA XL urine analyser intended for *in vitro* diagnostics. This control is also intended to evaluate microscopic examination test procedure. For professional use at clinical laboratories only.

## PRINCIPLE

After contact of sample with individual pads, diagnostics zones change their colour. The intensity and shade of zone coloration are proportional to the concentration of analyte. The URINORM XL N simulates normal levels of analytes and its reaction with diagnostic pad is negative. URINORM XL P simulates pathological state and provides positive reaction of all diagnostic pads.

## DESCRIPTION AND COMPOSITION

URINORM XL control urines are human urine-based containing constituents of human and animal origin as well as purified chemicals. URINORM XL N simulates normal values of individual analytes. URINORM XL P simulates pathological levels of each analyte and provide positive reaction for all diagnostic pads. URINORM XL control urines contain preservatives added to maintain product stability and avoid the microbial growth.

URINORM XL N – normal level 120 mL

URINORM XL P – pathological level 120 mL

## PREPARATION

URINORM XL control urines are ready-to use without any special preparation of sample. Do not add any preservatives to the control urine.

## MATERIAL REQUIRED BUT NOT PROVIDED

Diagnostic strips: DEKAPHAN AUTO; Cat. No. URP0030  
UNDEKAPHAN AUTO; Cat. No. URP0031

Automatic urine analyser: LAURA XL; Cat. No. INS00065.

## STABILITY AND STORAGE

URINORM XL control urines are stable until the expiration date unopened stored at 2–8 °C and protected from light and microbial contamination.

Once opened URINORM XL control urines are stable for 30 days at 20–25 °C and 18 months at 2–8 °C.

## METROLOGICAL TRACEABILITY

The QC reagents used with each analyser are traceable to international standards. Thus, the reference method, materials, and values are found to be specific to each instrument/analyser. Since the controls are value assigned using these reagents they are, therefore, traceable to the same reference standards. The value assignment of URINORM XL is done for every freshly produced lot by manufacturer.

## ASSAY PROCEDURE

Control urines URINORM XL should be treated the same as patient specimens.

1. Allow the control urines bottles to reach the room temperature.
2. Gently invert the bottles to homogenize the sample.
3. Remove the cap and pour 1–1.5 mL into the Eppendorf tube.
4. Perform the QC measurement on LAURA XL urine analyser – device automatically evaluates the results and compare them with declared values.
5. Close the bottle after use and store appropriately at 2–8 °C

## DECLARED VALUES

The following table summarizes declared values for each analyte for automatic urine analyser LAURA XL.

LOT	XXXXXXX
EXP	YYYY-MM

CHEMISTRY			
SCALE 1			
Analyte	URINORM XL N	URINORM XL P	Units
Specific gravity	X-X	X-X	
Nitrites	X	X	
pH	X-X	X-X	
Protein	X	X-X / X-X / X-X	g/L / mg/dL / ARB
Ascorbic acid	X	X-X / X-X / X-X	mmol/L / mg/dL / ARB
Glucose	X	X-X / X-X / X-X	mmol/L / mg/dL / ARB
Ketones	X	X-X / X-X / X-X	mmol/L / mg/dL / ARB
Urobilinogen	X	X-X / X-X / X-X	μmol/L / mg/dL / ARB
Bilirubin	X	X-X / X-X / X-X	μmol/L / mg/dL / ARB
Leucocytes	X	X-X / X-X / X-X	Leu / μL / ARB
Blood	X	X-X / X-X / X-X	Ery / μL / ARB
SEDIMENT			
Red blood cells	X-X / X-X	X-X / X-X	N/hpf / N/μL
White blood cells	X-X / X-X	X-X / X-X	N/hpf / N/μL

CHEMISTRY			
SCALE 2			
Analyte	URINORM XL N	URINORM XL P	Units
Specific gravity	X-X	X-X	
Nitrites	X	X	
pH	X-X	X-X	
Protein	X	X-X / X-X / X-X	g/L / mg/dL / ARB
Ascorbic acid	X	X-X / X-X / X-X	mmol/L / mg/dL / ARB
Glucose	X	X-X / X-X / X-X	mmol/L / mg/dL / ARB
Ketones	X	X-X / X-X / X-X	mmol/L / mg/dL / ARB
Urobilinogen	X	X-X / X-X / X-X	μmol/L / mg/dL / ARB
Bilirubin	X	X-X / X-X / X-X	μmol/L / mg/dL / ARB
Leucocytes	X	X-X / X-X / X-X	Leu / μL / ARB
Blood	X	X-X / X-X / X-X	Ery / μL / ARB
SEDIMENT			
Red blood cells	X-X / X-X	X-X / X-X	N/hpf / N/μL
White blood cells	X-X / X-X	X-X / X-X	N/hpf / N/μL

## EVALUATION

The result for each analyte is automatically evaluated by LAURA XL. The results are automatically compared with declared values set in LAURA XL user interface and the overall QC result is displayed.

## EXPECTED VALUES

Based on literature (whole population)<sup>1</sup>:

**Blood:** up to 5 Ery/μL

**Leucocytes:** up to 10 Leu/μL

**pH:** 5.5–7.0

**Urobilinogen:** up to 1 mg/dL

**Ketones:** up to 0.19 mmol/L

**Ascorbic acid:** not stated

**Glucose:** up to 1.4 mmol/L

**Nitrite:** 10<sup>5</sup> bacteria/mL

**Bilirubin:** up to 17 μmol/L

**Protein:** over 30 mg/24 hours

**Specific gravity:** 1.015–1.025

Each laboratory should verify that it is appropriate to apply these values to its patient population and if necessary, determine its own ranges.

## PERFORMANCE CHARACTERISTICS

Precision for URINORM XL control urines is defined as agreement with declared values in instruction for use for each lot.

Intra-assay: Two levels of control material were used as samples – 20 replicates in 2 consecutive runs for each level was performed in one day.

Inter-assay: Two levels of control material were used as samples – 2 replicates from each level were measured twice per day during 20 consecutive days.

Parameter	Intra-assay		Inter-assay	
	Level N	Level P	Level N	Level P
	Agreement (%)	Agreement (%)	Agreement (%)	Agreement (%)
SG	100	100	100	100
LEU	100	100	100	100
NIT	100	100	100	100
pH	100	100	100	100
ASC	100	100	100	100
PRO	100	100	100	100
GLU	100	100	100	100
KET	100	100	100	100
UBG	100	100	100	100
BIL	100	100	100	100
BLD	100	100	100	100
RBC	100	100	100	100
WBC	100	100	100	100

## LIMITATIONS

Control urines should not be used past the expiration date on the label. Do not use turbid control urines where there is an evidence microbial contamination.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Contains human sourced material. Because no known test method can offer complete assurance that products derived from human sources will not transmit infectious agents, the product should be handled as potentially infectious.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### URINORM XL N

UFI: EA5Y-MJYP-0E7C-JFT3

Warning

Control urines contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1).

### Hazard statement:

H317 May cause an allergic skin reaction.

### Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

### URINORM XL P

Reagent is not classified as hazardous.

## WASTE MANAGEMENT

Residues of control urines should be handled as potentially infectious material and should be disposed in accordance with internal regulations and by-laws. Paper and other packaging materials are to be discarded in the form of sorted waste (paper, glass, plastic).

## REFERENCES

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.

Date of revision: 9. 4. 2026

## USED SYMBOLS

Catalogue number

*In vitro* diagnostic medical device

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