PHAN® LAURA

| | | | | | | | | | | | | | |
|--|------------|---|----|-----|-----|----|-----|-----|-----|-----|-----|-----|----|
|  | | | | | | | | | | | | | |
| DIAGNOSTIC TEST STRIPS FOR URINALYSIS | | | | | | | | | | | | | |
| Product name | REF |  | SG | LEU | NIT | pH | PRO | GLU | KET | UBG | BIL | BLD | CP |
| DEKAPHAN® LAURA | URPH0028 | 100 | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ | ■ |
|  | | | | | | | | | | | | | |

INTENDED USE

Diagnostic strips for *in vitro* semiquantitative determination of specific gravity, leucocytes, nitrites, pH, protein, glucose, ketones, urobilinogen, bilirubin and blood in urine on readers LAURA and LAURA SMART. Intended for screening and monitoring of kidney function disturbances, leukocyturia, bacteriuria, acidosis / alkalosis, proteinuria, glucosuria, ketonuria, urobilinogenuria, bilirubinuria and haematouria. For professional use at clinical laboratories only. The semiquantitative analysis is not sufficient for the completing of diagnoses.

DESCRIPTION

For objective evaluation there are available 10 indication zones (specific gravity, leucocytes, nitrites, pH, protein, glucose, ketones, urobilinogen, bilirubin and blood) that are completed in various combinations as polyfunctional strips – all types of available products are mentioned in the table above. The diagnostic strips are intended for single use only.

SPECIMEN COLLECTION AND HANDLING

For specimen collection and preparation use only suitable clean and dry tubes or containers without traces of detergents and disinfectants.

Do not add the preservatives to the urine.

Use freshly voided urine, well mixed that has not been centrifuged. The urine should not be more than 2 hours old when tested. In case of longer standing, mix before use.

Do not expose the sample into the direct sunlight as this may cause the degradation of bilirubin and urobilinogen and this may turn in lower or false negative results.

Mind that woman menstruation blood may affect the results.

ASSAY PROCEDURE:

1. Remove only as many test strip required and reseal the tube immediately after use. Do not touch test pads of the strip. Do not divide test strips!
2. Immerse completely all reagent pads in urine specimen (for 1–2 seconds).
3. Run edge of the strip against rim of urine container to remove excess urine and wipe off the excess urine from the strip to the cellulose cotton wool as well. Hold strip in horizontal position.
4. Place the diagnostic strip to any LAURA reader. The strip is evaluated automatically.

EVALUATION

After immersion, while in contact with individual urine components, diagnostics zones change their colour. The intensity and shade of zone coloration are proportional to the concentration of analyte. While evaluated using the LAURA or LAURA SMART reader the intensity of reflected light is measured. The processing unit converts the reflected light intensity to the analytical value. The result is presented on the display and printed by the built-in thermal printer.

WORKING REAGENT CONCENTRATIONS

Specific gravity: poly(methylvinylether/maleic acid) 32 %; bromthymol blue 5.1 %

Leucocytes: indoxyl ester 0.43 %; diazonium salt 0.05 %

Nitrite: sulphaniamide 5.1 %; tetrahydrobenzo-[h]-quinoline 5.8 %

pH: methyl red 0.71 %; bromthymol blue 12.1 %

Protein: tetrabromphenolphthalein ester 0.21 %; tetrabromphenol blue 0.35 %

PRINCIPLE

Specific gravity (SG) – The test is based on the principle of ion exchange, which runs between polyelectrolyte and ions present in urine. Its result is colour change of acid-base indicator from blue-green colour in urine with low concentration of ions, through green and yellow-green in urine with increased concentration of ions, to amber yellow colour.

Leucocytes (LEU) – The test is based on enzymatic reaction. The test pad contains an indoxyl ester that is cleaved by granulocyte esterases. The released indoxyl reacts with a diazonium salt and pink or violet coloration is formed.

Nitrites (NIT) – The test is based on conversion of nitrate to nitrite by the action of certain species of bacteria contained in urine. The colour test is based on the principle of the Griess test. Any degree of pink coloration should be interpreted as a positive nitrite test.

pH – The test is based on the double acid-base indicator principle and gives a range of colours from orange through yellow and green to blue.

Protein (PRO) – The test is based on colour change of acid-base indicator, which is caused by presence of proteins. It is particularly sensitive to albumin, but is much less sensitive to globulin, mucoprotein, haemoglobin and Bence-Jones protein.

Glucose (GLU) – The test is based on the specific glucose oxidase/peroxidase reaction. The test is specific only for D-glucose.

Ketones (KET) – The test is based on the principle of Legal's test and is more susceptible to acetoacetic acid than to acetone. Test does not react with β-hydroxybutyric acid.

Urobilinogen (UBG) – The test is based on the coupling of urobilinogen with stabilized reagent. The test is specific for urobilinogen and stercobilinogen and is not susceptible to the interfering factors known at Ehrlich's test.

Bilirubin (BIL) – The test is based on the coupling of bilirubin with stabilized reagent. Even slightly orange coloration of the diagnostic pad means pathologic result.

Blood (BLD) – The test is based on the peroxidase activity of haemoglobin which catalyzes the oxidation of the indicator due to the presence of the organic hydroperoxide contained in the diagnostic pad. The label contains two colour scales; for detection of intact erythrocytes and free haemoglobin. The test is highly sensitive to free haemoglobin.

Compensation zone – The strip of a filter paper which is not impregnated by any agents and serves for suppression of impact of dark urines during objective assessment.

LIMITATIONS

Specific gravity – This reaction is affected by pH values of urine over 6.5, this can cause shift in colour response towards lower values of specific gravity.

Leucocytes – In case a urine sample is more markedly coloured (e.g. increased bilirubin content), the resulting colour could be affected by a sample coloration. The intensity of colour reaction is increased by alkaline pH and higher urine density.

Nitrite – Before testing the patient should intake vegetable-rich meals and discontinue antibiotic therapy and vitamin C for 3 days prior to the test. Sensitivity of this test declines with high specific gravity of urine. Increased diuresis can cause false negative results. Limited fluid intake prior testing can prevent from the excessive dilution of urine. Test can be applied only at fresh urine. Inaccurate results may occur at stale urines, in which nitrite can be formed by contamination of the specimen.

Protein – In strongly alkaline urines (pH >8) from patients on medication with quinine or quinoline containing drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quarternary ammonium groups. On the other hand, in the presence of non-ionic or anionic detergents, false-negative results may occur. Do not take the colour of the dry pad into consideration.

Glucose – The reaction is independent on pH and presence of ketone bodies.

Ketones – Drugs and diagnostics on the basis of phenolphthalein or sulphophthalein may turn red to purple because of alkaline reaction of the pad.

Urobilinogen – The reaction is not affected by pH of urine. The presence of bilirubin gives yellow colour. This colour which turns slowly to greenish-blue does not interfere with urobilinogen determination provided the reading is made 1 minute after wetting. The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of urobilinogen and thus leads to artificially low or false negative readings.

Bilirubin – The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of bilirubin and thus leads to artificially low or false negative readings. High concentrations of urobilinogen (above 100 µmol/L) interfere with the test. Also the red substances or the substances, which are turning red in low pH may interfere (e.g. Phenazopyridine). The reaction is not affected by pH of urine.

Blood – Microbial peroxidase associated with urinary tract infection may cause a false positive reaction. The sensitivity of the test is influenced by specific gravity or by inhibitors of pharmacological origin.

All diagnostics pads do not interfere with the common concentration of the ascorbic acid.

PLEASE NOTE

Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after discontinuing a drug. The sensitivity depends upon the variability of urines. The semiquantitative analysis is not sufficient for the completing of diagnoses.

MATERIAL REQUIRED BUT NOT PROVIDED

Control urines URINORM (Cat.No. REG00053), or Bio-Rad Liquichek, LAURA (Cat.No. INS00063), or LAURA SMART reader (Cat.No. INS00064).

QUALITY CONTROL

For quality control of precision and accuracy URINORM control urines, Cat.No. REG00053 are recommended. Obtained results are compared with the Target values in the instruction for use of URINORM. Bio-Rad Liquichek Urinalysis Control can be also used for quality control. It is recommended to perform QC measurements according to your local laboratory guidelines. More information about quality control can be found in LAURA and LAURA SMART user manuals.

STABILITY AND STORAGE

Keep diagnostic test strips in tightly closed original tubes in a dry and dark place at 2–30 °C. The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory. When stored under these conditions, test strips are stable to the expiry date given on the pack. Once opened, the diagnostic strips are stable for at least 3 months.

ANALYTICAL PERFORMANCE

Limit of detection is defined as the lowest concentration of the analyte which leads to positive results in ≥90 % cases. Limit of detection is not applicable for SG and pH.

Method comparison of PHAN LAURA diagnostic strips using LAURA SMART and LAURA urine reader was performed against COBAS U411 analyzer and Combur 10 test M. All concentration ranges were covered.

| Parameter | Result Values | Limit of detection | Method Comparison LAURA SMART | Method Comparison LAURA |
|------------|---|--------------------|--|---|
| SG | 1.000, 1.005, 1.010, 1.015, 1.020, 1.025, 1.030 | NA | 89 %* | 87 %* |
| LEU | neg., 25, 75, 500 LEU/µL neg., 1+, 2+, 3+ | 20/25 LEU/µL** | neg. 99 % pos. 98 % | neg. 83 % pos. 98 % |
| NIT | neg., pos. | 0.8 mg/L | neg. 100 % pos. 100 % | neg. 99 % pos. 100 % |
| pH | 5, 6, 6.5, 7, 8, 9 | NA | 95 %* | 98 %* |
| PRO | neg., 30, 100, 500 mg/dL neg., 0.30, 1, 5 g/L neg., 1+, 2+, 3+ | 0.22 g/L | neg. 96 % pos. 100 % | neg. 98 % pos. 100 % |
| GLU | neg., 50, 100, 300, 1000 mg/dL neg., 2.8, 5.5, 17, 55 mmol/L neg., 1+, 2+, 3+, 4+ | 1.7 mmol/L | neg. 99 % pos. 91 % | neg. 99 % pos. 93 % |
| KET | neg., 5.2, 16, 52, 156 mg/dL neg., 0.5, 1.5, 5, 15 mmol/L neg., ±, 1+, 2+, 3+ | 0.3 mmol/L | neg. 100 % pos. 89 % | neg. 99 % pos. 85 % |
| UBG | norm., 1, 3, 6, 12 mg/dL norm., 17, 51, 102, 203 µmol/L norm., 1+, 2+, 3+, 4+ | 0.8 mg/dL | neg. 100 % pos. 76 % | neg. 100 % pos. 84 % |
| BIL | neg., 1, 3, 6 mg/dL neg., 17, 51, 102 µmol/L neg., 1+, 2+, 3+ | 0.8 mg/dL | neg. 100 % pos. 94 % | neg. 99 % pos. 100 % |
| BLD | neg., 10, 50, 250 ERY/µL neg., 1+, 2+, 3+ | 7 ERY/µL | neg. 99 % pos. 98 % | neg. 99 % pos. 98 % |

* The acceptance criteria for SG and pH are that the satisfactory results may lay in ±1 concentration field.

** LOD for leucocytes using LAURA is 20 Leu/µL and LOD for leucocytes using LAURA SMART is 25 Leu/µL.

Precision

Precision testing consists of Intra-assay and Inter-assay measurements using control urines material.

Intra-assay - Two levels of control urine material was used as a sample. Each level was measured 20 times in one day.

Inter-assay - Two levels of control urine material was used as a sample. Each level was measured 2 times in one day and this testing was repeated in 20 consecutive days.

For both assays control urines BIORAD Liquichek was used as a control.

| | | Precision Intra-assay | | | | Precision Inter-assay | | | |
|---------------------|----------------|-----------------------|-------|------------------------|------------------------|-----------------------|-------|------------------------|------------------------|
| | | Result | | Agreement | | Result | | Agreement | |
| Parameter | Control | Laura SMART | Laura | Laura SMART | Laura | Laura SMART | Laura | Laura SMART | Laura |
| SG | Level 1 | 1.015 | 1.010 | 100 %* | 100 %* | 1.015 | 1.010 | 100 %* | 100 %* |
| | Level 2 | 1.010 | 1.005 | 100 %* | 100 %* | 1.010 | 1.005 | 100 %* | 100 %* |
| LEU (Leu/µL) | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | 500 | 500 | 100 % | 100 % | 500 | 500 | 100 % | 100 % |
| NIT | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | pos. | pos. | 100 % | 100 % | pos. | pos. | 100 % | 100 % |
| pH | Level 1 | 6 | 6 | 91 % | 84 % | 6 | 6 | 90 % | 100 % |
| | Level 2 | 7 | 7 | 100 % | 96 % | 7 | 7 | 96 % | 100 % |
| PRO (mg/dL) | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | 100 | 100 | 100 % | 100 % | 100 | 100 | 100 % | 100 % |
| GLU (mg/dL) | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | 300 | 300 | 100 % | 95 % | 300 | 300 | 100 % | 100 % |
| KET (mg/dL) | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | 16 | 16 | 100 % | 100 % | 16 | 16 | 100 % | 100 % |
| UBG (mg/dL) | Level 1 | norm. | neg. | 100 % | 100 % | norm. | norm. | 100 % | 100 % |
| | Level 2 | 6 | 6 | 100 %** | 100 %** | 6 | 6 | 100 %** | 100 %** |
| BIL | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | 6 | 6 | 100 % | 100 % | 6 | 6 | 100 % | 100 % |
| BLD (Ery/µL) | Level 1 | neg. | neg. | 100 % | 100 % | neg. | neg. | 100 % | 100 % |
| | Level 2 | 250 | 250 | 100 % | 100 % | 250 | 250 | 100 % | 100 % |

* The acceptance criteria for SG are that the satisfactory results may lay in ±1 concentration field.

** The acceptance criteria for UBG are that the satisfactory results may lay in ±1 concentration field due to different coloration of the pad than on the tube.

EXPECTED VALUES

Based on literature.^[1]

Blood: up to 5 Ery/µL

Leucocytes: up to 10 Leu/µL

pH: 5.5–7.0

Urobilinogen: up to 1 mg/dL

Ketones: up to 0.19 mmol/L

Glucose: up to 1.4 mmol/L

Nitrite: 10⁸ bacteria/mL

Bilirubin: up to 17 µmol/L

Protein: over 30 mg/24 hours

Specific gravity: 1.015–1.025

Each laboratory should verify that it is appropriate to apply these values to its patient population and if necessary determine its own ranges.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Do not touch test pads of the strip.

The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The product is not classified as dangerous.

WASTE MANAGEMENT

Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

REFERENCES


1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.

USED SYMBOLS

 Catalogue number

 Manufacturer

 Do not re-use

 Lot number

 Temperature limit

 Expiry date


 eIFU:

www.erba.com

 *In vitro* diagnostic medical device

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