MICROALBUPHAN[®] LAURA

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INTENDED USE

Diagnostic strips for *in vitro* semiquantitative determination of albumin and creatinine in urine on readers LAURA and LAURA SMART. Intended for screening and monitoring of risk of the beginning kidney damage. For professional use at clinical laboratories only. The semi-quantitative analysis is not sufficient for the completing of diagnoses.

DESCRIPTION

Diagnostic strips for objective *in vitro* semiquantitative determination of albumin and creatinine and albumin-creatinine ratio. The diagnostic strips are intended for single use only.

SPECIMEN COLLECTION AND HANDLING

For specimen collection and preparation use only suitable clean and dry tubes or containers without traces of detergents and disinfectants.

Do not add the preservatives to the urine.

Use freshly voided urine, well mixed that has not been centrifuged. The urine should not be more than 2 hours old when tested. In case of longer standing, mix before use.

Do not expose the sample into the direct sunlight as this may cause the degradation of bilirubin and urobilinogen and this may turn in lower or false negative results.

Mind that woman menstruation blood may affect the results.

ASSAY PROCEDURE

For objective evaluation of test strips, please, use the user manual for urine analyser LAURA or LAURA SMART.

- Remove only as many test strips as required and reseal the tube immediately after use. Do not divide test strips! Do not touch test pads of strips.
- Immerse completely all reagent pads in specimen (for 1–2 seconds).
- Run edge of the strip against rim of urine container to remove excess urine and wipe off the excess urine from the strip to the cellulose cotton wool as well. Hold strip in horizontal position.
- Evaluate the result using the reader LAURA or LAURA SMART.
- For the correct interpretation of the results take into account the ratio of albumin / creatinine.

EVALUATION

After immersion, while in contact with individual urine components, diagnostics zones change their colour. The intensity and shade of zone coloration are proportional to the concentration of analyte. While evaluated using the LAURA or LAURA SMART reader the intensity of reflected light is measured. The processing unit converts the reflected light intensity to the analytical value. The result is presented on the display and printed by the built-in thermal printer.

CALIBRATION

For more accurate results we recommend to calibrate for both diagnostic pads. Six-digit calibration code which is entered into the urine analyser LAURA or LAURA SMART is noted on the label, always under the lot number.

A detailed calibration procedure is mentioned in user manual for urine analyser LAURA or LAURA SMART.

WORKING REAGENT CONCENTRATIONS

Creatinine: 3,5-dinitrobenzoic acid 4.4 %, Albumin: tetrabromphenol blue 0.17 %

PRINCIPLE

Creatinine – The test is based on reaction of creatinine with 3,5-dinitrobenzoic acid in alkaline medium (Benedict-Behre reaction).

Albumin – The test is based on colour change of acid basic indicator, which is caused by presence of proteins (so-called proteins error of indicator). It is particularly sensitive to albumin but is much less sensitive to other urinary proteins.

LIMITATIONS

Creatinine: In a low positive sample, blood can cause falsely elevated results.

PLEASE NOTE

Knowledge of the effects drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after discontinuing a drug. The sensitivity depends upon variability of urines. The semi-quantitative analysis isn't sufficient for the completing of diagnoses.

MATERIAL REQUIRED BUT NOT PROVIDED

Control urines URINORM (Cat. No. REG00053), or Bio-Rad Liquichek, LAURA (Cat. No. INS00063), or LAURA SMART reader (Cat. No. INS00064).

QUALITY CONTROL

For quality control of precision and accuracy URINORM control urines, Cat.No. REG00053 are recommended. Obtained results are compared with the Target values in the instruction for use of URINORM. Bio-Rad Liquichek Urinalysis Control can be also used for quality control. It is recommended to perform QC measurements according to your local laboratory guidelines. More information about quality control can be found in LAURA and LAURA SMART user manuals.

STABILITY AND STORAGE

The recommended storage temperature is between (+2 to +25) °C. After first opening, the tube with remaining urine test strips has to be closed tightly immediately after removal of corresponding number of strips for evaluation. When stored under these conditions, test strips are stable to the expiry date given on the pack. Once opened, the diagnostic strips are stable for at least 3 months.

ANALYTICAL PERFORMANCE

Method comparison of MICROALBUPHAN diagnostic strips using LAURA SMART and LAURA urine reader was performed against Clinitek Status by Siemens. Agreement ±1 field.

Precision – Precision testing consists of Intra-assay and Inter-assay measurements using control urines material.

Repeatability (Intra-assay) – Two levels of control urine material were used as samples. Each level was measured 20 times in one day.

Reproducibility (Inter-assay) – Two levels of control urine material were used as samples. Each level was measured 2 times in one day and this testing was repeated in 20 consecutive days.

For both assays control urines BIORAD Liquichek was used as a control.

| Parameter | Range of assay | Method comparison | Repeatability | Reproducibility |
|------------|-----------------------------------|---------------------|-----------------------|----------------------|
| Albumin | 10–5000 mg/L (0,01–5 g/L) | 99 % | 100 % | 100 % |
| Creatinine | 0,1 →> 3 g/L (0,9 →> 26,5 mmol/L) | 90 % | 99,6 % | 100 % |

| Parameter | Control | Precision Intra-assay | | | | Precision Inter-assay | | | |
|------------------|----------------|-----------------------|-----------|-----------------------|----------------------|-----------------------|-----------|----------------------|-----------------------|
| | | Result | Agreement | | | Result | Agreement | | |
| Creatinine [g/L] | Level 1 | 0,25–2,00 | 0,25–2,00 | 99,6 % | 98 % | 0,25–1,0 | 0,25–2,0 | 85 % | 97,7 % |
| | Level 2 | 3 | 3 | 100 % | 100 % | 2–3 | > 3 | 100 % | 100 % |
| | Level 1 | 10 | 10 | 100 % | 100 % | 10–30 | 10–30 | 100 % | 100 % |
| Albumin [mg/L] | Level 1 | 1000 | 1000–5000 | 100 % | 100 % | 300–5000 | 300–5000 | 100 % | 100 % |
| | Level 2 | 1000 | 1000–5000 | 100 % | 100 % | 300–5000 | 300–5000 | 100 % | 100 % |

EXPECTED VALUES

Based on literature:

Albumin¹: <0.02 g/L (2 mg/dL)
Creatinine²: 2,2–26,5 mmol/L (0,25–3,0 g/L)

Albumin-to-Creatinine Ratio³: Albumin is normally present in urine at concentrations of less than 3.4 mg albumin / mmol creatinine (0.03 g albumin/g creatinine). Microalbuminuria is indicated at ratio results of 3.4–33.9 mg/mmol (30–300 mg/g), high values of ratio are indicated clinical albuminuria (proteinuria).

Albumin/Creatinine (mg/mmol) =

Albumin
(
g
/
L
)
×
1000

Creatinine
(
m
m
o
l
/
L
)

=

Albumin
(
m
g
/
L
)

Creatinine
(
m
m
o
l
/
L
)

Albumin/Creatinine (mg/g) =

Albumin
(
g
/
L
)
×
1000

Creatinine
(
g
/
L
)

=

Albumin
(
m
g
/
L
)

Creatinine
(
g
/
L
)

Albumin / Creatinine ratio is calculated by urine analyser LAURA or LAURA SMART automatically for each measured urine sample.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Do not touch test pads of the strip.

The strips must be kept away from moisture, direct sunlight, elevated temperature, and chemical fumes in the laboratory.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The product is not classified as dangerous.

WASTE MANAGEMENT

Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe hand- ling of such materials. Let waste recycle or put it to municipal waste.

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