
PHAN [®] AUTO													
			DIAGNOSTIC STRIPS FOR URINALYSIS FOR LAURA XL										
DIAGNOSTICKÉ PROUŽKY K VÝŠETŘENÍ MOČE ANALYZÁTORY LAURA XL													
		SG	LEU	NIT	pH	ASC	PRO	GLU	KET	UBG	BIL	BLD	CP
DEKAPHAN [®] AUTO	URPH0030	100	■	■	■	■	■	■	■	■	■	■	■
UNDEKAPHAN [®] AUTO	URPH0031	100	■	■	■	■	■	■	■	■	■	■	■
RFID The kind of strips applicable for automatic urine analyser LAURA XL only. / Proužky určené pouze pro přístroj LAURA XL [®] .													
													

INTENDED USE
Diagnostic strips for *in vitro* semiquantitative determination of specific gravity, leucocytes, nitrites, pH, ascorbic acid, protein, glucose, ketones, urobilinogen, bilirubin and blood in urine on automatic analyser LAURA XL. Intended for screening and monitoring of kidney function disturbances: leucocyturia, bacteriuria, acidosi/alkalosis, presence of ascorbic acid, proteinurie, glucosuria, keturia, urobilinogenuria, bilirubinuria and haematuria. Ascorbic acid in urine may affect the determination of other parameters (blood, glucose, nitrite). For professional use at clinical laboratories only. The semiquantitative analysis is not sufficient for the completing of diagnoses.

DESCRIPTION
For objective evaluation there are available 10 or 11 (depends on strip type) indication zones (specific gravity, leucocytes, nitrites, pH, ascorbic acid, protein, glucose, ketones, urobilinogen, bilirubin and blood) that are completed in various combinations as polyfunctional strips – all types of available products are mentioned in the table above. The diagnostic strips are intended for single use only.

SPECIMEN COLLECTION AND HANDLING
For specimen collection and preparation use only suitable clean and dry tubes or containers without traces of detergents and disinfectants.

Do not add the preservatives to the urine.

Use freshly voided urine, well mixed that has not been centrifuged. The urine should not be more than 2 hours old when tested. In case of longer standing, mix before use.

Do not expose the sample into the direct sunlight as this may cause the degradation of bilirubin and urobilinogen and this may turn in lower or false negative results. Mind that woman menstruation blood may affect the results.

ASSAY PROCEDURE:

- For insertion of strips into analyser follow the instruction of automatic urine analyser LAURA XL. Select MAINTENANCE → FEED STRIPS → chose the option of strips management (new strips vs. already used strips). Follow the pop-up window instructions.
- For more information about strips management refer to LAURA XL user manual.
- Log the tube into system through RFID tag placed at the bottom of the tube.
- Open the bottle and insert the whole volume of tube into feeder of automatic urine analyser LAURA XL. Do not divide test strips!
- Close the feeder and load the strips.
- After analysis by LAURA XL, the strip is evaluated automatically, and the result will be reported in LAURA XL user interface.

EVALUATION
After contact of sample with individual pads, diagnostics zones change their colour. The intensity and shade of zone coloration are proportional to the concentration of analyte. After analysis, the result will be automatically reported by LAURA XL analyser.

WORKING REAGENT CONCENTRATIONS

Specific gravity: poly(methylvinylether/maleic acid) 32 %; bromthymol blue 5.1 %

Leucocytes: indoxyl ester 0.43 %; diazonium salt 0.05 %

Nitrites: sulphaniamidé 5.1 %; tetrahydrobenzo-*h*-quinoline 5.8 %

pH: methyl red 0.71 %, bromthymol blue 12.1 %

Ascorbic Acid: phosphomolybdic acid 26 %

Protein: tetrabromphenolftalein ester 0.21 %; tetrabromphenol blue 0.35 %

PRINCIPLE
Specific gravity – The test is based on the principle of ion exchange, which runs between polyelectrolyte and ions present in the urine. Its result is a colour change of an acid-base indicator from the blue-green colour in the urine with low concentration of ions, through green and yellow-green in the urine with increased concentration of ions, to amber yellow colour.

Leucocytes – The test is based on the enzymatic reaction. The test pad contains an indoxyl ester, that is cleaved by the granulocyte esterases. The released indoxyl reacts with a diazonium salt and violet coloration is formed. The colour intensity is proportional to the leucocytes amount in a sample of tested urine.

Nitrites – The test is based on conversion of nitrate to nitrite by the action of certain species of bacteria contained in the urine. The colour test is based on the principle of the Griess test. Some degree of pink coloration should be interpreted as a positive nitrite test suggestive of 10³ or more organisms in 1 mL of the urine specimen.

pH – The test is based on the double indicator principle and gives a range of colours from orange through yellow and green to blue and permits differentiation to within 0.5 pH unit in the range pH 5 to 9.

Ascorbic Acid – The test is based on the reaction of phosphomolybdic acid which is reduced by ascorbic acid to molybdenum blue. The test is not specific for ascorbic acid because the green to greyish-blue colour of the test pad is exhibited also by other strongly reducing substances present in urine, such as gentisic acid and other acetylsalicylic acid metabolites. We recommend carrying out determination of ascorbic acid in urine especially in cases, in which ascorbic acid may disturb the tests for other urine constituents, such as glucose, blood and nitrite.

Protein – The test is based on colour change of acid-base indicator, which is caused by presence of proteins. It is particularly sensitive to albumin, but is much less sensitive to globulin, mucoprotein, haemoglobin and Bence-Jones protein.

Glucose – The test is based on the specific glucose oxidase/oxidase reaction and is specific for D-glucose. The reagent pad does not react with other sugars, it reacts with presence of D-glucose by green to the dark green coloration.

Ketones – The test is based on the principle of Legal's test and is more susceptible to the acetoacetic acid than to acetone. Test does not react with the β- hydroxybutyric acid. The colour scale is calibrated for the acetoacetic acid.

Urobilinogen – The test is based on the coupling of urobilinogen with a stabilized reagent. The test is specific for urobilinogen and sterobilinogen and is not susceptible to the interfering factors known at Ehrlich's test.

Bilirubin – The test is based on the coupling of bilirubin with stabilized reagent. The reaction is not affected by pH of the urine.

Blood – The test is based on the peroxidase activity of haemoglobin which catalyses the oxidation of the indicator due to the presence of the organic hydroperoxide contained in the diagnostic pad. The label contains two colour scales: for detection of intact erythrocytes and free haemoglobin. The test is highly sensitive to free haemoglobin and may detect its presence from concentrations corresponding approx. to 5 Ery/μL.

Compensation field – Pad, which isn't impregnated with any reagents. The compensation field is used for suppressing of the dark colour of the urine sample, because the dark colour could have the effect on the evaluation of reagent pads.

LIMITATIONS
Specific gravity – This reaction is affected by pH values of urine over 6.5, this can cause shift in colour response towards lower values of specific gravity. This phenomenon could be automatically corrected by LAURA XL (SG compensation).

Leucocytes – In case when the urine sample is more markedly coloured (e.g. increased bilirubin content), the resulting colour could be affected by a sample coloration. The intensity of the colour reaction is increased by alkaline pH and higher urine density.

Nitrites – Before testing the patient should intake vegetable-rich meals and discontinue antibiotic therapy for 3 days prior to the test. Sensitivity of this test declines with the high specific gravity of the urine. Increased diuresis can cause the false negative results. Limited fluid intake prior testing can prevent from the excessive dilution of the urine. The test can be applied only at the fresh urine. Inaccurate results may occur at the stale urine, in which nitrite can be formed by contamination of the specimen.

Protein – In strongly alkaline urines (pH >8) from patients on medication with guanine or quinoline containing drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quaternary ammonium groups. On the other hand, in the presence of non-ionic or anionic detergents, false-negative results may occur. Do not take the colour of the dry pad into consideration.

Glucose – The reaction is independent on pH and presence of ketone bodies.

Ketones – Drugs and diagnostics based on phenolphthalein or sulphophthalein may turn red to purple because of alkaline reaction of the pad.

Urobilinogen – The reaction is not affected by pH of urine. The presence of bilirubin gives yellow colour. This colour which turns slowly to greenish-blue does not interfere with urobilinogen determination provided the reading is made 1 minute after wetting. The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of urobilinogen and thus leads to artificially low or false negative readings.

Bilirubin – The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of bilirubin and thus leads to artificially low or false negative readings. High concentrations of urobilinogen (above 100 μmol/L) interfere with the test. Also the red substances or the substances, which are turning red in low pH may interfere (e.g. Phenazopyridine).

Blood – Microbial peroxidase associated with urinary tract infection may cause a false positive reaction. The sensitivity of the test is influenced by specific gravity or by inhibitors of pharmacological origin.

PLEASE NOTE

All diagnostic pads are protected against the common concentration of the ascorbic acid. Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after discontinuing a drug. The sensitivity depends upon the variability of urines. The semiquantitative analysis is not sufficient for the completing of diagnoses.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Control urines: URINORM XL, Cat. No. REG00060, or Bio-Rad Liquichek Urinalysis Control

Automatic urine analyser: LAURA XL, Cat. No. INS00065

QUALITY CONTROL

For quality control of precision and accuracy URINORM XL control urines (Cat. No. REG00060) are recommended. Obtained results are compared with the declared values in the instruction for use of URINORM XL. Bio-Rad Liquichek Urinalysis Control can be also used as third-party control material. It is recommended to perform QC measurements according to local laboratory guidelines. More information about quality control can be found in LAURA XL user manual.

STABILITY AND STORAGE

Keep diagnostic test strips in tightly closed original tubes in a dry and dark place at 2–30 °C. The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory. When stored under these conditions, test strips are stable to the expiry date given on the pack. The on-board stability is 5 days in feeder of LAURA XL automatic urine analyser.

ANALYTICAL PERFORMANCE

Limit of detection is defined as the lowest concentration of the analyte which leads to positive results in ≥ 90% cases. Limit of detection is not applicable for SG and pH.

The range of assay is defined as an interval between Limit of detection and the last positive level for each analyte.

Parameter	Limit of detection	Range of Assay
SG	1.003 kg/m ³	1.000–1.030 kg/m ³
LEU	15 Leu/μL	Neg.–ca. 500 Leu/μL
NIT	0.90 mg/L	Neg.–Pos.
pH	5.5*	5–9
ASC	0.60 mmol/L	Neg.–3.4 mmol/L
PRO	0.30 g/L	Neg.–5 g/L
GLU	1.70 mmol/L	Neg.–55 mmol/L
KET	0.50 mol/L	Neg.–15 mmol/L
UBG	17 μmol/L	Neg.–203 μmol/L
BIL	15 μmol/L	Neg.–102 μmol/L
BLD	4 Ery/μL	Neg.–ca. 250 Ery/μL

* The Limit of Detection for SG and pH zones is not applicable therefore the first cut-off value is presented.

EXPECTED VALUES

Based on literature (whole population),^[1]

Blood: up to 5 Ery/μL
Glucose: up to 1.4 mmol/L
Nitrite: up to 10 Leu/μL
10³ bacteria/mL
Leucocytes: up to 1 mg/dL
Bilirubin: up to 17 μmol/L
Urobilinogen: up to 1 mg/dL
Protein: over 30 mg/24 hours
Ketones: up to 0.19 mmol/L
Specific gravity: 1.015–1.025
Ascorbic acid: not stated

Each laboratory should verify that it is appropriate to apply these values to its patient population and if necessary, determine its own ranges.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Do not touch test pads of the strip.

The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The product is not classified as dangerous.

WASTE MANAGEMENT

Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

REFERENCES

- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.

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