BILIRUBIN DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00009	BIL D 200	R1: 4 × 40 mL, R2: 1 × 40 mL instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of bilirubin direct in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases. Important for monitoring in neonates. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Red blood cells at the end of their circulating lives are broken down in the reticuloendothelial system, mainly the spleen. The resulting heme is converted to bilirubin upon removal of iron. This process accounts for about 80 % of the 500 µmol (292 mg) of bilirubin formed daily. Other sources of bilirubin include the breakdown of myoglobin and cytochromes and the catabolism of immature red blood cells in the bone marrow.


Once formed, bilirubin is transported to the liver bound to albumin. This fraction of bilirubin is referred to as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid (mono- and diglucuronides) by the enzyme uridyl diphosphate glucuronyl transferase to form conjugated bilirubin. Conjugated bilirubin or direct bilirubin is excreted via the biliary system into the intestine, where it is metabolized by bacteria to a group of products known collectively as sterco-bilinogen. Elimination is almost complete and serum levels are normally negligible.

Total bilirubin is the sum of the unconjugated and conjugated fractions. Total bilirubin is elevated in hepatitis, cirrhosis, hemolytic disorders, several inherited enzyme deficiencies, and conditions causing hepatic obstruction.

TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

PRINCIPLE

Bilirubin reacts with diazotized sulphanilic acid in a strongly acidic medium and form an intensely coloured diazo dye – azobilirubin^{1,2,3,4,5}.

Absorbance of the Azobilirubin measured at 546 nm is proportional to the concentration of the direct bilirubin in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	R2
Sulphanilic acid 28.9 mmol/L	Sodium nitrite 2.9 mmol/L
HCl 23.0 mmol/L	

COMPOSITION OF REACTION MIXTURE

Sulphanilic acid	22.0 mmol/L
HCl	17.5 mmol/L
Sodium nitrite	0.55 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Prepare working reagent by mixing of 4 portion of reagent R1 with 1 portion of reagent R2.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 546 nm may be used, general laboratory equipment.

XL MULTICAL 4x3, Cat. No. XSYS0034

XL MULTICAL 10x3, Cat. No. XSYS0122

ERBA NORM 4x5, Cat. No. BLT00080

ERBA NORM 10x5, Cat. No. XSYS0123

ERBA PATH 4x5, Cat. No. BLT00081

ERBA PATH 10x5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

The working reagent is stable for 7 days at 2–8 °C, when protected from contamination and light. It is recommended to prepare fresh working solution before assay is performed.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma^{6,7}: 2 days at 15–25 °C
7 days at 4–8 °C
6 months at -20 °C

Protect sample from light.

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2-point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

• after reagent lot change

• as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the Doumas method¹¹.

ASSAY PROCEDURE

Wavelength: 546 / 670 nm

Cuvette: 1 cm

	Reagent blank	Calibrator	Sample
Working reagent	1.00 mL	1.00 mL	1.00 mL
Sample	–	–	0.05 mL
Calibrator	–	0.05 mL	–
Distilled water	0.05 mL	–	–

Mix and incubate 5 min. at 37 °C. Measure absorbance of the sample A_{sam} and calibrator A_{cal} against reagent blank.

CALCULATION

$$1. \text{ Direct Bilirubin (mg/dL)} = \frac{A_{\text{sam}}}{A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

2. Using factor (f):

$$\text{Direct Bilirubin (mg/dL)} = A_{\text{sam}} \times f \quad f = 16$$

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End point	Normal Low (mg/dL)	0
Wavelength 1 (nm)	546	Normal High (mg/dL)	0.2
Wavelength 2 (nm)	670	Linearity Low (mg/dL)	0.04
Sample Volume (µL)	25/50	Linearity High (mg/dL)	20.5
Working Reagent Volume (µL)	500/1000	Factor	16
Incubation time (min.)	5	Blank with	Reagent
Reaction temperature (°C)	37	Absorbance limit (max.)	0.15
Reaction direction	Increasing	Units	mg/dL

UNIT CONVERSION

mg/dL × 17.1 = µmol/L

EXPECTED VALUES¹²

In serum:

Adults and infants: 0–0.2 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.04 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

Linearity: 20.5 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	0.59	0.005	0.80
Sample 2	1.56	0.013	0.82

Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	0.55	0.012	2.27
Sample 2	1.62	0.033	2.01

Accuracy

Two different validated control materials were used. Determined bias is and 6.1 % at the target value 0.35 mg/dL and 7.9 % at the target value 0.99 mg/dL.

Comparison

A comparison between XL-640 automatic system BILIRUBIN DIRECT (y) and a commercially available test (x) using 120 samples gave following results:

Linear regression:

$$y = 0.922x + 0.154 \text{ mg/dL} \quad r = 0.992$$

Passing-Bablok¹³:

$$y = 0.993x + 0.134 \text{ mg/dL} \quad r = 0.882$$

Interferences

Criterion: Recovery within ±10 % of initial value of direct bilirubin concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 3 g/L, triglycerides up to 850 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels, except positive interferences of Rifampicin¹⁴.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 546 nm against the distilled water is 0.15.

- High concentration of haemoglobin and triglycerides in sample can interfere with determination of direct bilirubin. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety and Performance is available at Eudamed see:

<https://ec.europa.eu/tools/eudamed>.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: U61F-WEUA-9H7R-X1MU

Danger

Hazard statements:

H314 Causes severe skin burns and eye damage.

Precautionary statements:

P260 Do not breathe vapours.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Supplemental information

EUH 208 Contains sulphanilic acid. May produce an allergic reaction.

R2

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Pearlman FC and Lee RTY. Detection and measurement of total bilirubin in serum, with use of surfactants as solubilizing agents. Clin Chem 20: 447-453, 1974.
2. Zoppi F, Peracino A, Fenili D, Marcovina S and Ramella C. Metodo per la determinazione della bilirubina totale e coniugata. Uso di un tensioattivo cationico come agente solubilizzante. Giorn It Chim Cl 1: 343-359, 1976.
3. Malloy HT, Evelyn KA: J. Biol. Chem. 119: 481-490, 1937.
4. Walters MI and Gerarde HW, An ultramicromethod for the determination of conjugated and total bilirubin in serum or plasma, Microchemical Journal. 15: 231-243, 1970.
5. Henry, RJ, Editor, Clinical Chemistry, Principles and Technics, p.1058, Harper and Row, Publishers, Hagerstown, Maryland, 1974.
6. Quality of Diagnostic Samples, Recommendations of the Working Group on Preamalytical Quality of the German Society for Clinical Chemistry and Laboratory medicine, 3rd completely revised ed. 2010.
7. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
8. Balistreri WF, Shaw LM. Liver function. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders: 729-761, 1987.
9. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 192-202, 1998.
10. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1125-77, 1999.
11. Dumas BT, Perry BW, Jendrzczak B, et al. Pitfalls in the American Monitor Kit Methods for Determination of Total and Direct Bilirubin. Clin Chem 28 (11), 2305-2308, 1982.
12. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
13. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
14. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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