GLUCOSE

Cat. No.	Pack Name	Packaging (Content)
BLT00025	GLU 500	R1: 2 × 250 mL, R2 standard: 2 × 5 mL instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of glucose in human serum, plasma and urine on various automatic systems. Intended for screening, monitoring and diagnosis of diabetes, hypoglycemia. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Glucose is the major carbohydrate present in the peripheral blood. Oxidation of glucose is the major source of cellular energy in the body. Glucose derived from dietary sources is converted to glycogen for storage in the liver or to fatty acids for storage in adipose tissue. The concentration of glucose in blood is controlled within narrow limits by many hormones, the most important of which are produced by the pancreas. Accurate measurement of glucose in body fluid is important in diagnosis and management of diabetes, hypoglycemia, adrenal dysfunction and various other conditions. High levels of serum glucose may be seen in case of Diabetes mellitus, in patients receiving glucose containing fluids intravenously, during severe stress and in cerebrovascular accidents. Decreased levels of glucose can be due to insulin administration, as a result of insulinoma, inborn errors of carbohydrate metabolism or fasting. Glucose measurement in urine is used as a diabetes screening procedure and to aid in the evaluation of glycosuria, to detect renal tubular defects, and in the management of diabetes mellitus.

PRINCIPLE

Glucose in the sample is oxidised to yield gluconic acid and hydrogen peroxide in the presence of glucose oxidase. The enzyme peroxidase catalyses the oxidative coupling of 4-aminoantipyrine (4AAP) with phenol to yield a coloured quinoneimine complex (Tinder's reaction¹):

Absorbance of coloured quinoneimine complex measured at 490–550 nm is proportional to the concentration of glucose in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2 standard	
Phosphate buffer (pH 7.15)	250 mmol/L	Glucose	See bottle label
Glucose oxidase	>25 KU/L		
Peroxidase	>2 KU/L		
Phenol	5 mmol/L		
4-aminoantipyrine	0.5 mmol/L		
Sodium azide	0.5 g/L		

COMPOSITION OF REACTION MIXTURE

Phosphate buffer (pH 7.15)	248 mmol/L
Glucose oxidase	>25 KU/L
Peroxidase	>2 KU/L
Phenol	5 mmol/L
4-aminoantipyrine	0.5 mmol/L
Sodium azide	0.5 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 490–550 nm may be used, general laboratory equipment.
 XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagent is stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum:
 Plasma: Li-heparin and EDTA plasma.
 Separate at the latest 1h after blood collection from cellular contents.

Urine:
 The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.
 See the limitations and interferences section for details about possible sample interferences.

Stability in plasma after addition of a glycolytic inhibitor (Fluoride, monoiodacetate, mannose) ² :		
	2 days at	20–25 °C
	7 days at	4–8 °C
	1 day at	-20 °C

Discard contaminated specimens.

Stability in serum (separated from cellular contents, hemolysis free) without adding a glycolytic inhibitor ^{3,4} :		
	8 hours at	25 °C
	72 hours at	4 °C

Stability in urine ⁵ :		
	24 hours at	4–8 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank
 Calibration frequency: it is recommended to do a calibration
 • after reagent lot change
 • as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.
 The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL, R2 standard and controls ERBA NORM and ERBA PATH have been standardized against ID/MS.

ASSAY PROCEDURE

Wavelength: 500 (490–550) nm
 Cuvette: 1 cm

	Reagent blank	Calibrator (Standard)	Sample
Reagent 1	1000 mL	1000 mL	1000 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and incubate 5–10 min. at 37 °C. Measure absorbance of the sample A_{sam} and calibrator (standard) A_{cal} against reagent blank.

CALCULATION

$$\text{Glucose (mg/dL)} = \frac{A_{\text{sam}}}{A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator (standard) concentration}$$

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point	Normal High (mg/dL)	100
Wavelength 1 (nm)	505	Linearity Low (mg/dL)	0.96
Wavelength 2 (nm)	670	Linearity High (mg/dL)	550
Sample Volume (µL)	5/10	Concentration of Standard	See bottle label
Working Reagent Volume (µL)	500/1000	Blank with	Reagent
Incubation time (min.)	5	Absorbance limit (max.)	0.3
Reaction temperature (°C)	37	Units	mg/dL
Reaction direction	Increasing		
Normal Low (mg/dL)	74		

UNIT CONVERSION

mg/dL × 0.056 = mmol/L

EXPECTED VALUES^{5,6}

At 37 °C

Glucose Fasting:

Cord:	45–96 mg/dL
Newborn, 1 d:	40–60 mg/dL
Newborn, >1 d:	50–80 mg/dL
Child:	60–100 mg/dL
Adult:	74–100 mg/dL
>60 y:	82–115 mg/dL
>90 y:	75–121 mg/dL

Random Glucose:

≥1 y:	70–140 mg/dL
Glucose 2 h Postprandial:	<120 mg/dL
WB (Hep) Adult:	65–95 mg/dL
Urine:	1–15 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

0.96 mg/dL
 Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV < 20% (n = 30).

Linearity:

550 mg/dL
 Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	81.9	1.16	1.41
Sample 2	245.8	1.30	0.53

Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	108.6	1.36	1.26
Sample 2	340.9	1.02	0.30

Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	81.7	1.86	2.28
Sample 2	246.1	3.94	1.60

Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	109.0	2.19	2.01
Sample 2	321.3	6.45	2.01

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is 2.6 % at the target value 69.5 mg/dL, 0.9 % at the target value 211.4 mg/dL for serum, 2.1 % at the target value 27.7 mg/dL and 4.1 % at the target value 294.6 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system GLUCOSE (y) and a commercially available test (x) using 146 samples (serum) gave following results:
 Linear regression:
 $y = 1.027x - 7.499$ mg/dL $r = 0.982$
 Passing-Bablok⁷:
 $y = 1.025x - 7.560$ mg/dL $r = 0.963$

Interferences

Criterion: Recovery within ±10 % of initial value of glucose concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs:
 N-acetylcysteine, Metamizole and Acetaminofen (including metabolite N-acetyl-p-benzoquinone imine) may decrease the measured results in serum⁸.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 505 nm against the distilled water is 0.3.
 - High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of glucose. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
 The Summary of Safety and Performance is available at Eudamed see: <https://ec.europa.eu/tools/eudamed>.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1
 UFI: WY5U-AJ46-VE76-0WED

Warning

Hazard statement:
 H315 Causes skin irritation.
 H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.
 P302 + P352 IF ON SKIN: Wash with plenty of water.
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

R2

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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


1. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. *Analyst* 1972; 97: 142-5.
2. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; 30-31.
3. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. 750-808.
4. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002; 48: 436-72.
5. Tietz NW, (Ed.), *Textbook of Clinical Chemistry*. Burtis CA and Ashwood ER, Fifth Edition, 2012.
6. *Tietz Textbook of Clinical Chemistry*. 5th ed. Edited by Burtis CA, Ashwood ER, Bruns DE. WB Saunders Company, Philadelphia 2006, 25: 837-907.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11): 783-790.
8. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, *Clin Biochem* 49, 100-104, 2016.

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LITERATÚRA

1. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. *Analyst* 1972; 97: 142-5.
2. Guder WG, Zawta B et al. *The Quality of Diagnostic Samples*. 1st ed. Darmstadt: GIT Verlag; 2001; 30-31.
3. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. 750-808.
4. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. *Clin Chem* 2002; 48: 436-72.
5. Tietz NW, (Ed.), *Textbook of Clinical Chemistry*. Burtis CA and Ashwood ER, Fifth Edition, 2012.
6. *Tietz Textbook of Clinical Chemistry*. 5th ed. Edited by Burtis CA, Ashwood ER, Bruns DE. WB Saunders Company, Philadelphia 2006, 25: 837-907.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11): 783-790.
8. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, *Clin Biochem* 49, 100-104, 2016.

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