LDL DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00041	LDL 80	R1: 2 × 30 mL, R2: 2 × 10 mL, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of LDL cholesterol in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of coronary heart disease and lipoprotein metabolism disorders. For professional use in clinical laboratory only.


CLINICAL SIGNIFICANCE

Low Density Lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on triglyceride-rich Very Low Density Lipoproteins (VLDLs). Specific LDL receptors exist to facilitate the elimination of LDL from plasma by liver parenchymal cells. It has been shown that most of the cholesterol stored in atherosclerotic plaques originates from LDL. For this reason the LDL Cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.

Accurate measurement of LDL Cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

PRINCIPLE

The assay is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol-methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL, and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHO) and cholesterol esterase (CHE), whereas HDL reacts with the enzymes. Addition of R2 containing a specific detergent releases LDL from the PVS/PEGME complex. The released LDL reacts with the enzymes to produce hydrogen peroxide (H₂O₂) which is quantified by the Trinder's reaction^{1,2,3}.

The absorbance of the produced quinone dye at 600 nm is proportional to the LDL concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
MES buffer (pH 6.5)	50 mmol/L
Polyvinyl sulfonic acid	50 mg/L
Polyethylene-glycol-methyl ether	30 ml/L
4-aminoantipyrine, 4-AAP	0.9 g/L
Cholesterol esterase	5 kU/L
Cholesterol oxidase	20 kU/L
Peroxidase	5 kU/L
Detergent	

R2	
MES buffer (pH 6.5)	50 mmol/L
N-Bis(4-sulfobutyl)-3-methylaniline, TODB	3 mmol/L
Detergent	

COMPOSITION OF REACTION MIXTURE

MES buffer (pH 6.5)	49.6 mmol/L
Polyvinyl sulfonic acid	37 mg/L
Polyethylene-glycol-methyl ether	22 ml/L
4-aminoantipyrine, 4-AAP	0.7 g/L
N-Bis(4-sulfobutyl)-3-methylaniline, TODB	0.7 mmol/L
Cholesterol esterase	3.7 kU/L
Cholesterol oxidase	14.9 kU/L
Peroxidase	3.7 kU/L
Detergent	

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance at 600/700 nm may be used, general laboratory equipment.

HDL/LDL Calibrator, Cat. No. XSYS0061
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma.

Fasting and non-fasting samples can be used¹.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ⁴ :	1 day at	20–25 °C
	7 days at	4–8 °C
	3 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with HDL/LDL Calibrator is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, HDL/LDL calibrator and controls ERBA NORM and ERBA PATH have been standardized against the reference material NIST SRM 1951.

ASSAY PROCEDURE

Wavelength: 600/700 nm

Cuvette: 1 cm

	Reagent blank	Calibrator	Sample
Reagent 1	0.375 mL	0.375 mL	0.375 mL
Sample	–	–	0.003 mL
Calibrator	–	0.003 mL	–
Distilled water	0.003 mL	–	–

Mix and after 5 min. incubation read the initial absorbance for blank A_{bl} , sample A_{sam} and calibrator A_{cal} . Then add:

Reagent 2	0.125 mL	0.125 mL	0.125 mL
-----------	----------	----------	----------

Mix and after 5 min. incubation read the final absorbance for blank A_{bl} , sample A_{sam} and calibrator A_{cal} . Calculate resulting absorbance like the difference between the final and initial absorbance $A = (A_{\text{FINAL}} - A_{\text{INITIAL}})$.

CALCULATION

$$\text{LDL cholesterol (mg/dL)} = \frac{A_{\text{sam}} - A_{\text{bl}}}{A_{\text{cal}} - A_{\text{bl}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	1-Point End	Normal Low (mg/dL)	0
Wavelength 1 (nm)	600	Normal High (mg/dL)	130
Wavelength 2 (nm)	670	Linearity Low (mg/dL)	0.5
Sample Volume (µL)	3	Linearity High (mg/dL)	413
Reagent 1 Volume (µL)	375	Concentration of Standard	See bottle label
Reagent 2 Volume (µL)	125	Blank with	Reagent
Incubation time (min.)	5	Absorbance limit (max.)	0.3
Reaction temperature (°C)	37	Units	mg/dL
Reaction direction	Increasing		

UNIT CONVERSION

mg/dL × 0.026 = mmol/L

EXPECTED VALUES¹⁶

In serum:	Male	Female
5–9 y	63–129	68–140 mg/dL
10–14 y	64–133	68–136 mg/dL
15–19 y	62–130	59–137 mg/dL
20–24 y	66–147	57–159 mg/dL
25–29 y	70–165	71–164 mg/dL
30–34 y	78–185	70–156 mg/dL
35–39 y	81–189	75–172 mg/dL
40–44 y	87–186	74–174 mg/dL
45–49 y	97–202	79–186 mg/dL
50–54 y	89–197	88–201 mg/dL
55–59 y	88–203	89–210 mg/dL
60–64 y	83–210	100–224 mg/dL
65–69 y	98–210	92–221 mg/dL
>69 y	88–186	96–206 mg/dL

Coronary heart disease risk:

Optimal	<100 mg/dL
Near/above optimal	100–129 mg/dL
Borderline high	130–159 mg/dL
High	160–189 mg/dL
Very high	>189 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.46 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV < 20 % (n = 30).

Linearity: 413 mg/dL

Linearity is the highest measured activity with recovery within ± 10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	73.7	0.59	0.80	Sample 1	73.3	1.59	2.17
Sample 2	132.2	0.72	0.54	Sample 2	136.9	3.34	2.44

Accuracy

Two different validated control materials were used. Determined bias is 9.9 % at the target value 128.8 mg/dL and 5.5 % at the target value 188.5 mg/dL.

Comparison

A comparison between XL-640 automatic system LDL DIRECT (y) and a commercially available test (x) using 140 samples gave following results:

Linear regression:

$$y = 0.920x + 6.591 \text{ mg/dL} \quad r = 0.986$$

Passing-Bablok¹⁷:

$$y = 0.922x + 6.161 \text{ mg/dL} \quad r = 0.988$$

Interferences

Criterion: Recovery within ± 10 % of initial value of LDL cholesterol concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 650 mg/dL.

N-acetylcysteine, Metamizole and Acetaminofen (Paracetamol) including its metabolite N-acetyl-p-benzoquinone imine may cause false-negative results^{18,19}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 600 nm against the distilled water is 0.3.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of LDL cholesterol. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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

1. Pisani T, Gebiski CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Barham D, Trinder P, An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
3. Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl Lab 31: 325-330, 1985.
4. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
5. Rifai N, Warnick GR, McNamara JR, et al. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 38: 150-160, 1992.
6. Naito HK, Strong JP, Scott MG, et al. Atherogenesis: current topics on etiology and risk factors. Clin Chem 41: 132-133, 1995.
7. "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285: 2486, 2001.
8. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26: 566, 1985.
9. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
10. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
11. Badimon JJ, Badimon L, Fuester V, Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85: 1234-41, 1990.
12. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
13. Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
14. Williams P, et al, High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
15. Friedwald, W. T., Levy, R. I. and Frederickson, D. S.: Estimation of the concentration of low density lipoprotein cholesterol in plasma without use of ultracentrifuge. Clin. Chem. 18, 449-502, 1972.
16. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
17. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
18. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
19. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100-104, 2016.

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LITERATÚRA

1. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Barham D, Trinder P, An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
3. Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl Lab 31: 325-330, 1985.
4. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
5. Rifai N, Warnick GR, McNamara JR, et al. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 38: 150-160, 1992.
6. Naito HK, Strong JP, Scott MG, et al. Atherogenesis: current topics on etiology and risk factors. Clin Chem 41: 132-133, 1995.
7. "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285: 2486, 2001.
8. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26: 566, 1985.
9. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
10. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
11. Badimon JJ, Badimon L, Fuester V, Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85: 1234-41, 1990.
12. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
13. Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
14. Williams P, et al, High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
15. Friedwald, W. T., Levy, R. I. and Frederickson, D. S.: Estimation of the concentration of low density lipoprotein cholesterol in plasma without use of ultracentrifuge. Clin. Chem. 18, 449-502, 1972.
16. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
17. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
18. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
19. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100-104, 2016.

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REF

Katalógové číslo

LOT

Číslo šarže

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eFU:
www.erba.com**IVD**Diagnostický zdravotnícky prostriedok *in vitro*

Výrobca

Obmedzenie teploty

CONT

Obsah

QUALITY SYSTEM CERTIFIED
ISO 13485, IVDRErba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com

CC/IFU/075/26/A

Dátum revízie: 27. 4. 2026