# UREA

Cat. No.	Pack Name	Packaging (Content)
BLT00060	UREA 1000	R1: 4 × 200 mL, R2: 1 × 200 mL instruction for use
BLT00061	UREA 250	R1: 4 × 50 mL, R2: 1 × 50 mL, R3: 1 × 5 mL instruction for use

## INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of urea in human serum, plasma and urine on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases, kidney function. For professional use in clinical laboratories only.

## CLINICAL SIGNIFICANCE

Urea is the major end product of protein nitrogen metabolism. It is synthesized by the urea cycle in the liver from ammonia which is produced by amino acid deamination. Urea is excreted mostly by the kidneys but minimal amounts are also excreted in sweat and degraded in the intestines by bacterial action. Determination of blood urea nitrogen (BUN) is the most widely used screening test for renal function. When used in conjunction with serum creatinine determinations it can aid in the differential diagnosis of the three types of azotemia: prerenal, renal and postrenal. Elevations in blood urea nitrogen concentration are seen in inadequate renal perfusion, shock, diminished blood volume (prerenal causes), chronic nephritis, nephrosclerosis, tubular necrosis, glomerular nephritis (renal causes) and urinary tract obstruction (postrenal causes). Transient elevations may also be seen during periods of high protein intake. Unpredictable levels occur with liver diseases.

## PRINCIPLE

The enzyme methodology employed in this reagent is based on the reaction first described by Talke and Schubert. To shorten and simplify the assay, the calculations are based on the discovery of Tiffany et al. that urea concentration is proportional to absorbance change over a fixed time interval<sup>1,2,3,4</sup>. Urea is hydrolyzed by urease to form ammonium and carbonate. In the second reaction 2-oxoglutarate reacts with ammonium in the presence of glutamate dehydrogenase (GLDH) and the coenzyme NADH to produce L-glutamate.

The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due the oxidation of NADH. The rate of oxidation of NADH is proportional to the urea concentration in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Tris buffer	100 mmol/L	NADH	1.66 mmol/L
2-oxoglutarate	5.49 mmol/L	Sodium azide	0.9 g/L
Urease (Jack Bean)	≥10 kU/L		
GLDH (Microorganism)	≥3.8 kU/L		
Sodium azide	0.9 g/L		

**R3 standard**  
Urea See bottle label

## COMPOSITION OF REACTION MIXTURE

Tris buffer	79 mmol/L
2-oxoglutarate	4.35 mmol/L
Urease (Jack Bean)	≥7.9 kU/L
GLDH (Microorganism)	≥3.0 kU/L
NADH	0.33 mmol/L
Sodium azide	0.89 g/L

## REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 340 nm may be used, general laboratory equipment.  
XL MULTICAL 4×3, Cat. No. XSYS0034  
XL MULTICAL 10×3, Cat. No. XSYS0122  
ERBA NORM 4×5, Cat. No. BLT00080  
ERBA NORM 10×5, Cat. No. XSYS0123  
ERBA PATH 4×5, Cat. No. BLT00081  
ERBA PATH 10×5, Cat. No. XSYS0124

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

### Two reagents method – substrate start

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

### Monoreagent method

Stability of working reagent: 5 days at 15–25 °C in dark  
4 weeks at 2–8 °C in dark

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.  
For specimen collection and preparation only use suitable tubes or collection containers.  
Only the specimens listed below were tested and found acceptable.

### Serum.

Plasma: Li-heparin and K<sub>2</sub>-EDTA plasma. Do not use ammonium heparin.

Urine: Dilute urine samples in 1+99 ratio with distilled water and multiply results by 100. Bacterial growth in the specimen and high atmospheric ammonia concentrations as well as contamination by ammonium ions may cause erroneously elevated results.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma <sup>6</sup> :		
7 days at	15–25 °C	
7 days at	2–8 °C	
1 year at	-20 °C	

Stability in urine <sup>6</sup> :		
2 days at	15–25 °C	
7 days at	2–8 °C	
1 month at	-20 °C	

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.  
2 point calibration (blank and calibrator); distilled water is recommended as blank  
Calibration frequency: it is recommended to do a calibration  
• after reagent lot change  
• as required by internal quality control procedures

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, calibrator XL MULTICAL, R3 standard and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS.

## ASSAY PROCEDURE

Wavelength: 340 nm

Cuvette: 1 cm

## Two reagents method

	Reagent blank	Calibrator (Standard)	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and after 1 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix and measure the initial absorbance after 30 sec (A<sub>1</sub>), start timer simultaneously and read again exactly after 1 min (A<sub>2</sub>). Measure against reagent blank. Calculate absorbance change ΔA = (A<sub>2</sub>–A<sub>1</sub>)/min.

## Monoreagent method

	Reagent blank	Calibrator (Standard)	Sample
Working reagent	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and measure the initial absorbance after 30 sec (A<sub>1</sub>), start timer simultaneously and read again exactly after 1 min (A<sub>2</sub>). Measure against reagent blank. Calculate absorbance change ΔA = (A<sub>2</sub>–A<sub>1</sub>)/min.

## CALCULATION

$$\text{Urea (mg/dL)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator (standard) concentration}$$

## ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic	Normal Low mg/dL	17
Wavelength (nm)	340	Normal High mg/dL	43
Sample Volume (μL)	50/100	Linearity Low mg/dL	1.28
Working Reagent Volume (μL)	500/1000	Linearity High mg/dL	219
Lag time (sec.)	20	Blank with	Reagent
Kinetic interval (sec.)	60	Absorbance limit (min.)	1.1
No. of readings	1	Concentration of Standard	See bottle label
Reaction temperature (°C)	37	Units	mg/dL
Reaction direction	Decreasing		

## UNIT CONVERSION

mg/dL × 0.1665 = mmol/L  
Urea (mg/dL) × 0.467 = BUN (mg/dL)  
BUN (mg/dL) × 2.14 = Urea (mg/dL)

## EXPECTED VALUES

### In Serum / Plasma<sup>7</sup>

Adults	17–43 mg/dL	Children	1–3 years	11–36 mg/dL
Global		4–13 years	15–40 mg/dL	15–36 mg/dL
Women <50 years	15–40 mg/dL	14–19 years	21–43 mg/dL	18–45 mg/dL
Women >50 years	21–43 mg/dL			
Men <50 years	19–44 mg/dL			
Men >50 years	18–55 mg/dL			

Urea / Creatinine ratio<sup>7</sup>: 20–35 (mg/dL)/(mg/dL)

Urea in Urine<sup>8</sup>: 26–43 g/24 h

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

### Limit of quantification:

Serum / plasma: 1.28 mg/dL  
Urine: 124 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

### Linearity:

Serum / plasma: 219 mg/dL  
Urine: 21900 mg/dL  
Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

### Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	36.6	0.36	0.98	Sample 1	32.8	1.07	3.27
Sample 2	93.9	0.75	0.79	Sample 2	87.1	2.27	2.61

Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	906	19.6	2.16	Sample 1	1080	38.7	3.58
Sample 2	1656	46.8	2.83	Sample 2	2170	89.5	4.12

## Accuracy

Two different validated control materials for serum and urine were used. Determined bias is -7.0 % at the target value 124.6 mg/dL, -8.1 % at the target value 192.0 mg/dL for serum, 9.2 % at the target value 835 mg/dL and 9.0 % at the target value 1772 mg/dL for urine.

## Comparison

A comparison between XL-640 automatic system UREA (y) and a commercially available test (x) using 147 samples (serum) gave following results:

Linear regression:

$$y = 1.001x + 5.611 \text{ mg/dL} \quad r = 0.996$$

Passing-Bablok<sup>9</sup>:

$$y = 1.040x + 4.166 \text{ mg/dL} \quad r = 0.994$$

## Interferences

Criterion: Recovery within ±10 % of initial value of urea concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs: No interference in serum or urine was found at therapeutic concentrations using common drug panels<sup>10</sup>.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Minimal allowable absorbance of the reagent blank measured at 340 nm against the distilled water is 1.1.  
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of urea. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### R1, R2, R3

Reagents are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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


1. Richterich R, Colombo JP. *Klinische Chemie*. 4th ed. Basel: Karger S 1978:319-324.
2. Talke H, Schubert GA. Enzymatische Harnstoffbestimmung in Blut und Serum im optischen Test nach Warburg. *Klin Wochenschr* 1965;43:174.
3. Tiffany TO, Jansen JM, Burtis CA, et al. Enzymatic kinetic rate and end-point analyses of substrate, by use of a GeMSAEC Fast Analyzer. *Clin Chem* 1972;18:829-840.
4. Sampson EJ, Baired MA, Burtis CA, et al. A coupled-enzyme equilibrium method for measuring urea in serum: Optimization and evaluation of the AACC study group on urea candidate reference method. *Clin Chem* 1980;26:816-826.
5. Rock RC, Walker WG, Jennings CD. Nitrogen metabolites and renal function. In: Tietz NW, ed. *Fundamentals of Clinical Chemistry*. 3rd ed. Philadelphia: WB Saunders 1987;669-704.
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
7. Thomas L. *Clinical Laboratory Diagnostics*. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 374-7.
8. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1838.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov; 26(11): 783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.

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## LITERATÚRA

1. Richterich R, Colombo JP. Klinische Chemie. 4th ed. Basel: Karger S 1978:319-324.
2. Talke H, Schubert GA. Enzymatische Harnstoffbestimmung in Blut und Serum im optischen Test nach Warburg. Klin Wochenschr 1965;43:174.
3. Tiffany TO, Jansen JM, Burtis CA, et al. Enzymatic kinetic rate and end-point analyses of substrate, by use of a GeMSAEC Fast Analyzer. Clin Chem 1972;18:829-840.
4. Sampson EJ, Baird MA, Burtis CA, et al. A coupled-enzyme equilibrium method for measuring urea in serum: Optimization and evaluation of the AACC study group on urea candidate reference method. Clin Chem 1980;26:816-826.
5. Rock RC, Walker WG, Jennings CD. Nitrogen metabolites and renal function. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders 1987;669-704.
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
7. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 374-7.
8. Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1838.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov; 26(11): 783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.

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**REF**

Katalógové číslo

**LOT**

Číslo šarže

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eFU:  
[www.erba.com](http://www.erba.com)**IVD**Diagnostický zdravotnícky prostriedok *in vitro*

Výrobca

Obmedzenie teploty

**CONT**

Obsah

QUALITY SYSTEM CERTIFIED  
ISO 13485, IVDRErba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: [diagnostics@erba.com](mailto:diagnostics@erba.com), [www.erba.com](http://www.erba.com)

CC/IFU/025/26/A

Dátum revízie: 25. 5. 2026