CRP CONTROL

Cat. No.	Pack Name	Packaging (Content)
BLT20014	CRP CON L	1 × 1 mL, instruction for use

INTENDED USE

CRP CONTROL is control material for *in vitro* quantitative determination of C-reactive protein in human serum on automatic systems ERBA XL by turbidimetric immunoassay. Intended for monitoring accuracy and precision for the quantitative methods as specified in the value sheet. For professional use in clinical laboratory only.

PRINCIPLE

This human serum control is to be used in the same way as an unknown specimen that would be used in accordance with the direction of the assay procedure. The results obtained for the control are to be compared with the assigned values given in the "ASSAYED VALUES" section of the insert and an evaluation made by standard statistical techniques to determine if the procedure is within the control limits.

DESCRIPTION AND COMPOSITION

The control is diluted human serum containing low levels of CRP with phosphate buffered saline. Liquid stabilized. Preservative: 0.095 % sodium azide.

ASSAYED VALUES

Lot	Expiry	Parameter	Value	Interval	1 SD	Unit
XXXXXXXXXX	YYYY-MM	CRP	0.00	0.00–0.00	0.000	mg/dL
			00.0	00.0–00.0	0.00	mg/L

Assayed values for XL-Erba analysers with calibration using automatic dilution – see the end of the document.

UNIT CONVERSION

mg/dL × 10 = mg/L

CONTROL PREPARATION

Control is liquid, ready to use.

STORAGE AND HANDLING

The unopened CRP CON L stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. After first opening the vial, control can be used for 6 weeks if stored tightly closed at 2–8 °C after use. Do not freeze.

MATERIALS REQUIRED, BUT NOT PROVIDED

CRP, Cat. No. XSYS0047
CRP CAL SH, Cat. No. XSYS053
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

ASSAY PROCEDURE

Pour the required volume into a sample cup and analyze in the same way as patient samples. The controls should be run daily in parallel with patient samples and after every calibration. Control intervals must be adapted to individual laboratory's requirements. Follow the applicable government regulations and local guidelines for quality control.

EVALUATION

The XL-analysers perform the evaluation automatically.

EXPECTED RESULTS

The target value stated in the "ASSAYED VALUES" section. Determination was performed under strictly standardized conditions on Erba analysers using Erba system reagents and the calibrator CRP CAL SH.

The target value specified is the mean of all values obtained. The corresponding control range is calculated as the target value ±3 standard deviations. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material. Results should be within the defined ranges.

For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range.

TRACEABILITY

This control has been standardized against the reference material ERM-DA474/IFCC.

PERFORMANCE CHARACTERISTICS

Performance characteristics of individual analytes on the Erba XL automatic systems are available on websites www.erba.com.

LIMITATIONS:

Erroneous results can occur from technical errors associated with the assay procedure.

Improper storage or handling of the control can also affect the results. Do not use the control if there is visible evidence of microbial growth in the vial.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Health protection

Serum used in the manufacture of this product has been tested and found non-reactive for hepatitis B surface antigen (HbsAg), antibody to Hepatitis C (HCV), HIV-1 p24 antigen, and antibody to HIV-1/2. The test procedures do not guarantee that all infectious agents will be detected. Because no test method can offer complete assurance that Hepatitis B virus Hepatitis C virus and HIV-1/2 or other infectious agents are absent, the material should be handled as potentially infectious. During the work all principles of general safety as well as the principles of work with biological materials are to be maintained. Above all observe the principles of personal hygiene, use personal protection aids (clothes, shoes, etc.), do not eat, drink and smoke. After finishing the work and before a meal wash your hands with soap and warm water and apply a suitable regenerating cream. When the bottle is broken, the contaminated surface must be disinfected with Persteril or Chloramine.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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