ALPHA AMYLASE (SINGLE REAGENT)

Cat. No.	Pack Name	Packaging (Content)
XSYS0003	AMY 110	R1: 5 × 22 mL, RFID tag, instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of α -amylase in human serum, plasma and urine on automatic systems ERBA XL. In combination with other parameters it is intended mainly for monitoring and diagnosis of acute pancreatitis and acute abdominal disorders. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

The α -amylases are hydrolytic enzymes, which break down starch into maltose. In the human body α -amylases originate from various organs: the pancreatic amylase is produced by the pancreas and released into the intestinal tract, the salivary amylase is synthesized in the salivary glands and secreted into saliva. The amylase present in the blood is eliminated through the kidney and excreted into the urine. Therefore, elevation of serum activity is reflected in a rise of urinary amylase activity. Because of the sparsity of specific clinical symptoms of pancreatic diseases, α -amylase determinations are of considerable importance in pancreatic diagnostics. They are mainly used in the diagnosis and monitoring of acute pancreatitis. Hyperamylasemia does not, however, only occur with acute pancreatitis or in the inflammatory phase of chronic pancreatitis, but also in renal failure (reduced glomerular filtration), pulmonary inflammation, diseases of the salivary gland, diabetic ketoacidosis, cerebral trauma, surgical interventions or in the case of macroamylasemia. To confirm pancreatic specificity, it is recommended that an additional pancreas-specific enzyme, lipase or pancreatic α -amylase, also be determined.

PRINCIPLE

The direct amylase assay involves the use of a chromogenic substrate 2-chloro-4-nitrophenyl- β -1-4-galactopyranosylmaltotriose (CNP3)¹. The α -amylase hydrolyzes the substrate CNP3, releasing 2-chloro-4-nitrophenol (CNP) and forming 2-chloro-4-nitrophenyl- β -1-4-galactopyranosylmaltotriose (CNP2), maltotriose (G3) and glucose (G). Reaction scheme is following:

The rate of CNP formation is proportional to the activity of α -amylase present in the sample and can be measured kinetically at 405 or 415 nm.

REAGENT DESCRIPTION AND COMPOSITION

R1	
MES buffer	50 mmol/L
Calcium chloride	3.81 mmol/L
Sodium chloride	300 mmol/L
Potassium thiocyanate	450 mmol/L
Sodium azide	0.90 g/L
CNP3	0.91 mmol/L

COMPOSITION OF REACTION MIXTURE

MES buffer	49 mmol/L
Calcium chloride	3.74 mmol/L
Sodium chloride	294 mmol/L
Potassium thiocyanate	441 mmol/L
Sodium azide	0.88 g/L
CNP3	0.89 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Urine: Collect urine without additives. α -amylase is unstable in acid urine. Assay promptly or adjust pH to alkaline range (just above pH 7) before storage⁸.

Stability in serum / plasma⁹:	7 days at 15–25 °C
	7 days at 2–8 °C
	1 year at -20 °C

Stability in urine⁹:	2 days at 15–25 °C
	10 days at 2–8 °C
	3 weeks at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method¹⁰.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	Rate A
Curve type	Linear
Wavelength (prim. / sec.)	405 or 415 / 700 nm
Reading time	100–200 s after adding of R1
Reaction direction	Increase
Unit	U/L (μkat/L)

Reagent volumes

R1 200 μL

Sample 4 μL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

UNIT CONVERSION

U/L × 0.0167 = μkat/L

EXPECTED VALUES¹¹

At 37 °C

Serum: up to 80 U/L

Urine: up to 500 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 4.18 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 1549 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	56.0	0.38	0.67
Sample 2	153.1	0.91	0.60

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	56.2	1.57	2.80
Sample 2	153.5	3.68	2.39

Accuracy

Two different validated control materials were used. Determined bias is -7.1 % at the target value 305.0 U/L and -5.6 % at the target value 499.4 U/L.

Comparison

A comparison between XL-640 automatic system ALPHAAMYLASE (SINGLE REAGENT) (y) and a commercially available test (x) using 40 samples gave following results:

Linear regression:

$y = 0.979x - 5.58 \text{ U/L}$ $r = 0.986$

Passing-Bablok¹²:

$y = 0.781x + 2.53 \text{ U/L}$ $r = 0.920$

Interferences

Criterion: Recovery within ±10 % of initial value of α -amylase activity in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs:

Serum: No interference was found at therapeutic concentrations using common drug panels. Exception: Icodextrin-based drugs may lead to decreased amylase results¹³.

Urine: No interference was found at therapeutic concentrations using common drug panels. Exception: Acetylcysteine, Ascorbic acid, Cefoxitin, Gentamicin sulphate, Levodopa, Methylidopa and Ofloxacin may lead to decreased amylase results¹³.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 405 nm against the distilled water is 0.14.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of α -amylase. Some drugs can also interfere. See paragraph Interferences.
- Saliva and skin contain α -amylase therefore never pipette reagents by mouth and avoid contamination of samples and reagents. Even trace contamination can affected results.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety and Performance is available at Eudamed see:

<https://ec.europa.eu/tools/eudamed>.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: 6RYE-AE3K-3H7A-Q7QR

Danger

Contains: potassium thiocyanate

Hazard statement:

H318 Causes serious eye damage.

Precautionary statement:

P280 Wear eye protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Chavez RG, et al, U.S. Patent 4,963,479.
2. Winn-Deen ES, David H, Sigler H, Chavez R, Development of a direct assay for α -amylase, Clin. Chem. 34 (1988), p. 2005-2008
3. J. F. Ziva, and P. R. Pannall, "Plasma Enzymes in Diagnosis" in Clinical Chemistry in Diagnosis and Treatment. Lloyd London 1979 : Chapter XV : 341-2.
4. Foo, Y.A. and Brosalki, S.B. Ann. Clin. Biochem. 1986; 23: 624-37.
5. Bais, R. Am. Jnl. of Clin. Path. 1982; 78 : 184-8.
6. Lorentz K. α -Amylase. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft;1998.p.192-202.
7. Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company;1999.p.689-98.
8. Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia PA: WB Saunders Company 1995;46-51.
9. Junge W, Wortmann W, Wilke B, Waldenstroem J et al. Development and evaluation of assays for determination of total and pancreatic amylase at 37°C according to the principle recommended by the IFCC. Clin Biochem 2001; 34: 607-15.
10. Lorentz K. Approved recommendation on EFCC methods for the measurement of catalytic concentration of enzymes. Part 9. IFCC method for α -amylase (1,4- α -D-Glucan-4-glucanohydrolase, EC 3.2.1.1). Clin Chem Lab Med 1998; 38:185-203.
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
12. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
13. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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