CHLORIDE

Cat. No.	Pack Name	Packaging (Content)
XSYS0008	CL 120	R1: 10 × 12 mL, RFID tag, instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of chloride in human serum, plasma and urine on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis in case of suspicion of hyperchloridemia or hypochloridemia. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Chloride is the most important anion in serum besides bicarbonate. Together with sodium it is an essential osmotically active component in plasma which is involved in maintenance of water distribution and anion-cation-balance. Serum concentrations of chloride behave parallel to sodium levels and reciprocally to bicarbonate. Increased chloride values occur in dehydration, metabolic acidosis related with prolonged diarrhoea and bicarbonate loss, renal insufficiencies and endocrinological disorders as reduced or increased adrenal function. Decreased values are observed in metabolic acidosis with increased production of organic acids, salt-losing nephritis and excessive sweating.

PRINCIPLE

When chloride is mixed with a solution of undissociated mercuric thiocyanate, the chloride preferentially combines with mercury forming mercuric chloride. The thiocyanate that is released then combines with ferric ions present in the solution forming strongly coloured ferric thiocyanate¹.

The resulting absorbance at 505 nm is proportional to the chloride concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	Concentration
Mercuric thiocyanate	1.01 mmol/L
Ferric nitrate	37.1 mmol/L
Nitric acid	62.4 mmol/L
Mercuric nitrate	0.105 mmol/L

COMPOSITION OF REACTION MIXTURE

Mercuric thiocyanate	1.00 mmol/L
Ferric nitrate	36.7 mmol/L
Zinc sulfate	61.8 mmol/L
Nitric acid	0.104 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.
 Plasma: heparin
 Urine.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.
 See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma / urine⁴: 7 days at 20–25 °C
 7 days at 4–8 °C
 at least one year at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank.
 Calibration frequency: 22 days

Calibration is needed:
 • after reagent lot change
 • as required by internal quality control procedures
 • calibration interval may be extended based on acceptable verification of calibration by the laboratory.

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.
 The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against coulometry.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	505/660 nm
Reading time	90–150 s after adding of R1
Reaction direction	Increase
Unit	mmol/L (mEq/L)
Reagent volumes	
R1	200 µL
Sample	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

UNIT CONVERSION

mmol/L = mEq/L

EXPECTED VALUES⁵

Serum:

Adult	98–107 mmol/L
>90 y	98–111 mmol/L
Cord	96–104 mmol/L
Premature	95–110 mmol/L
0–30 d	98–113 mmol/L

Urine, 24 h:

Adult	110–250 mmol/24 h
>60 y	95–195 mmol/L
Infant	2–10 mmol/24 h
<6 y	15–40 mmol/24 h
6–10 y	
Female	18–74 mmol/24 h
Male	36–110 mmol/24 h
10–14 y	
Female	36–173 mmol/24 h
Male	64–176 mmol/24 h

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 0.43 mmol/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 200 mmol/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mmol/L)	SD (mmol/L)	CV (%)	Intermediate precision (serum)	Mean (mmol/L)	SD (mmol/L)	CV (%)
Sample 1	110.1	0.79	0.71	Sample 1	109.6	2.64	2.41
Sample 2	119.6	0.87	0.73	Sample 2	119.5	2.35	1.96

Repeatability (urine)	Mean (mmol/L)	SD (mmol/L)	CV (%)	Intermediate precision (urine)	Mean (mmol/L)	SD (mmol/L)	CV (%)
Sample 1	91.3	1.06	1.16	Sample 1	122.1	2.07	1.70
Sample 2	120.1	1.09	0.91	Sample 2	142.4	3.17	2.23

Accuracy

Two different validated control materials for serum were used. Determined bias is 4.1 % at the target value 112 mmol/L and -2.4 % at the target value 131 mmol/L.

Comparison

A comparison between XL-640 automatic system CHLORIDE (y) and a commercially available test (x) using 145 samples (serum) gave following results:

Linear regression:
 $y = 0.912x + 9.555$ mmol/L $r = 0.968$

Passing-Bablok⁶:
 $y = 0.956x + 5.156$ mmol/L $r = 0.964$

Interferences

Criterion: Recovery within ±10 % of initial value of chloride concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 12 mg/dL, triglycerides up to 850 mg/dL.

Drugs: aspirin and bicarbonate may increase results in serum, bromides and iodides may increase results in serum and urine⁷.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of chlorides. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1
 UFI: 058F-WERG-3HT7-HPU2

Danger

Contains: methanol, mercury dithiocyanate, mercury dinitrate

Hazard statement:

H302 Harmful if swallowed.
 H314 Causes severe skin burns and eye damage.
 H371 May cause damage to organs.

Precautionary statement:

P260 Do not breathe vapours.
 P280 Wear protective gloves/protective clothing/eye protection.
 P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
 P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P308 + P313 IF exposed or concerned: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Levinson, S.S.: Direct determination of serum chloride with a semiautomated discrete analyzer. Clin. Chem. 22: 273–274, 1976.
2. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 295–298, 1998.
3. Scott GS, Heusel JW, LeGrys VA, Siggard-Andersen O. Electrolytes and blood gases. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company: 1056–1094, 1999.
4. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag, 22–23, 2001.
5. Tietz Textbook of Clinical Chemistry and Molecular diagnostics. Burtis, C.A., Ashwood, E. R., Bruns, D.E.; 5th edition, WB Saunders Company, 2012.
6. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov; 26(11): 783–790.
7. Tietz Clinical Guide to Laboratory Tests 4th edition, Alan HB Wu, WB Saunders Company, 2006.

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