ALKALINE PHOSPHATASE

Cat. No.	Pack Name	Packaging (Content)
XSYS0002	ALP 110	R1: 2 × 44 mL, R2: 2 × 11 mL, RFID tag, instruction for use, chimneys: 2 pcs

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of alkaline phosphatase in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver disease and bone disease. For professional use in clinical laboratories only.


CLINICAL SIGNIFICANCE

Human ALP consists of a group of enzymes which hydrolyse phosphates at an alkaline pH. Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from the germ cells. It occurs in osteoblasts, hepatocytes, leukocytes, the kidneys, spleen, placenta, prostate and the small intestine. The liver-bone-kidney type is particularly important.

A rise in the alkaline phosphatase occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia or fractures. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth.

PRINCIPLE

The products are based on the principle of photometric determination according to IFCC recommendation¹. This method utilises 4-nitrophenyl phosphate as the substrate. Under optimised conditions ALP present in the sample catalyses the following reaction:

At the pH of the reaction, 4-nitrophenol has an intense yellow colour. The reagent also contains a metal ion buffer system to ensure that optimal concentrations of zinc and magnesium are maintained. The metal ion buffer can also chelate other potentially inhibitory ions which may be present. The reaction is monitored by measuring the rate of increase in absorbance at 405 or 415 nm which is proportional to the activity of ALP in the serum.

REAGENT DESCRIPTION AND COMPOSITION

R1	R2
AMP buffer, pH 10.5	AMP buffer, pH 9.5
Magnesium acetate	Magnesium acetate
Zinc sulfate	Zinc sulfate
HEDTA	HEDTA
Sodium azide	p-nitrophenyl phosphate
	Sodium azide

COMPOSITION OF REACTION MIXTURE

AMP buffer	350 mmol/L
Magnesium acetate	2.00 mmol/L
Zinc sulfate	1.00 mmol/L
HEDTA	2.00 mmol/L
p-nitrophenyl phosphate	16.0 mmol/L
Sodium azide	0.9 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 32 days if refrigerated (2–10 °C) and not contaminated and if chimneys are used.

Absorption of atmospheric CO₂ in open reagent bottle R1 leads to the reduced stability of reagents. To minimize the rate of absorption of CO₂ the use of "chimneys" provided with the kit is recommended. Use of chimneys has shown improvement in calibration stability, provided the storage and use of reagents is as per instructions. Chimneys provided with the kit are to be placed directly into the reagent bottle immediately after opening the cap. Chimneys can be re-used again for the bottle of the same kit. Before using the chimney in a new bottle, it is recommended to wash chimneys well with distilled water and dry them thoroughly.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases.

When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma²: 7 days at 20–25 °C
 7 days at 4–8 °C
 2 months at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 12 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to original method IFCC¹.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	Rate A
Curve type	Linear
Wavelength (prim. / sec.)	405 (415) / 700 nm
Reading time	60–240 s after adding of R2
Reaction direction	Increase
Unit	U/L (µkat/L)
Reagent volumes	
R1	160 µL
R2	40 µL
Sample	4 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION

U/L × 0.0167 = µkat/L

EXPECTED VALUES³

At 37 °C

Females: 4–15 years:	54–369 U/L	Males: 1–12 years:	54–369 U/L
20–50 years:	42–98 U/L	20–50 years:	53–128 U/L
≥60 years:	53–141 U/L	≥60 years:	56–119 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section are representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com

Limit of quantification: 4.94 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 1617 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	125	0.55	0.44
Sample 2	374	1.41	0.38

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	130	3.77	2.90
Sample 2	392	7.86	2.01

Accuracy

Two different validated control materials were used. Determined bias is -4.1 % at the target value 163.8 U/L and -1.6 % at the target value 392.1 U/L.

Comparison

A comparison between XL-640 automatic system ALP (y) and a commercially available test (x) using 47 samples gave following results:

Linear regression:
 $y = 1.049x - 0.18 \text{ U/L}$ $r = 1.000$
 Passing-Bablok⁷:
 $y = 1.048x - 0.03 \text{ U/L}$ $r = 1.000$

Interferences

Criterion: Recovery within ±10 % of initial value of an alkaline phosphatase activity in sample without interfering substance. Following substances do not interfere: haemoglobin up to 9 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of ALP. See paragraph Interferences.
- Drugs: No interference was found at therapeutic concentrations using common drug panels⁸.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

The reagent is not classified as dangerous.

R2

UFI: 5FAQ-SJAM-HE7A-6Q8K

Warning

Hazard statement:

H315 Causes skin irritation.
 H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.
 P302 + P352 IF ON SKIN: Wash with plenty of water.
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
 e-mail: diagnostics@erba.com, www.erba.com

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1. IFCC method for the measurement of ALP J. Clin. Chem. Clin. Biochem. 1983; 21: 731-48.
2. Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
3. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
4. Zilva JF, Pannall PR, "Plasma Enzymes in Diagnosis" in Clinical Chemistry in Diagnosis and Treatment. Lloyd London 1979: Chapter 15: 343.
5. Young DS. Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990: 3: 19-25.
6. Kaplan and Pesce (Eds.) Clinical Chemistry, Theory analysis and correlation. Second Edition. CV Mosby Co. 1989.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
8. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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