TOTAL PROTEIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0018	TP 440	R1: 10 × 44 mL, RFID tag, instruction for use, chimney: 1pc

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of total protein in human serum and plasma on automatic systems ERBA XL. Intended for monitoring of gross changes in protein levels caused by various disease states. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE


Plasma proteins derive primarily from synthesis in the liver, plasma cells, lymph nodes, spleen, and bone marrow. In disease states both the total plasma protein level and the ratio of the individual fractions may be dramatically altered from their normal values.

Hypoproteinemia may be caused by such conditions as nephrotic syndrome, extensive bleeding, sprue (deficient protein absorption), severe burns, salt retention syndromes, and Kwashiorkor (acute protein starvation).

Hyperproteinemia may be observed in cases of severe dehydration and disease states such as multiple myeloma. Changes in the proportions of the plasma proteins may occur in one or several of the protein fractions and often without alterations in the quantity of the total protein. The A/G ratio has commonly been used as an index of the distribution between the albumin and globulin fractions. This ratio can be significantly altered in such conditions as cirrhosis of the liver, glomerulonephritis, nephrotic syndrome, acute hepatitis, lupus erythematosus, and in some acute and chronic infections.

PRINCIPLE

Biuret method^{1,2,3,4,5}. The peptide bonds of protein react with copper (II) ions in alkaline solution to form a blue violet ion complex, the so called biuret reaction, each copper ion complexing with 5 or 6 peptide bonds. Tartrate is added as a stabiliser whilst iodide is used to prevent auto-reduction of the alkaline copper complex.

Absorbance of the blue violet complex measured at 546 nm is proportional to the protein concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
Copper sulphate	12 mmol/L
Potassium sodium tartrate	31.9 mmol/L
Potassium iodide	30.1 mmol/L
Sodium hydroxide	600 mmol/L

COMPOSITION OF REACTION MIXTURE

Copper sulphate	11.8 mmol/L
Potassium sodium tartrate	31.3 mmol/L
Potassium iodide	29.5 mmol/L
Sodium hydroxide	588 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated and if chimneys are used.

Absorption of atmospheric CO₂ in open reagent bottle R1 leads to the reduced stability of reagents. To minimize the rate of absorption of CO₂, the use of "chimneys" provided with the kit is recommended. Use of chimneys has shown improvement in calibration stability, provided the storage and use of reagents is as per instructions. Chimneys provided with the kit are to be placed directly into the reagent bottle immediately after opening the cap. Chimneys can be re-used again for the bottle of the same kit. Before using the chimney in a new bottle, it is recommended to wash chimneys well with distilled water and dry them thoroughly.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.
 For specimen collection and preparation only use suitable tubes or collection containers.
 Only the specimens listed below were tested and found acceptable.
 Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.
 See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma:	6 days at	20–25 °C
	4 weeks at	4–8 °C
	1 year at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank.
 Calibration frequency: 30 days
 Calibration is needed:
 • after reagent lot change
 • as required by internal quality control procedures
 • calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the SRM 927c reference material.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	546 nm
Reading time	10 min. after adding of R1
Reaction direction	Increase
Unit	g/dL (g/L)
Reagent volumes	
R1	200 µL
Sample	4 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

UNIT CONVERSION

g/dL × 10 = g/L

EXPECTED VALUES*

In Serum:

Adults	6.4–8.3 g/dL
Premature	3.6–6.0 g/dL
Newborn	4.6–7.0 g/dL
1 week	4.4–7.6 g/dL
7–12 months	5.1–7.3 g/dL
1–2 years	5.6–7.5 g/dL
>2 years	6.0–8.0 g/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 0.26 g/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

Linearity: 14.0 g/dL

Linearity is the highest measured activity with recovery within ±10% from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	6.06	0.056	0.93
Sample 2	9.39	0.041	0.43

Intermediate precision	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	5.99	0.114	1.91
Sample 2	9.35	0.158	1.69

Accuracy

Two different validated control materials were used. Determined bias is -2.4% at the target value 7.21 g/dL and 1.7% at the target value 8.43 g/dL.

Comparison

A comparison between XL-640 automatic system TOTAL PROTEIN (y) and a commercially available test (x) using 150 samples gave following results:

Linear regression:
 $y = 0.955x + 0.183 \text{ g/dL}$ $r = 0.987$
 Passing-Bablok⁶:
 $y = 0.964x + 0.135 \text{ g/dL}$ $r = 0.978$

Interferences

Criterion: Recovery within ±10% of initial value of total protein concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 3 g/L, bilirubin up to 28 mg/dL, triglycerides up to 320 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels¹⁰.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of total protein. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.
 Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1
 UFI: ARHK-VJ1C-KE7U-9PC4

Danger

Contains: sodium hydroxide

Hazard statement:

H314 Causes severe skin burns and eye damage.
 H412 Harmful to aquatic life with long lasting effects.

Precautionary statement:

P260 Do not breathe vapours.
 P280 Wear protective gloves / protective clothing / eye protection / face protection.
 P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
 P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
 P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Weichselbaum TE., An accurate and rapid method for the determination of proteins in small amounts of blood serum and plasma. Am J Clin Pathol, 10, 40-49 (1946)
2. Cornall, A. G., Bardawill, C. J., David, M. M., Determination of serum proteins by means of the biuret reaction, J. Biol. Chem. 177, 751, 1949.
3. Dumas, B. T., Bayse, D. D., Carter, R. J., Peters, T., Jr., and Schaffer, R., A candidate Reference Method for determination of total protein in serum. I. Development and validation. Clin. Chem. 27, 1642-1650 (1981).
4. Chromy, V., and Fischer, J., Photometric determination of total protein in lipemic sera. Clin. Chem. 23, 754-756 (1977).
5. Chromý, V., Fischer, J., Vozniček, Assay of total serum protein by the biuret method, reaction rate and lipemia, J. Z. Med. Labor.-Diagn. 21, 333-337, 1980.
6. Brobeck JR, ed. Physiological Basis of Medical Practice, 9th ed. Baltimore, MD: Wilkins and Wilkins 1973;4-7.
7. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
8. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11): 783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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LITERATÚRA

1. Weichselbaum TE., An accurate and rapid method for the determination of proteins in small amounts of blood serum and plasma. Am J Clin Pathol,10, 40-49 (1946)
2. Cornall, A. G., Bardawill, C. J., David, M. M., Determination of serum proteins by means of the biuret reaction, J. Biol. Chem. 177, 751, 1949.
3. Dumas, B. T., Bayse, D. D., Carter, R. J., Peters, T., Jr., and Schaffer, R., A candidate Reference Method for determination of total protein in serum. I. Development and validation. Clin. Chem. 27, 1642-1650 (1981).
4. Chromý, V., and Fischer, J., Photometric determination of total protein in lipemic sera. Clin. Chem. 23, 754-756 (1977).
5. Chromý, V., Fischer, J., Vozniček, Assay of total serum protein by the biuret method, reaction rate and lipemia, J. Z. Med. Labor.-Diagn. 21, 333-337, 1980.
6. Brobeck JR, ed. Physiological Basis of Medical Practice, 9th ed. Baltimore, MD: Wilkins and Wilkins 1973;4-7.
7. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
8. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11): 783-790.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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