# MICROPROTEIN

Cat. No.	Pack Name	Packaging (Content)
XYS0027	MP 120	R1: 10 × 12 mL, R2 standard: 1 × 5 mL, RFID tag, instruction for use

Reaction direction	Increase
Unit	mg/dL (mg/L)
Reagent volumes	
R1	200 µL
Sample volumes	4 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:sample does not change.

## UNIT CONVERSION

mg/dL × 10 = mg/L

## EXPECTED VALUES<sup>9</sup>

### Urine, 24 hours:

Adult 1–14 mg/dL

### Excretion

Adult <100 mg/day

Pregnancy <150 mg/day

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on [www.erba.com](http://www.erba.com).

### Limit of quantification:

1.57 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30).

### Linearity:

210 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

### Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision	SD (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	23.1	0.41	1.78	Sample 1	21.5	0.77	3.57
Sample 2	66.8	1.03	1.54	Sample 2	63.2	1.57	2.49

### Accuracy

Two different validated control materials were used. Determined bias is 8.9 % at the target value 18.9 mg/dL and 6.4 % at the target value 67.9 mg/dL.

### Comparison

A comparison between XL-640 automatic system MICROPROTEIN (y) and a commercially available test (x) using 120 samples gave following results:

Linear regression:

$$y = 0.908x + 4.882 \text{ mg/dL} \quad r = 0.997$$

Passing-Bablok<sup>10</sup>:

$$y = 0.917x + 4.138 \text{ mg/dL} \quad r = 0.997$$

### Interferences

In this method, some kinds of surface active agents may affect the color. Cationic surfactants, in general, give the same color as proteins. Because anions inhibit the color reaction, wash the equipment thoroughly, using distilled water, until no surface active agent remains. Then dry equipment completely before using it.

Haemoglobin interferes<sup>11</sup>. Neomycin and Chlorpromazine may cause false-positive results<sup>12</sup>.

Small amounts of protein attached to the cuvette wall after measurement of certain other tests will cause an erroneously high measured value when the test solution is transferred to the cuvette. If this should occur, wash the cuvette completely and measure again.

### Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- Surfactants, hemoglobin or some drugs in the sample can interfere with MP determination. See paragraph interference.

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

#### R1

UFI: GA8F-XE48-QH7T-UD06

#### Warning

Contains: methanol

#### Hazard Statements:

H371 May cause damage to organs.

#### Precautionary Statements:

P260 Do not breathe vapours/spray.

P264 Wash hands thoroughly after handling.

P308 + P311 IF exposed or concerned: Call a POISON CENTER or doctor.

#### R2

Reagent is not classified as dangerous.

### WASTE MANAGEMENT

Please refer to local legal requirements.

## INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of protein in urine on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of kidney function disorders. For professional use in clinical laboratory only.

## CLINICAL SIGNIFICANCE

The role of the renal system in the conversion of plasma proteins has been recognized for some time. Under normal physiological conditions small molecular weight proteins such as insulin pass through the glomeruli in relatively large amounts. Intermediate size proteins such as Transferrin and Albumin also pass through but only in relatively small amounts. Most of these proteins are reabsorbed in the renal tubules such that normal urine contains less than 150 mg of protein per day. This also includes the protein of non-serum origin normally secreted by the distal tubule (mucoprotein) and collecting ducts. Increased levels of urinary protein, (proteinuria) usually more than 0.15 g per 24 hours (150 mg/24 hours), almost always indicates disease. Proteinuria may be classified as renal proteinurea or proteinuria with normal renal function. Renal proteinuria may be further classified as Glomerular or tubular proteinuria.

Glomerular proteinuria is due to increased glomerular permeability (nephrotic syndrome) and may be seen in glomerular nephritis or secondary to other diseases such as diabetic nephropathy. Albumin is usually the predominant protein in the urine. Tubular proteinuria may be due to renal tubular damage from any cause especially pyelonephritis. Tubular proteinuria results in modest increases in the low molecular weight proteins. Increased protein excretion is seen during normal pregnancy, after strenuous exercise or following prolonged maintenance of an upright posture. Increase in low molecular weight proteins may be due to the production of Bence Jones protein, haemoglobinuria as a result of severe haemolysis and myoglobinuria as a result of severe muscle damage.

## PRINCIPLE

Direct, colorimetric method with pyrogallol red. Pyrogallol red is combined with molybdate to form a red complex with a maximum absorbance at 470 nm. The assay is based on the shift in absorbance that occurs when the pyrogallol red-molybdate complex binds basic amino groups of protein molecules. Under the conditions of the test in the presence of protein, a blue-purple complex is formed with a maximum absorbance at 600 nm<sup>1,2,3,4</sup>.

Protein + Pyrogallol red-molybdate complex → Blue-purple complex

The absorbance of blue-purple complex measured at 600 nm is directly proportional to the protein concentration in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

### R1

Succinate buffer (pH 2.5) 15 mmol/L

Pyrogallol red 60 µmol/L

Sodium molybdate 43 µmol/L

### R2 standard

Protein See bottle label

## COMPOSITION OF REACTION MIXTURE

Succinate buffer (pH 2.5) 14.7 mmol/L

Pyrogallol red 5.8 mmol/L

Sodium molybdate 42.2 µmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

## MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

All control solutions with MP (e.g. Bio-Rad Liquicheck, Level 1 and Level 2).

Erba XL analysers: XL-200, Cat. No. INS00002

XL-640, Cat. No. INS00008

XL-1000, Cat. No. INS00010

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Urine: Use random or 24-hour urine specimens. Use no preservatives. Refrigerate specimen during collection.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

### Stability in urine<sup>5,6</sup>:

1 day at 15–25 °C

7 days at 4–8 °C

1 month at -20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with R2 standard is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 7 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

## QUALITY CONTROL

For quality control all control solutions with protein total values determined by this method can be used. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method and R2 standard have been standardized against the reference material NIST SRM 927.

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
1. Tietz NW. Textbook of clinical chemistry, 2nd ed. Philadelphia: WB Saunders Company, 1994.
2. Fujita Y, Mori I, Kitano S. Color reaction between pyrogallol red molybdenum (VI) complex and protein. Bunseki Kagaku 32: 379-386, 1983.
3. Watanabe N, Kamei S, Ohkubo A, Yamanaka M, Ohsawa S, Makino K, Tokuda K. Urinary protein as measured with pyrogallol red-molybdate complex manually and in Hitachi 726 automated analyzer. Clin Chem 32: 1551-1544, 1986.
4. Orsonneau JL, Douet P, Massoubre C, Lustenberger P, Bernard S. An improved pyrogallol red-molybdate method for determining total urinary protein. Clin Chem 35: 2233-2236, 1989.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 52-53, 2001.
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
7. Johnson AM, Rohlfis EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 477-540, 1999.
8. Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 1308-1326, 1998.
9. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
10. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
11. Yilmaz FM, Yücel D. Effect of Addition of Hemolysate on Urine and Cerebrospinal Fluid Assays for Protein. Clin Chem 52:152-153, 2006.
12. Fujita Y, Mori I, Kitano S, Color Reaction Between Pyrogallol Red - Molybdenum(VI) Complex and Protein, Bunseki Kagaku 32, E379-E386, 1983.

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## LITERATÚRA

1. Tietz NW. Textbook of clinical chemistry, 2nd ed. Philadelphia: WB Saunders Company, 1994.
2. Fujita Y, Mori I, Kitano S. Color reaction between pyrogallol red molybdenum (VI) complex and protein. Bunseki Kagaku 32: 379-386, 1983.
3. Watanabe N, Kamei S, Ohkubo A, Yamanaka M, Ohsawa S, Makino K, Tokuda K. Urinary protein as measured with pyrogallol red-molybdate complex manually and in Hitachi 726 automated analyzer. Clin Chem 32: 1551-1544, 1986.
4. Orsonneau JL, Douet P, Massoubre C, Lustenberger P, Bernard S. An improved pyrogallol red-molybdate method for determining total urinary protein. Clin Chem 35: 2233-2236, 1989.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 52-53, 2001.
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
7. Johnson AM, Rohlfis EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 477-540, 1999.
8. Felgenhauer K. Laboratory diagnosis of neurological diseases. In: Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 1308-1326, 1998.
9. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
10. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
11. Yilmaz FM, Yücel D. Effect of Addition of Hemolysate on Urine and Cerebrospinal Fluid Assays for Protein. Clin Chem 52:152-153, 2006.
12. Fujita Y, Mori I, Kitano S, Color Reaction Between Pyrogallol Red - Molybdenum(VI) Complex and Protein, Bunseki Kagaku 32, E379-E386, 1983.

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