LDL DIRECT

Cat. No.	Pack Name	Packaging (Content)
XSYS0044	LDL C 80	R1: 2 × 30 mL, R2: 2 × 10 mL, RFID tag, instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of LDL cholesterol in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of coronary heart disease and lipoprotein metabolism disorders. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Low Density Lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on triglyceride-rich Very Low Density Lipoproteins (VLDLs). Specific LDL receptors exist to facilitate the elimination of LDL from plasma by liver parenchymal cells. It has been shown that most of the cholesterol stored in atherosclerotic plaques originates from LDL. For this reason the LDL Cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis. Accurate measurement of LDL Cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture.

PRINCIPLE

The assay is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol-methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL, and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHO) and cholesterol esterase (CHE), whereas HDL reacts with the enzymes. Addition of R2 containing a specific detergent releases LDL from the PVS/PEGME complex. The released LDL reacts with the enzymes to produce hydrogen peroxide (H₂O₂) which is quantified by the Trinder's reaction^{12,21}.

The absorbance of the produced quinone dye at 600 nm is proportional to the LDL concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
MES buffer (pH 6.5)	50 mmol/L
Polyvinyl sulfonic acid	50 mg/L
Polyethylene-glycol-methyl ether	30 mL
4-aminoantipyrine, 4-AAP	0.9 g/L
Cholesterol esterase	5 KU/L
Cholesterol oxidase	20 KU/L
Peroxidase	5 KU/L
Detergent	

R2	
MES buffer (pH 6.5)	50 mmol/L
N-Bis(4-sulfobutyl)-3-methylaniline, TODB	5 mmol/L
Detergent	

COMPOSITION OF REACTION MIXTURE

MES buffer (pH 6.5)	49.6 mmol/L
Polyvinyl sulfonic acid	37 mmol/L
Polyethylene-glycol-methyl ether	22 g/L
4-aminoantipyrine, 4-AAP	0.7 mg/L
N-Bis(4-sulfobutyl)-3-methylaniline, TODB	0.7 mL
Cholesterol esterase	3.7 mmol/L
Cholesterol oxidase	14.9 KU/L
Peroxidase	3.7 KU/L
Detergent	

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

HDL/LDL Calibrator, Cat. No. XSYS0061
 ERBA NORM 4x5, Cat. No. BLT00080
 ERBA NORM 10x5, Cat. No. XSYS0123
 ERBA PATH 4x5, Cat. No. BLT00081
 ERBA PATH 10x5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma.

Fasting and non-fasting samples can be used¹.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay. See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ² :	1 day at	20–25 °C
	7 days at	4–8 °C
	3 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with HDL/LDL Calibrator is recommended. 2 point calibration (blank and calibrator); distilled water is recommended as blank. Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, HDL/LDL calibrator and controls ERBA NORM and ERBA PATH have been standardized against the reference material NIST SRM 1951.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	2-Point
Curve type	Linear
Wavelength (prim. / sec.)	600/700 nm
Reading time 1	just before adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (mmol/L)
Reagent volumes	
R1	180 µL
R2	60 µL
Sample volumes	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION

mg/dL × 0.026 = mmol/L

EXPECTED VALUES¹⁶

In serum:	Male	Female
5–9 y	63–129	68–140 mg/dL
10–14 y	64–133	68–136 mg/dL
15–19 y	62–130	59–137 mg/dL
20–24 y	66–147	57–159 mg/dL
25–29 y	70–165	71–164 mg/dL
30–34 y	78–185	70–156 mg/dL
35–39 y	81–189	75–172 mg/dL
40–44 y	87–186	74–174 mg/dL
45–49 y	97–202	79–186 mg/dL
50–54 y	89–197	88–201 mg/dL
55–59 y	88–203	89–210 mg/dL
60–64 y	83–210	100–224 mg/dL
65–69 y	98–210	92–221 mg/dL
>69 y	88–186	96–206 mg/dL

Coronary heart disease risk:

Optimal	<100 mg/dL
Near/above optimal	100–129 mg/dL
Borderline high	130–159 mg/dL
High	160–189 mg/dL
Very high	>189 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 0.46 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 413 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	73.7	0.59	0.80	Sample 1	73.3	1.59	2.17
Sample 2	132.2	0.72	0.54	Sample 2	136.9	3.34	2.44

Accuracy

Two different validated control materials were used. Determined bias is 9.9 % at the target value 128.8 mg/dL and 5.5 % at the target value 188.5 mg/dL.

Comparison

A comparison between XL-640 automatic system LDL DIRECT (y) and a commercially available test (x) using 140 samples gave following results:

Linear regression: y = 0.920x + 6.591 mg/dL r = 0.986

Passing-Bablok¹⁷: y = 0.922x + 6.161 mg/dL r = 0.988

Interferences

Criterion: Recovery within ±10 % of initial value of LDL cholesterol concentration in the sample without interfering substance. Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 650 mg/dL, N-acetylcysteine, Metamizole and Acetaminofen (Paracetamol) including its metabolite N-acetyl-p-benzoquinone imine may cause false-negative results^{18,19}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of LDL cholesterol. Some drugs can also interfere. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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

- Pisani T, Gebiski CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
- Barham D, Trinder P, An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
- Wiebe DA, Warnick GR. Measurement of high-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 127-144, 1997.
- Nauck M, Maerz W, Wieland H. New immunoseparationbased homogenous assay for HDL-cholesterol compared with three homogenous and two heterogeneous methods for HDLcholesterol. Clin Chem 44: 1443-1451, 1998.
- Sidhu D, Naugler C. Fasting time and lipid levels in a community-based population: A cross sectional study. Arch. Intern. Med. Dec 10, 172(22): 1707-1710, 2012.
- Ontario Community Laboratory Guideline for Adult Lipid Testing (CLP017) 2013.
- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
- Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
- Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
- Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
- Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
- National Institutes of Health publication No. 93-3095, September, 1993.
- Williams P, et al., High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
- Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
- Anonymous, Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), 285, 2486-2497, JAMA 2001.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
- Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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LITERATÚRA

1. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Barham D, Trinder P, An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
3. Wiebe DA, Warnick GR. Measurement of high-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 127-144, 1997.
4. Nauck M, Maerz W, Wieland H. New immunoseparationbased homogenous assay for HDL-cholesterol compared with three homogenous and two heterogeneous methods for HDLcholesterol. Clin Chem 44: 1443-1451, 1998.
5. Sidhu D, Naugler C. Fasting time and lipid levels in a community-based population: A cross sectional study. Arch. Intern. Med. Dec 10, 172(22): 1707-1710, 2012.
6. Ontario Community Laboratory Guideline for Adult Lipid Testing (CLP017) 2013.
7. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
8. Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
9. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
10. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62: 707, 1977.
11. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
12. National Institutes of Health publication No. 93-3095, September, 1993.
13. Williams P, et al., High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
14. Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
15. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
16. Anonymous, Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), 285, 2486-2497, JAMA 2001.
17. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
18. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
19. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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