
IRON

Cat. No.	Pack Name	Packaging (Content)
XSYS0049	FE 125	R1: 4 × 25 mL, R2: 4 × 6.5 mL, R3 standard: 1 × 4 mL, RFID tag, instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of iron in human serum and plasma on automatic systems ERBA XL. Intended for screening, monitoring and diagnosis of iron deficiency anemia, hemochromatosis. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Ingested iron is mainly absorbed in the form of Fe²⁺ in the duodenum and upper jejunum. The trivalent form and the heme-bound Fe³⁺ component of iron in food has to be reduced by vitamin C. About 1 mg of iron is assimilated daily. Upon reaching the mucosal cells, Fe²⁺ ions become bound to transport substances. Before passing into the plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bound to transferrin in this form. The transport of Fe ions in blood plasma takes place via transferrin-iron complexes. A maximum of 2 Fe³⁺ ions per protein molecule can be transported. Serum iron is almost completely bound to transferrin. Iron (non-heme) measurements are used in the diagnosis and treatment of diseases such as iron deficiency anemia, hemochromatosis (a disease associated with widespread deposit in the tissue of the two iron-containing pigments, hemosiderin and hemofuscin, and characterized by pigmentation of the skin), and chronic renal disease. Iron determinations are performed for the diagnosis and monitoring of microcytic anemia (e.g. due to iron metabolism disorders and hemoglobinopathy), macrocytic anemia (e.g. due to vitamin B12 deficiency, folic acid deficiency and drug-induced metabolic disorders of unknown origin) as well as normocytic anemias such as renal anemia (erythropoietin deficiency), hemolytic anemia, hemoglobinopathy, bone marrow disease and toxic bone marrow damage.

PRINCIPLE

Transferrin-bound iron is released at an acidic pH and reduced from ferric to ferrous ions by ascorbic acid and thiourea. These ions react with Ferene S to form a violet coloured complex^{1,2,3,4}.

The resulting change in absorbance measured at 600 nm is proportional to the iron concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1
Acetate buffer (pH 4.5) 2.2 mol/L
Thiourea 180 mmol/L

R2
Thiourea 120 mmol/L
Ascorbic acid 240 mmol/L
Ferene S 3.0 mmol/L

R3 standard
Iron 500 µg/dL (89.5 µmol/L)

COMPOSITION OF REACTION MIXTURE
Acetate buffer 1.6 mol/L
Thiourea 150 mmol/L
Ascorbic acid 42.9 mmol/L
Ferene S 0.5 mmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.
Serum.

Plasma: Li-heparin plasma. Do not use EDTA or oxalate plasma. Separate serum or plasma from the clot or cells within 1 hour.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay. See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma^{5,6}:

7 days	at 15–25 °C
3 weeks	at 2–8 °C
1 year	at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with R3 standard is recommended. 2 point calibration (blank and calibrator); distilled water is recommended as blank. Calibration frequency: 14 days
Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, R3 standard and controls ERBA NORM and ERBA PATH have been standardized against to AAS standard.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	2-Point
Curve type	Linear
Wavelength (prim. / sec.)	600 / 700 nm
Reading time 1	just before adding of R2
Reading time 2	10 min. after adding of R1
Reaction direction	Increase
Unit	µg/dL (µmol/L)
Reagent volumes	
R1	180 µL
R2	45 µL
Sample volumes	27 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

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Iron Saturation:

20–55 %

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems is available on www.erba.com.

Limit of quantification:

1.5 µg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity:

1000 µg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (µg/dL)	SD (µg/dL)	CV (%)
Sample 1	123	1.41	1.15
Sample 2	217	1.75	0.81

Intermediate precision	Mean (µg/dL)	SD (µg/dL)	CV (%)
Sample 1	127	3.70	2.91
Sample 2	222	4.72	2.13

Accuracy

Two different validated control materials were used. Determined bias is -1.6 % at the target value 106 µg/dL and 0.5 % at the target value 188 µg/dL.

Comparison

A comparison between XL-640 automatic system IRON (y) and a commercially available test (x) using 101 samples gave following results:

Linear regression:
y = 1.004x - 4.586 µg/dL r = 0.996
Passing-Bablok¹¹:
y = 1.011x - 4.638 µg/dL r = 0.996

Interferences

Criterion: Recovery within ±10 % of initial value of iron concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin interferes (haemoglobin may lead to artificially increased values due to contamination of the sample with haemoglobin-bound iron), bilirubin min. up to 50 mg/dL, triglycerides up to 850 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels¹².

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of iron. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: 072Y-EJ4C-YE7J-R5TS

Warning

Contains: thiourea

Hazard statements:
H351 Suspected of causing cancer.

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection.
P308+P313 IF exposed or concerned: Get medical advice/attention.

R2

The reagent is not classified as dangerous.

R3

UFI: USJS-S9VF-UM7A-TRHT

Warning

Contains: hydroxylammonium chloride

Hazard statements:
H317 May cause an allergic skin reaction.
H351 Suspected of causing cancer.

Precautionary statements

P261 Avoid breathing vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection.
P302+P352 IF ON SKIN: Wash with plenty of water.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Stookey LL. Ferrozine-A new spectrophotometric reagent for iron. *Anal Chem* 42: 779-781, 1970.
2. Itano M. Serum Iron Survey. *Am J Clin Pathol* 70: 516-522, 1978.
3. Artiss JD, Vinogradov S, Zak B. Spectrophotometric study of several sensitive reagents for serum iron. *Clin Biochem* 14: 311-315, 1981.
4. Higgins T. Novel chromogen for serum iron determinations. *Clin Chem* 27(9):1619, 1981.
5. Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT-Verlag 1996.
6. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
7. Weissman, N., Pileggi, V.J., in *Clinical Chemistry: Principles and Technics*, 2nd Ed., R. J. Henry et al, editors, Hagerstown (MD), Harper & Row, 692-693, 1974.
8. Wick M. Iron metabolism and its disorders. In: Thomas L, editor. *Clinical laboratory diagnostics*. 1st ed. Frankfurt: THBooks Verlagsgesellschaft: 268-273, 1998.
9. Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company: 1642-1710, 1999.
10. Tietz, N. W.: *Textbook Of Clin. Chem.*, 3rd edition, W. B. Saunders Co., Philadelphia, 1999.
11. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem*, Nov;26(11): 783-790, 1988.
12. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001; 38:376-385.

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