# UIBC

Cat. No.	Pack Name	Packaging (Content)
XSYS0050	UIBC 125	R1: 4 x 25 mL, R2: 4 x 6.5 mL, R3 standard: 1 x 4 mL, RFID tag, instruction for use

## INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of unsaturated iron-binding capacity in human serum and plasma on automatic systems ERBA XL. Intended for screening, monitoring and diagnosis of iron deficiency anemia, hemochromatosis. For professional use in clinical laboratory only.

## CLINICAL SIGNIFICANCE

The prosthetic group of haemoglobin is the iron complex of protoporphyrin IX (heme) in which the centrally located iron atom acts as a stabilizer of oxyhemoglobin. Numerous enzymes and coenzymes require iron, e.g. peroxidases, catalases, cytochromes (which are also heme proteins), many of the enzymes of the Krebs cycle, and monoamine oxidase (which is involved in neurotransmission).


The total iron content of the body is about 3 to 3.5 g. Of this amount about 2.5 g is contained in erythrocytes or their precursors in the bone marrow.

Plasma contains only about 2.5 mg of iron. Iron is transported as Fe(III) bound to the plasma protein apotransferrin. The apotransferrin-Fe(III) complex is called transferrin. Normally only about one third of the iron-binding sites of transferrin are occupied by Fe(III). The additional amount of iron that can be bound is the unsaturated (or latent) iron-binding capacity (UIBC). The sum of the serum iron and UIBC represents total iron-binding capacity (TIBC). TIBC is a measurement for the maximum iron concentration that transferrin can bind.

The serum TIBC varies in disorders of iron metabolism. In iron-deficiency anemia the TIBC is elevated and the transferrin saturation is lowered to 15 % or less. Low serum iron associated with low TIBC is characteristic of the anemia of chronic disorders, malignant tumors, and infections.

## PRINCIPLE

Direct, colorimetric method with Ferrozine. In an alkaline environment known ferrous iron concentration incubated with serum, binds specifically with transferrin at unsaturated iron binding sites. Remaining unbound ferrous ions are measured with the chromogen reaction. The difference between the amount of excess iron and the total amount added to the serum is equivalent to the quantity bound to transferrin. This is the UIBC of the sample. The TIBC is equal to the serum iron concentration plus the UIBC<sup>1,2,3,4</sup>.

The resulting change in absorbance measured at 570 or 578 nm is directly proportional to the unbound excess iron concentration and indirectly proportional to the unsaturated iron-binding capacity in the sample.

## REAGENT DESCRIPTION AND COMPOSITION

R1	R2	R3 standard
Tris buffer (pH 8.45)	220 mmol/L Hydroxylamine hydrochloride	220 mmol/L
Ferrous ammonium sulfate	12.1 μmol/L Ferrozine	3.4 mmol/L
Hydroxylamine hydrochloride	100 mmol/L	
Iron	500 μg/dL (89.5 μmol/L)	

## COMPOSITION OF REACTION MIXTURE

Tris buffer	166 mmol/L
Ferrous ammonium sulfate	9.1 μmol/L
Hydroxylamine hydrochloride	117 mmol/L
Ferrozine	0.6 mmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

## MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

ERBA NORM 4x5, Cat. No. BLT00080  
 ERBA NORM 10x5, Cat. No. XSYS0123  
 ERBA PATH 4x5, Cat. No. BLT00081  
 ERBA PATH 10x5, Cat. No. XSYS0124  
 Erba XL analysers: XL-200, Cat. No. INS00002  
 XL-640, Cat. No. INS00008  
 XL-1000, Cat. No. INS00010

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.  
 Plasma: Li-heparin plasma. Li-heparin plasma values are approximately 6 % lower than serum values. The binding of iron to transferrin is strongly influenced by anions, specifically bicarbonate<sup>5</sup>. In order to avoid drift of UIBC recovery in samples over time the environmental concentration of CO<sub>2</sub> in the laboratory should be kept as constant as possible.

Specimens should be collected in the morning to avoid low results due to diurnal variation<sup>6</sup>. The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay. See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma <sup>7,8</sup> :	4 days at	15–25 °C
	7 days at	2–8 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with R3 standard is recommended. 2 point calibration (blank and calibrator); distilled water is recommended as blank. Calibration frequency: 7 days. Calibration is needed:  
 • after reagent lot change  
 • as required by internal quality control procedures  
 • calibration interval may be extended based on acceptable verification of calibration by the laboratory

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, R3 standard and controls ERBA NORM and ERBA PATH have been standardized against to AAS standard.

## ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

## ASSAY PARAMETERS FOR ERBA XL AUTOMATIC SYSTEMS

Assay type	2-Point
Curve type	Linear
Wavelength (prim. / sec.)	570 (578) / 700 nm
Reading time 1	just before adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	μg/dL (μmol/L)
Reagent volumes	
R1	200 μL
R2	50 μL
Sample volumes	15 μL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

## UNIT CONVERSION

μg/dL × 0.179 = μmol/L

## EXPECTED VALUES<sup>12,13</sup>

### UIBC:

Woman: 80–375 μg/dL  
 Man: 75–360 μg/dL

### TIBC:

Adult: 250–425 μg/dL  
 Infant: 100–400 μg/dL

TIBC = Iron + UIBC

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

**Limit of quantification:** 10.0 μg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

**Linearity:** 830 μg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

### Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (μg/dL)	SD (μg/dL)	CV (%)
Sample 1	153	3.4	2.22
Sample 2	216	4.7	2.17

Intermediate precision	Mean (μg/dL)	SD (μg/dL)	CV (%)
Sample 1	155	5.8	3.76
Sample 2	214	6.8	3.20

## Accuracy

Two different validated control materials were used. Determined bias is 16.8 % at the target value 55.7 μg/dL and 4.2 % at the target value 131 μg/dL.

## Comparison

A comparison between XL-640 automatic system UIBC (y) and a commercially available test (x) using 119 samples gave following results:

Linear regression:  
 $y = 1.045x + 10.827 \mu\text{g/dL}$   $r = 0.993$

Passing-Bablok<sup>14</sup>:  
 $y = 1.048x + 10.189 \mu\text{g/dL}$   $r = 0.991$

## Interferences

Criterion: Recovery within ±10 % of initial value of UIBC concentration in the sample without interfering substance. Following substances do not interfere: haemoglobin interferes (higher haemoglobin concentrations lead to artificially decreased values due to contamination of the sample with haemoglobin-bound iron), bilirubin min. up to 50 mg/dL, triglycerides min. up to 2600 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels<sup>15</sup>.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of UIBC. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### R1

UF1: M6RS-P9CG-4M7H-G87V

### Warning

Contains: thiourea  
**Hazard statements:**  
 H351 Suspected of causing cancer.

**Precautionary statements:**  
 P280 Wear protective gloves/protective clothing/eye protection.  
 P308 + P313 IF exposed or concerned: Get medical advice/attention.

### Supplemental information:

EUH208 Contains hydroxylamine hydrochloride. May produce an allergic reaction.

### R2

UF1: 9V6Y-RJDT-HE78-FM98

### Warning

Contains: hydroxylamine hydrochloride  
**Hazard statements:**

H317 May cause an allergic skin reaction.  
 H351 Suspected of causing cancer.

**Precautionary statements:**  
 P261 Avoid breathing mist/vapours/spray.  
 P280 Wear protective gloves/protective clothing/eye protection.  
 P308 + P313 IF exposed or concerned: Get medical advice/attention.  
 P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

### R3

UF1: USJS-S9VF-UM7A-TRHT

### Warning

Contains: hydroxylamine hydrochloride  
**Hazard statements:**

H317 May cause an allergic skin reaction.  
 H351 Suspected of causing cancer.

**Precautionary statements:**  
 P261 Avoid breathing vapours/spray.  
 P280 Wear protective gloves/protective clothing/eye protection.  
 P302 + P352 IF ON SKIN: Wash with plenty of water.  
 P308 + P313 IF exposed or concerned: Get medical advice/attention.  
 P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Stookey LL. Ferrozine-A new spectrophotometric reagent for iron. *Anal Chem* 42: 779-781, 1970.
2. Persijn JP, Van der Slik W, Riethorst A. Determination of serum iron and latent iron-binding capacity (LIBC). *Clin Chim Acta* 35: 91-98, 1971.
3. Artiss JD, Vinogradov S, Zak B. Spectrophotometric study of several sensitive reagents for serum iron. *Clin Biochem* 14: 311-315, 1981.
4. Higgins T. Novel chromogen for serum iron determinations. *Clin Chem* 27(9):1619, 1981.
5. Harris WR. Thermodynamics of Anion Binding to Human Serum Transferrin. *Biochemistry* 24: 7412-7418, 1985.
6. Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Tietz NW, ed. *Fundamentals of Clinical Chemistry*. 3rd ed. Philadelphia: WB Saunders: 789-824, 1987.
7. Weissman N, Pileggi VJ. Inorganic ions. In: Henry RJ, Cannon DC, Winkelman W, eds. *Clinical Chemistry, Principles and Techniques*. 2nd ed. Hagerstown: Harper & Row 4: 639-754, 1974.
8. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
9. Fairbanks VF, Klee GG. Biochemical aspects of hematology. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company: 1642-1710, 1999.
10. Wick M, Pingerra W, Lehmann P. Clinical aspects and laboratory. Iron metabolism, anemias. 5th ed. Wien, New York: Springer; 2003.
11. Bauer JD. Hemoglobin, porphyrin, and iron metabolism. In: Kaplan LA, Pesce AJ, eds. *Clinical Chemistry, theory, analysis, and correlation*. St. Louis: Mosby Company: 611-655, 1984.
12. Tietz, N. W.: *Textbook of Clin. Chem.*, 3rd edition, W. B. Saunders Co., Philadelphia, 1999.
13. Burtis CA, Bruns DE. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, 306, 2014.
14. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem*, Nov;26(11): 783-790, 1988.
15. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385.

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