CRP CALIBRATOR

Cat. No.	Pack Name	Packaging (Content)
XSYS0053	CRP CAL SH	1 × 1 mL, instruction for use

INTENDED USE

CRP CALIBRATOR is calibration material for *in vitro* quantitative determination of C-reactive protein in human serum on automatic systems ERBA XL by turbidimetric immunoassay. For professional use in clinical laboratory only.

PRINCIPLE

CRP CAL SH is a liquid calibrator based on human plasma and pleural fluid. The concentration of CRP is suitable for calibrations of clinical chemistry assays on automated clinical chemistry analysers Erba XL.

DESCRIPTION AND COMPOSITION

The calibrator is diluted human plasma and pleural fluid containing high levels of CRP with phosphate buffered saline. Liquid stabilized.

CERTIFIED VALUES

Calibrator Lot	Expiry	Parameter	Concentration	
XXXXXXXXXX	YYYY-MM	CRP	00.00 mg/dL	000.0 mg/L

UNIT CONVERSION

mg/dL × 10 = mg/L

CALIBRATOR PREPARATION

Calibrator is liquid, ready to use. For calibration generate a 6-point calibration curve by diluting the calibrator 1:32, 1:16, 1:8, 1:4, 1:2 and undiluted in saline.

STORAGE AND HANDLING

The unopened CRP CAL SH stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

After first opening the vial, calibrator can be used for 6 weeks if stored tightly closed at 2–8 °C after use. Do not freeze.

MATERIALS REQUIRED, BUT NOT PROVIDED

CRP, Cat. No. XSYS0047
CRP CON L, Cat. No. BLT20014
CRP CON H, Cat. No. BLT20013
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

ASSAY PROCEDURE

Use calibrator as specified in the relevant Instruction For Use for the CRP.

EVALUATION

The XL-analysers perform the evaluation automatically.

EXPECTED RESULTS

The calibrator value was determined using the immunoturbidimetric method. Determinations were performed under strictly standardized conditions on XL analysers using Erba system reagents and the standard reference materials.

The calibrator values were obtained via single determinations performed in different laboratories in several independent series. The calibrator value specified is the mean of all values obtained. Calibrator was tested within the functional testing of all diagnostic kits.

TRACEABILITY

This calibrator has been standardized against the reference material ERM-DA474/IFCC.

PERFORMANCE CHARACTERISTICS

Performance characteristics of individual analytes on the Erba XL automatic systems are available on websites www.erba.com.

LIMITATIONS:

Erroneous results can occur from inaccuracies at generate of calibration curve and technical errors associated with the assay procedure.

Improper storage or handling of the calibrator can also affect the results. Do not use the calibrator if there is visible evidence of microbial growth in the vial.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Health protection

Serum used in the manufacture of this product has been tested and found non-reactive for hepatitis B surface antigen (HbsAg), antibody to Hepatitis C (HCV), HIV-1 p24 antigen, and antibody to HIV-1/2. The test procedures do not guarantee that all infectious agents will be detected. Because no test method can offer complete assurance that Hepatitis B virus Hepatitis C virus and HIV-1/2 or other infectious agents are absent, the material should be handled as potentially infectious. During the work all principles of general safety as well as the principles of work with biological materials are to be maintained. Above all observe the principles of personal hygiene, use personal protection aids (clothes, shoes, etc.), do not eat, drink and smoke. After finishing the work and before a meal wash your hands with soap and warm water and apply a suitable regenerating cream. When the bottle is broken, the contaminated surface must be disinfected with Persteril or Chloramine.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The calibrator is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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