# LIPASE

Cat. No.	Pack Name	Packaging (Content)
XSYS0081	LIP 110	R1: 2 × 44 mL, R2: 2 × 11 mL, RFID tag, instruction for use

## INTENDED USE



The kit is intended for *in vitro* photometric quantitative determination of lipase in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of pancreatic disorders such as pancreatitis and obstruction of pancreatic duct. For professional use in clinical laboratories only.

## CLINICAL SIGNIFICANCE

Lipases (LIP) are enzymes which hydrolyze glycerol ester of long fatty acids. The enzyme and its co-factor colipase is produced in the pancreas, lipase being also secreted in small amounts by the salivary glands as well as by gastric, pulmonary and intestinal mucosa. Bile acids and colipase form micellar complexes with the lipids and bind lipase on the substrate / water interface. Determination of lipase is used for investigation of pancreatic disorders. In acute pancreatitis the lipase concentrations rise to 2–50 fold to upper reference limit within 4–8 hours after begin of abdominal pain peaking at 24 hours and decreasing within 8 to 14 days. Elevated lipase values can also be observed in chronic pancreatitis and obstruction of the pancreatic duct.

## PRINCIPLE

Enzymatic colorimetric assay with 1,2-O-dilauryl-rac-glycerol-3-glutaric-acid-(6-methylresorufin) ester as substrate<sup>1,2,3,4</sup>. The chromogenic lipase substrate 1,2-O-dilauryl-rac-glycerol-3-glutaric-acid-(6-methylresorufin) ester is cleaved by the catalytic action of alkaline lipase solution to form 1,2-O-dilauryl-rac-glycerol and an unstable intermediate, glutaric acid-(6-methylresorufin) ester. This decomposes spontaneously in alkaline solution to form glutaric acid and methylresorufin. Addition of detergent and colipase increases the specificity of the assay for pancreatic lipase.

The rate of colour formation is proportional to the activity of LIP present in the sample and can be measured kinetically at 570 or 578 nm.

## REAGENT DESCRIPTION AND COMPOSITION

R1	R2
Good's Buffer pH 8.0	Tartrate Buffer pH 4.0
Taurodesoxycholate	Lipase Color Substrate
Desoxycholate	
Calcium ions	
Colipase	

## COMPOSITION OF REACTION MIXTURE

Good's / Tartrate Buffer	
Taurodesoxycholate	≥0.8 mmol/L
Desoxycholate	≥0.8 mmol/L
Calcium ions	≥0.8 mmol/L
Colipase	≥0.8 mg/L
Lipase Color Substrate	≥0.02 mmol/L

## REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

## MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034  
 XL MULTICAL 10×3, Cat. No. XSYS0122  
 ERBA NORM 4×5, Cat. No. BLT00080  
 ERBA NORM 10×5, Cat. No. XSYS0123  
 ERBA PATH 4×5, Cat. No. BLT00081  
 ERBA PATH 10×5, Cat. No. XSYS0124  
 Erba XL analysers: XL-200, Cat. No. INS00002  
 XL-640, Cat. No. INS00008  
 XL-1000, Cat. No. INS00010

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

## SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma <sup>5</sup> :	7 days at	20–25 °C
	3 weeks at	4–8 °C
	12 months at	-20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 8 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

## QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

## TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to enzymatic colorimetric assay<sup>1,2,3,4</sup>.

## ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see [www.erba.com](http://www.erba.com).

## Assay parameters for ERBA XL automatic systems

Assay type	Rate A
Curve type	Linear
Wavelength (prim. / sec.)	570 (578) / 660 nm
Reading time	90–195 s after adding of R2
Reaction direction	Increase
Unit	U/L (µkat/L)
Reagent volumes	
R1	160 µL
R2	40 µL
Sample	3.3 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

## UNIT CONVERSION

U/L × 0.0167 = µkat/L

## EXPECTED VALUES<sup>7</sup>

At 37 °C ≤60 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

## ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on [www.erba.com](http://www.erba.com).

**Limit of quantification:** 6.88 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

**Linearity:** 678 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

## Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	48.2	0.36	0.75
Sample 2	86.0	0.54	0.63

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	47.3	1.17	2.48
Sample 2	85.0	1.50	1.77

## Accuracy

Two different validated control materials were used. Determined bias is -7.4 % at the target value 91.2 U/L and -12.3 % at the target value 122.4 U/L.

## Comparison

A comparison between XL-640 automatic system LIP (y) and a commercially available test (x) using 53 samples gave following results:

Linear regression:

$$y = 1.017x - 2.76 \text{ U/L} \quad r = 0.981$$

Passing-Bablok<sup>6</sup>:

$$y = 1.071x - 4.03 \text{ U/L} \quad r = 0.976$$

## Interferences

Criterion: Recovery within ±10 % of initial value of LIP activity in the sample without interfering substance. Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

## Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of LIP. See paragraph Interferences.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

### R1

The reagent is not classified as dangerous.

### R2

UFI: MKEV-5WC0-KJ50-42R5

Danger

Contains: propan-1-ol

## Hazard statement:

H318 Causes serious eye damage.

## Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Neumann U, Junius M, Batz HG, et al. New substrates for the optical determination of lipase. EP 207252 1987.
2. Borgström B. The action of bile salts and other detergents on pancreatic lipase and the interaction with colipase. Biochimica et Biophysica Acta 1977;488:381-391.
3. Gargouri Y, Julien R, Bois A, et al. Studies on the detergent inhibition of pancreatic lipase activity. J of Lipid Research 1983;24:1336-1342.
4. Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv Clin Enzymol 1986;4:60-67.
5. Lorentz K Lipase. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH- Books Verlagsgesellschaft;1998. p. 95-7.
6. Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 689-708.
7. Tietz N, Shuey DF. Lipase in serum – the elusive enzyme: an overview. Clin Chem 1993; 39:746- 56.
8. Lott J, Patel ST, Sawhney AK, Kazmierczak SC, LoveJE. Assays of serum lipase: analytical and clinical considerations. Clin Chem 1986;32:1290-1302.
9. Guder W, Fonseca-Wollheim W, et al, Quality of Diagnostic Samples, German Society for Clinical Chemistry and Laboratory Medicine, 3rd completely revised edition 2010.
10. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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## LITERATÚRA

1. Neumann U, Junius M, Batz HG, et al. New substrates for the optical determination of lipase. EP 207252 1987.
2. Borgström B. The action of bile salts and other detergents on pancreatic lipase and the interaction with colipase. *Biochimica et Biophysica Acta* 1977;488:381-391.
3. Gargouri Y, Julien R, Bois A, et al. Studies on the detergent inhibition of pancreatic lipase activity. *J of Lipid Research* 1983;24:1336-1342.
4. Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. *Adv Clin Enzymol* 1986;4:60-67.
5. Lorentz K Lipase. In: Thomas L, editor. *Clinical laboratory diagnostics*. 1st ed. Frankfurt: TH- Books Verlagsgesellschaft;1998. p. 95-7.
6. Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 689-708.
7. Tietz N, Shuey DF. Lipase in serum – the elusive enzyme: an overview. *Clin Chem* 1993; 39:746- 56.
8. Lott J, Patel ST, Sawhney AK, Kazmierczak SC, LoveJE. Assays of serum lipase: analytical and clinical considerations. *Clin Chem* 1986;32:1290-1302.
9. Guder W, Fonseca-Wollheim W, et al, *Quality of Diagnostic Samples*, German Society for Clinical Chemistry and Laboratory Medicine, 3rd completely revised edition 2010.
10. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

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**REF**

Katalógové číslo

**LOT**

Číslo šarže

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eIFU:  
[www.erba.com](http://www.erba.com)**IVD**Diagnostický zdravotnícky prostriedok *in vitro*

Výrobca

Obmedzenie teploty

**CONT**

Obsah

QUALITY SYSTEM CERTIFIED  
ISO 13485, IVDRErba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: [diagnostics@erba.com](mailto:diagnostics@erba.com), [www.erba.com](http://www.erba.com)

CC/IFU/021/26/A

Dátum revízie: 3. 3. 2026