GLUCOSE HK

Cat. No.	Pack Name	Packaging (Content)
XSYS0095	GLU HK 330	R1: 6 × 44 mL, R2: 6 × 11 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of glucose in human serum, plasma and urine on automatic systems ERBA XL. Intended for screening, monitoring and diagnosis of diabetes, hypoglycemia. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Glucose is the major carbohydrate present in the peripheral blood. Oxidation of glucose is the major source of cellular energy in the body. Glucose derived from dietary sources is converted to glycogen for storage in the liver or to fatty acids for storage in adipose tissue. The concentration of glucose in blood is controlled within narrow limits by many hormones, the most important of which are produced by the pancreas.


Accurate measurement of glucose in body fluid is important in diagnosis and management of diabetes, hypoglycemia, adrenal dysfunction and various other conditions. High levels of serum glucose may be seen in case of Diabetes mellitus, in patients receiving glucose containing fluids intravenously, during severe stress and in cerebrovascular accidents.

Decreased levels of glucose can be due to insulin administration, as a result of insulinoma, inborn errors of carbohydrate metabolism or fasting.

Glucose measurement in urine is used as a diabetes screening procedure and to aid in the evaluation of glycosuria, to detect renal tubular defects, and in the management of diabetes mellitus.

PRINCIPLE

Hexokinase (HK) catalyzes the phosphorylation of glucose to glucose-6-phosphate by ATP. Glucose-6-phosphate dehydrogenase (G-6-PDH) oxidizes glucose-6-phosphate in the presence of NAD⁺ to gluconate-6-phosphate. No other carbohydrate is oxidized^{1,2}.

The absorbance of NADH measured at 340 nm is directly proportional to the glucose concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Tris buffer, pH 7.40	100 mmol/L	MgCl ₂	25 mmol/L
ATP	2.5 mmol/L	Hexokinase	>10 kU/L
NAD ⁺	2.5 mmol/L	G-6-PDH	>10 kU/L

COMPOSITION OF REACTION MIXTURE

Tris buffer	79 mmol/L
ATP	2.0 mmol/L
NAD ⁺	2.0 mmol/L
MgCl ₂	5.0 mmol/L
Hexokinase	>2 kU/L
G-6-PDH	>2 kU/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 60 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and EDTA plasma.

Separate at the latest 1h after blood collection from cellular contents.

Urine.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in plasma after addition of a glycolytic inhibitor (Fluoride, monoiodacetate, mannose)³:

2 days at	20–25 °C
7 days at	4–8 °C
1 day at	-20 °C

Stability in serum (separated from cellular contents, hemolysis free) without adding a glycolytic inhibitor^{4,5}:

8 hours at	25 °C
72 hours at	4 °C

Stability in urine⁴

Discard contaminated specimens.

24 hours at	4–8 °C
-------------	--------

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits

should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	2-Point
Curve type	Linear
Wavelength (prim. / sec.)	340 / 405 (415) nm
Reading time 1	just before adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (mmol/L)

Reagent volumes

R1	180 µL (160 µL / XL-1000)
R2	45 µL (40 µL / XL-1000)

Sample

2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION

mg/dL × 0.056 = mmol/L

EXPECTED VALUES^{6,7}

At 37 °C

Glucose Fasting:

Cord:	45–96 mg/dL
Newborn, 1 d:	40–60 mg/dL
Newborn, >1 d:	50–80 mg/dL
Child:	60–100 mg/dL
Adult:	74–100 mg/dL
>60 y:	82–115 mg/dL
>90 y:	75–121 mg/dL

Random Glucose:

≥1 y: 70–140 mg/dL

Glucose 2 h Postprandial: <120 mg/dL

WB (Hep) Adult: 65–95 mg/dL

Urine: 1–15 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 1.47 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 1000 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	79.7	0.33	0.42
Sample 2	237.6	1.59	0.67

Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	115.8	1.04	0.90
Sample 2	346.6	1.28	0.37

Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	79.9	1.67	2.09
Sample 2	241.9	3.54	1.46

Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	105.4	1.48	1.41
Sample 2	312.4	8.83	2.83

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is -0.7 % at the target value 69.5 mg/dL, -1.4 % at the target value 211.4 mg/dL for serum, 3.7 % at the target value 27.7 mg/dL and 1.1 % at the target value 294.6 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system GLUCOSE HK (y) and a commercially available test (x) using 146 samples (serum) gave following results:

Linear regression:

$$y = 1.026x - 6.545 \text{ mg/dL} \quad r = 0.990$$

Passing-Bablok⁸:

$$y = 1.029x - 6.901 \text{ mg/dL} \quad r = 0.982$$

Interferences

Criterion: Recovery within ±10 % of initial value of glucose concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs:

Serum: No interference was found at therapeutic concentrations using common drug panels⁹.

Urine: No interference was found at therapeutic concentrations using common drug panels⁹.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.

- Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of glucose. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. The Summary of Safety and Performance is available at Eudamed see: <https://ec.europa.eu/tools/eudamed>.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Kunst A, Draeger B, Ziegenhorn J. In: Bergmeyer. Methods of Enzymatic Analysis, 3rd ed. Volume VI, Metabolites 1: Carbohydrates 1984;163-172.
2. Tietz NW, ed. Clinical Guide to Laboratory Tests, 4th ed. Philadelphia: WB Saunders Co 2006; 444-451.
3. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; 30-31.
4. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. 750-808.
5. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2002; 48: 436-72.
6. Tietz NW, (Ed.), Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
7. Tietz Textbook of Clinical Chemistry. 5th ed. Edited by Burtis CA, Ashwood ER, Bruns DE. WB Saunders Company, Philadelphia 2006, 25: 837-907.
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov; 26(11): 783-790.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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LITERATÚRA

1. Kunst A, Draeger B, Ziegenhorn J. In: Bergmeyer. Methods of Enzymatic Analysis, 3rd ed. Volume VI, Metabolites 1: Carbohydrates 1984;163-172.
2. Tietz NW, ed. Clinical Guide to Laboratory Tests, 4th ed. Philadelphia: WB Saunders Co 2006; 444-451.
3. Guder W, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; 30-31.
4. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. 750-808.
5. Sacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chem 2002; 48: 436-72.
6. Tietz NW, (Ed.), Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
7. Tietz Textbook of Clinical Chemistry. 5th ed. Edited by Burtis CA, Ashwood ER, Bruns DE. WB Saunders Company, Philadelphia 2006, 25: 837-907.
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov; 26(11): 783-790.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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