FERRITIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0101	FRTN	R1: 2 × 14.5 mL, R2: 2 × 7.7 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* quantitative determination of ferritin in human serum and plasma on automatic systems ERBA XL by turbidimetric immunoassay. Intended for screening, monitoring and diagnosis of iron deficiency of all population. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

The plasma Ferritin concentration declines very early in the development of iron deficiency. On the other hand, a large number of chronic diseases result in increased serum Ferritin concentrations. These diseases include chronic infections, chronic inflammatory disorders such as rheumatoid arthritis or renal disease, Gaucher's disease, and numerous types of malignancies, especially lymphomas, leukaemia, breast cancer and neuroblastoma. Increase in plasma Ferritin concentration also occurs in viral hepatitis or following toxic liver injury because of release of Ferritin from damaged liver cells. Plasma Ferritin concentration is also increased with increases of iron stores, as seen in patients with hemosiderosis or hemochromatosis. Besides the use of Ferritin as an iron metabolism parameter, Ferritin as also gained importance as a tumour marker for therapeutic drug monitoring and follow-up^{1,2,3}.

PRINCIPLE

Determination of the concentration of ferritin by photometric measurement of antigen-antibody reaction of antibodies to human ferritin bound to latex particles with ferritin present in the sample. The turbidity of immune complex, measured at 600 nm, is directly proportional to the ferritin sample concentration⁴.

REAGENT DESCRIPTION AND COMPOSITION

R1	
TRIS buffer, pH 7.2	<100 g/L
Proclin 300	<0.0025 %
Non reactant components	
R2	
Latex particles coated with rabbit F(ab) ₂	
antibodies against human Ferritin	<110 mg/L
Sodium azide	<0.1 %
2-chloroacetamide	<0.1 %
Non reactant components	

REAGENT PREPARATION

Reagents are liquid, ready to use. R2 needs to be mixed via inversion 5–10 times before being placed on the instrument and the stirring must be repeated on a weekly basis. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

FRTN CAL, Cat. No. XSYS0104	
FRTN CON L, Cat. No. BLT20042	
FRTN CON H, Cat. No. BLT20041	
Erba XL analysers:	
XL-200, Cat. No. INS00002	
XL-640, Cat. No. INS00008	
XL-1000, Cat. No. INS00010	

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: 60 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable⁵.

Serum, heparinised or EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum, plasma⁶:	7 days at	20–25 °C
	2 days at	4–8 °C
	1 year at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with FRTN CAL is recommended.

5-point calibration: generate a 5-point calibration curve by diluting the calibrator 1:16, 1:8, 1:4, 1:2

and undiluted in saline; saline is used as blank

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control FRTN CON L and FRTN CON H are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator FRTN CAL and controls FRTN CON L and FRTN CON H have been standardized against the reference material NIBSC 94/572 (WHO International Standard)⁶.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com

ASSAY PARAMETERS FOR ERBA XL AUTOMATIC SYSTEMS

Assay type	2-Point
Curve type	Cubic Spline
Wavelength (prim. / sec.)	600 / - nm
Reading time 1	40 s after adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	ng/mL (µg/L)

Reagent volumes

R1	120 µL
R2	60 µL
Sample volumes	20 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION

ng/mL = µg/L

EXPECTED VALUES²

Newborn	25–200 ng/mL
1 month	200–600 ng/mL
2–5 months	50–200 ng/mL
6 months–15 years	7–140 ng/mL
Adult	
Male	20–250 ng/mL
Female	10–120 ng/mL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com

Limit of quantification: 4.42 ng/mL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30)

Linearity: approx. 420 ng/mL

Linearity is given by concentration of the calibrator used for generating 5-point calibration curve.

Hook Effect: safety zone >5500 ng/mL

Precision

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (ng/mL)	SD (ng/mL)	CV (%)
Sample 1	55.2	2.13	3.87
Sample 2	260.0	11.36	4.37

Intermediate precision	Mean (ng/mL)	SD (ng/mL)	CV (%)
Sample 1	52.0	2.21	4.24
Sample 2	257.3	8.52	3.31

Accuracy

Two different validated control materials were used. Determined bias is and -5.1 % at the target value 160 ng/mL and 5.4 % at the target value 338 ng/mL.

Comparison

A comparison between XL-640 automatic system FRTN (y) and a commercially available test (x) using 150 samples gave following results:

Linear regression:

$$y = 1.024x + 1.454 \text{ ng/mL} \quad r = 0.994$$

Passing-Bablok⁷:

$$y = 1.061x - 0.534 \text{ ng/mL} \quad r = 0.996$$

Interferences

Criterion: Recovery within ±10 % of initial value of ferritin concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 13.3 g/L, bilirubin up to 35 mg/dL, triglycerides up to 1910 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁸.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of ferritin. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: CSDW-AW9R-DJ56-Y7UV

Warning

Contains: reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1)

Hazard statements:

H317 May cause an allergic skin reaction.

Precautionary statements:

P280 Wear protective gloves.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

R2

Supplemental information:

EUH 208 Contains 2-chloroacetamide. May produce an allergic reaction.

Health protection

Used animal-source materials are accompanied by an official state veterinarian health certificate. The animal components were collected from healthy animals under the care of a registered veterinarian and were apparently free from infectious and contagious diseases and injurious parasites. However, this material should be handled as thought capable of transmitting infectious diseases.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Tietz NW, ed., Clinical Guide to Laboratory Tests. 3rd. ed. Philadelphia: W.B. Saunders Company Ltd., 1995.
2. Burtis CA, Ashwood ER, Bruns DE, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 5th edition: W.B. Saunders Company Ltd., 2012.
3. Jacobs DJ, Demott WR, Grady HJ, Horvat RJ, Huestis DW, Kasten, eds., Laboratory Test Handbook. 4th. ed. Ohio, Hudson: Lexi-Comp Inc., 1996.
4. Bernard A, Lauwerys R, Turbidimetric latex immunoassay for serum ferritin, Journal of Immunological methods 71, 141-147, 1984.
5. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
6. Thorpe SJ, Walker D, Arosio P, Heath A, Cook JD, Worwood M, International collaborative study to evaluate a recombinant L ferritin preparation as an International Standard, Clin Chem 43 (9), 1582-1587, 1997.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
8. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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