
HDL DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00028	HDL 80	R1: 2 × 30 mL, R2: 2 × 10 mL, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of HDL cholesterol in human serum and plasma on various automatic systems. In combination with other parameters, it is intended for prognosis, prediction, and diagnosis of coronary heart disease of all population. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

High-density lipoproteins (HDL) compose one of the major classes of plasma lipoproteins. They are synthesized in liver as complexes of apolipoprotein and phospholipid and are capable of picking up cholesterol and carrying it from arteries to the liver, where the cholesterol is converted to bile acids and excreted into the intestine.

An inverse relationship between HDL Cholesterol (HDL C) levels in serum and the incidence/prevalence of coronary heart disease (CHD) has been demonstrated in a number of epidemiological studies. The importance of HDL C as a risk factor for CHD is now recognized.

Accurate measurement of HDL C is of vital importance when assessing patient's risk for CHD.

PRINCIPLE

The assay is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol-methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHO) and cholesterol esterase (CHE).

The enzymes selectively react with HDL to produce hydrogen peroxide (H₂O₂) which is detected through the Trinder's reaction^{1,2,3,4}.

The absorbance of the produced quinone dye at 600 nm is proportional to the HDL C concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
MES buffer (pH 6.5)	6.5 mmol/L
N-Bis(4-sulfobutyl)-3-methylaniline, TODB	3 mmol/L
Polyvinyl sulfonic acid	50 mg/L
Polyethylene-glycol-methyl ether	30 ml/L
Magnesium chloride	2 mmol/L

R2	
MES buffer (pH 6.5)	50 mmol/L
Cholesterol esterase	5 kU/L
Cholesterol oxidase	20 kU/L
Peroxidase	5 kU/L
4-aminoantipyrine, 4-AAP	0.9 g/L
Detergent	0.5 %

COMPOSITION OF REACTION MIXTURE

MES buffer (pH 6.5)	17.2 mmol/L
N-Bis(4-sulfobutyl)-3-methylaniline, TODB	2.0 mmol/L
4-aminoantipyrine, 4-AAP	0.2 g/L
Polyvinyl sulfonic acid	37 mg/L
Polyethylene-glycol-methyl ether	22 ml/L
Magnesium chloride	1.5 mmol/L
Cholesterol esterase	1.2 kU/L
Cholesterol oxidase	5.0 kU/L
Peroxidase	1.2 kU/L
Detergent	0.1 %

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 600/700 nm may be used, general laboratory equipment.

HDL/LDL Calibrator, Cat. No. XSYS0061
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagents are stable for 30 days at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.
 For specimen collection and preparation only use suitable tubes or collection containers.
 Only the specimens listed below were tested and found acceptable.

Serum.
 Plasma: Li-heparin plasma.
 Fasting and non-fasting samples can be used^{5,6}.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.
 See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ⁷ :	2 days at	20–25 °C
	7 days at	4–8 °C
	3 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with HDL/LDL Calibrator is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank
 Calibration frequency: it is recommended to do a calibration
 • after reagent lot change
 • as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM 4×5 or ERBA NORM 10×5 and ERBA PATH 4×5 or ERBA PATH 10×5 are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, HDL/LDL calibrator and controls ERBA NORM and ERBA PATH have been standardized against the reference material NIST SRM 1951.

ASSAY PROCEDURE

Wavelength: 600/700 nm
 Cuvette: 1 cm

	Reagent blank	Calibrator	Sample
Reagent 1	0.375 mL	0.375 mL	0.375 mL
Sample	–	–	0.005 mL
Calibrator	–	0.005 mL	–
Distilled water	0.005 mL	–	–

Mix and after 5 min. incubation read the initial absorbance for blank A_{bl} , sample A_{sam} and calibrator A_{cal} . Then add:

Reagent 2	0.125 mL	0.125 mL	0.125 mL
-----------	----------	----------	----------

Mix and after 5 min. incubation read the final absorbance for blank A_{bl} , sample A_{sam} and calibrator A_{cal} . Calculate resulting absorbance like the difference between the final and initial absorbance $A = (A_{\text{FINAL}} - A_{\text{INITIAL}})$.

CALCULATION

$$\text{HDL cholesterol (mg/dL)} = \frac{A_{\text{sam}} - A_{\text{bl}}}{A_{\text{cal}} - A_{\text{bl}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	1-Point End	Normal Low (mg/dL)	28
Wavelength 1 (nm)	600	Normal High (mg/dL)	96
Wavelength 2 (nm)	670	Linearity Low (mg/dL)	0.6
Sample Volume (µL)	5	Linearity High (mg/dL)	208
Reagent 1 Volume (µL)	375	Concentration of Standard	See bottle label
Reagent 2 Volume (µL)	125	Blank with	Reagent
Incubation time (min.)	5	Absorbance limit (max.)	0.3
Reaction temperature (°C)	37	Units	mg/dL
Reaction direction	Increasing		

UNIT CONVERSION

mg/dL × 0.026 = mmol/L

EXPECTED VALUES¹⁵

In serum:	Male	Female
5–9 y	38–75	36–73 mg/dL
10–14 y	37–74	37–70 mg/dL
15–19 y	30–63	35–74 mg/dL
20–24 y	30–63	33–79 mg/dL
25–29 y	31–63	37–83 mg/dL
30–34 y	28–63	36–77 mg/dL
35–39 y	29–62	34–82 mg/dL
40–44 y	27–67	34–88 mg/dL
45–49 y	30–64	34–87 mg/dL
50–54 y	28–63	37–92 mg/dL
55–59 y	28–71	37–91 mg/dL
60–64 y	30–74	38–92 mg/dL
65–69 y	30–75	35–96 mg/dL
>69 y	31–75	33–92 mg/dL

Adult Treatment Panel III classification¹⁶:

Low <40 mg/dL
 High >59 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.63 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 208 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	42.7	0.44	1.03	Sample 1	43.1	0.86	2.01
Sample 2	79.0	0.39	0.49	Sample 2	84.2	1.76	2.09

Accuracy

Two different validated control materials were used. Determined bias is -11.1 % at the target value 55.4 mg/dL and -14.1 % at the target value 83.8 mg/dL.

Comparison

A comparison between XL-640 automatic system HDL DIRECT (y) and a commercially available test (x) using 131 samples gave following results:

Linear regression:

$$y = 1.087x - 1.265 \text{ mg/dL} \quad r = 0.990$$

Passing-Bablok¹⁷:

$$y = 1.086x - 1.402 \text{ mg/dL} \quad r = 0.991$$

Interferences

Criterion: Recovery within ±10 % of initial value of cholesterol concentration in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 12.5 g/L, bilirubin up to 20 mg/dL, triglycerides up to 850 mg/dL.

N-acetylcysteine, Metamizole and Acetaminofen (Paracetamol) including its metabolite N-acetyl-p-benzoquinone imine may cause false-negative results^{18,19}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 600 nm against the distilled water is 0.3.
 - High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of HDL cholesterol. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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
1. Pisani T, Gebbski CP, Leary ET, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
3. Wiebe DA, Warnick GR. Measurement of high-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 127-144, 1997.
4. Nauck M, Maerz W, Wieland H. New immunoseparationbased homogenous assay for HDL-cholesterol compared with three homogenous and two heterogeneous methods for HDL-cholesterol. Clin Chem 44: 1443-1451, 1998.
5. Sidhu D, Naugler C. Fasting time and lipid levels in a community-based population: A cross sectional study. Arch. Intern. Med. Dec 10, 172(22): 1707-1710, 2012.
6. Ontario Community Laboratory Guideline for Adult Lipid Testing (CLP017) 2013.
7. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
8. Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
9. Barr DP, Russ EM, Eder HA. Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
10. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
11. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
12. National Institutes of Health publication No. 93-3095, September, 1993.
13. Williams P, et al., High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
14. Kannel WB, Castelli WP, Gordon T. Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
15. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
16. Anonymous, Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), 285, 2486-2497, JAMA 2001.
17. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
18. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
19. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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LITERATÚRA

1. Pisani T, Gebiski CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 119: 1127-1135, 1995.
2. Barham D, Trinder P. An improved color reagent for the determination of blood glucose by the oxidase system. Analyst 97: 142-145, 1972.
3. Wiebe DA, Warnick GR. Measurement of high-density lipoprotein cholesterol. In: Rifai N, Warnick GR, Dominiczak MH, eds. Handbook of lipoprotein testing. Washington: AACC Press: 127-144, 1997.
4. Nauck M, Maerz W, Wieland H. New immunoseparationbased homogenous assay for HDL-cholesterol compared with three homogenous and two heterogeneous methods for HDL-cholesterol. Clin Chem 44: 1443-1451, 1998.
5. Sidhu D, Naugler C. Fasting time and lipid levels in a community-based population: A cross sectional study. Arch. Intern. Med. Dec 10, 172(22): 1707-1710, 2012.
6. Ontario Community Laboratory Guideline for Adult Lipid Testing (CLP017) 2013.
7. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag: 22-23, 2001.
8. Dominiczak M, McNamara J. The system of Cardiovascular prevention. 103–125; Nauk M, Wiebe D, Warnick G. Measurement of High-Density-Lipoprotein Cholesterol. 221–244. In: Handbook of Lipoprotein Testing (eds. Rifai, Warnick and Dominiczak), 2nd edition.
9. Barr DP, Russ EM, Eder HA. Protein-lipid relationships in human plasma, Am. J. Med. 11: 480, 1951.
10. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med. 62; 707, 1977.
11. Castelli, WP et al., HDL Cholesterol and other lipids in coronary heart disease, Circulation, 55: 767, 1977.
12. National Institutes of Health publication No. 93-3095, September, 1993.
13. Williams P, et al., High density lipoprotein and coronary risk factor, Lancet, 1: 72, 1979.
14. Kannel WB, Castelli WP, Gordon T, Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med. 90: 85, 1979.
15. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
16. Anonymous, Executive summary of the third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), 285, 2486-2497, JAMA 2001.
17. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
18. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
19. Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski MD, N-acetylcysteine interference of Trinder-based assays, Clin Biochem 49, 100 -104, 2016.

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