CREATINE KINASE

Cat. No.	Pack Name	Packaging (Content)
XSYS0022	CK 110	R1: 2 x 44 mL, R2: 2 x 11 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of creatine kinase in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for monitoring and diagnosis of myocardial or skeletal muscle disease. For professional use in clinical laboratories only.


CLINICAL SIGNIFICANCE

Creatine kinase (CK) is a dimeric enzyme occurring in four different forms: a mitochondrial isoenzyme and the cytosolic isoenzymes CK-MM (muscle type), CK-BB (brain type) and CK-MB (myocardial type). The determination of CK and CK-isoenzyme activities is utilized in the diagnosis and monitoring of myocardial infarction and myopathies such as the progressive Duchenne muscular dystrophy. Following injury to the myocardium, as occurs with acute myocardial infarction, CK is released from the damaged myocardial cells. In early cases a rise in the CK activity can be found just 4 hours after an infarction, the CK-activities reaches a maximum after 12–24 hours and then falls back to the normal range after 3–4 days. Myocardial damage is very likely when the total CK activity is above 190 U/L, the CK-MB activity is above 24 U/L (37 °C) and the CK-MB activity fraction exceeds 6 % of the total.

PRINCIPLE

The assay method using creatine phosphate and ADP was first described by Oliver¹, modified by Rosalki² and further improved for optimal test conditions by Szasz et al³. CK is rapidly inactivated by oxidation of the sulfhydryl groups in the active center. The enzyme can be reactivated by the addition of N-acetylcysteine (NAC)³. Interference by adenylate kinase is prevented by the addition of diadenosine pentaphosphate⁴ and AMP^{3,4}.

Standardized methods for the determination of CK using the "reverse reaction" and activation by NAC were recommended by the German Society for Clinical Chemistry (DGKC)⁴ in 1977 and the International Federation of Clinical Chemistry (IFCC)⁵ in 1991. In 2002 the IFCC confirmed their recommendation and extended it to 37 °C^{6,7}. This assay meets the recommendations of the IFCC and DGKC.

The rate of formation of NADPH is proportional to the activity of CK present in the sample and can be measured kinetically at 340 nm.

REAGENT DESCRIPTION AND COMPOSITION

R1	R2		
Imidazole buffer, pH 6.1	Imidazole buffer, pH 8.9	125 mmol/L	125 mmol/L
Glucose	ADP	25 mmol/L	15.2 mmol/L
Magnesium acetate	D-glucose-6-phosphate-dehydrogenase	12.5 mmol/L	>8800 U/L
EDTA	Creatine phosphate	2 mmol/L	250 mmol/L
N-acetylcysteine	AMP	25 mmol/L	25 mmol/L
NADP ⁺	Diadenosine pentaphosphate	2.4 mmol/L	103 µmol/L
Hexokinase		>6800 U/L	

COMPOSITION OF REACTION MIXTURE

Imidazole buffer	Hexokinase	123 mmol/L	>5333 U/L
Glucose	ADP	19.6 mmol/L	3.0 mmol/L
Magnesium acetate	D-glucose-6-phosphate-dehydrogenase	9.8 mmol/L	>1725 U/L
EDTA	Creatine phosphate	1.6 mmol/L	49 mmol/L
N-acetylcysteine	AMP	19.6 mmol/L	4.9 mmol/L
NADP ⁺	Diadenosine pentaphosphate	1.9 mmol/L	20 µmol/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4x3, Cat. No. XSYS0034
 XL MULTICAL 10x3, Cat. No. XSYS0122
 ERBA NORM 4x5, Cat. No. BLT00080
 ERBA NORM 10x5, Cat. No. XSYS0123
 ERBA PATH 4x5, Cat. No. BLT00081
 ERBA PATH 10x5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
 On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.
 For specimen collection and preparation only use suitable tubes or collection containers.
 Only the specimens listed below were tested and found acceptable.
 Serum.

Plasma: Li-heparin and K₂-, K₃-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma ¹⁰ :		
	2 days at	15–25 °C
	7 days at	2–8 °C
	4 weeks at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank.

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method⁸.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	Rate A
Curve type	Linear
Wavelength (prim. / sec.)	340/405 (415) nm
Reading time	180–330 s after adding of R2
Reaction direction	Increase
Unit	U/L (µkat/L)
Reagent volumes	
R1	160 µL
R2	40 µL
Sample	4 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1 : R2 : sample does not change.

UNIT CONVERSION

U/L x 0.0167 = µkat/L

EXPECTED VALUES¹¹

At 37 °C

Females:	24–145 U/L	Children ¹² :	
Males:	46–171 U/L	Umbilical cord blood	175–402 U/L
		Newborns	468–1200 U/L
		≤5 days	195–700 U/L
		<6 months	41–330 U/L
		>6 months	24–229 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 5.12 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 2340 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	133.0	0.86	0.64
Sample 2	313.3	0.86	0.27

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	143.1	2.99	2.09
Sample 2	333.3	6.05	1.81

Accuracy

Two different validated control materials were used. Determined bias is 2.1 % at the target value 208.1 U/L and -2.5 % at the target value 513.2 U/L.

Comparison

A comparison between XL-640 automatic system CK (y) and a commercially available test (x) using 51 samples gave following results:

Linear regression:
 $y = 0.975x - 8.47$ U/L $r = 0.991$
 Passing-Bablok¹³:
 $y = 0.992x - 4.00$ U/L $r = 0.990$

Interferences

Criterion: Recovery within ±10 % of initial value of CK activity in the sample without interfering substance. Following substances do not interfere: haemoglobin up to 7 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels¹⁴.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of CK. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1
 UFI: JXTU-AWJV-3J5M-K095

R2
 UFI: 81UU-UW88-EJ53-7AV7

R1, R2

Danger

Contains: imidazole

Hazard statement:

H360D May damage the unborn child.

Precautionary statement:

P201 Obtain special instructions before use.

P280 Wear protective gloves / protective clothing / eye protection / face protection.

P308 + P313 IF exposed or concerned: Get medical advice / attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Oliver IT. A spectrophotometric method for the determination of creatine phosphokinase and myokinase. *Biochem J* 1955;61:116-122.
2. Rosalki SB. An improved procedure for serum creatine phosphokinase determination. *J Lab Clin Med* 1967;69:696-705.
3. Szasz G, Gruber W, Bernt E. Creatine kinase in serum: 1. Determination of optimum reaction conditions. *Clin Chem* 1976;22(5):650-656.
4. Standard method for the determination of creatine kinase activity. *J Clin Chem Clin Biochem* 1977;15:249-260.
5. Hørder M, Elser RC, Gerhardt M, et al. Approved Recommendation on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. Part 7. IFCC Method for Creatine Kinase. *Eur J Clin Chem Clin Biochem* 1991;29:435-456.
6. Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine Kinase. *Clin Chem Lab Med* 2002;40(6):635-642.
7. Klauke R, Schmidt E, Lorentz K. Recommendations for carrying out standard ECCLS procedures (1988) for the catalytic concentrations of creatine kinase, aspartate aminotransferase, alanine aminotransferase and γ -glutamyltransferase at 37 °C. *Eur J Clin Chem Clin Biochem* 1993;31:901-909.
8. Thomas L, ed. *Labor und Diagnose*, 8th ed. Bd 1:TH-Books Verlagsgesellschaft 2012.
9. Stein W. *Laboratory Diagnosis of Acute Myocardial Infarction*. Darmstadt: GIT Verlag 1988;34-37.
10. Guder WG, Narayanan S, Wissner H, et al. List of Analytes; Preanalytical Variables. Brochure in: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT-Verlag 1996.
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
12. Stein W. Creatine kinase (total activity), creatine kinase isoenzymes and variants. In: Thomas L, ed. *Clinical laboratory diagnostics*. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p.71-80
13. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
14. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385

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LITERATÚRA

1. Oliver IT. A spectrophotometric method for the determination of creatine phosphokinase and myokinase. *Biochem J* 1955;61:116-122.
2. Rosalki SB. An improved procedure for serum creatine phosphokinase determination. *J Lab Clin Med* 1967;69:696-705.
3. Szasz G, Gruber W, Bernt E. Creatine kinase in serum: 1. Determination of optimum reaction conditions. *Clin Chem* 1976;22(5):650-656.
4. Standard method for the determination of creatine kinase activity. *J Clin Chem Clin Biochem* 1977;15:249-260.
5. Hørder M, Elser RC, Gerhardt M, et al. Approved Recommendation on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. Part 7. IFCC Method for Creatine Kinase. *Eur J Clin Chem Clin Biochem* 1991;29:435-456.
6. Schumann G, Bonora R, Ceriotti F, et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C – Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine Kinase. *Clin Chem Lab Med* 2002;40(6):635-642.
7. Klauke R, Schmidt E, Lorentz K. Recommendations for carrying out standard ECCLS procedures (1988) for the catalytic concentrations of creatine kinase, aspartate aminotransferase, alanine aminotransferase and γ -glutamyltransferase at 37 °C. *Eur J Clin Chem Clin Biochem* 1993;31:901-909.
8. Thomas L, ed. *Labor und Diagnose*, 8th ed. Bd 1:TH-Books Verlagsgesellschaft 2012.
9. Stein W. *Laboratory Diagnosis of Acute Myocardial Infarction*. Darmstadt: GIT Verlag 1988;34-37.
10. Guder WG, Narayanan S, Wissler H, et al. List of Analytes; Preanalytical Variables. Brochure in: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT-Verlag 1996.
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
12. Stein W. Creatine kinase (total activity), creatine kinase isoenzymes and variants. In: Thomas L, ed. *Clinical laboratory diagnostics*. Frankfurt: TH-Books Verlagsgesellschaft;1998.p.71-80
13. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.
14. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. *Ann Clin Biochem* 2001;38:376-385

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