## PYRIDOXAL-5-PHOSPHATE 50

| Cat. No. | Pack Name | Packaging (Content)   |
|----------|-----------|---|
| 10003294 | PDP 50    | 1 bottle (contains 50 pcs of tablets),<br>instruction for use |

### INTENDED USE

The kit is intended as supplementary reagent for *in vitro* photometric quantitative determination of ALT/GPT and AST/GOT in human serum or plasma on various automatic systems. For professional use in clinical laboratory only.

### CLINICAL SIGNIFICANCE

In patients with vitamin B6 deficiency, serum aminotransferase activity may be decreased. The apparent reduction in aminotransferase activity may be related to decreased pyridoxal-5-phosphate (PDP), the prosthetic group for aminotransferases, resulting in an increase in the ratio of apoenzyme to holoenzyme. Addition of PDP, recommended by IFCC1, stabilizes the activity of ALT/GPT or AST/GOT and avoids falsely low values in samples containing insufficient endogenous PDP, e.g. form patients with myocardial infarction, liver disease and intensive care patients.

### PRINCIPLE

PDP functions as coenzyme in amino transfer reactions of alanine ALT/GPT and aspartate AST/GPT amino transferases. Enzymes with bound coenzyme as well as coenzyme deficient apoenzymes may be present in serum. Therefore, addition of PDP usually produces a marked increase in amino transferase activity. The increase ranges from zero to 3 or 4 fold and averages approximately 20 % for ALT/GPT and approximately 50 % for AST/GOT.

### REAGENT DESCRIPTION AND COMPOSITION

PDP 6 µmol/tablet

### REAGENT PREPARATION

With the following kits:

ALT/GPT 250 Cat. No. BLT00052  
ALT/GPT 330 Cat. No. XSYS0017  
AST/GOT 250 Cat. No. BLT00050  
AST/GOT 330 Cat. No. XSYS0016

Dissolve one tablet of PDP 50 in one bottle of the reagent R1.

With the following kits:

ALT/GPT 500 Cat. No. BLT00053  
AST/GOT 500 Cat. No. BLT00051

Dissolve two tablets of PDP 50 in one bottle of the reagent R1.

### STABILITY AND STORAGE

The unopened bottle with PDP tablets is stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. After opening PDP tablets are stable until expiry date at 2–8 °C if stored at original bottle, closed carefully, protected from light and without any contamination. After dissolving of PDP with reagent R1 (ALT/GPT or AST/GOT) stability 7 days at 2–8 °C.

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

The reagent is not classified as dangerous.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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### REFERENCES / LITERATURA / LITERATÚRA

- Schumann G, Bonora R, Cieriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 40: 725-33, 2002.
- Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 55-65, 1998.
- Klaue R, Schmidt L, Lorentz K. Recommendations for carrying out standard ECCLS procedures (1988) for the catalytic concentrations of creatine kinase, aspartate aminotransferase, alanine aminotransferase and gamma-glutamyltransferase at 37 °C. Standardization Committee of the German Society for Clinical Chemistry, Enzyme Working Group of the German Society for Clinical Chemistry. Eur J Clin Chem Clin Biochem 31: 901-9, 1993.

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Katalógové číslo

Číslo šarže

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Čítajte návod k použitiu

In vitro Diagnostikum

Výrobca

Teplota skladovania

Obsah

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

CC/IFU/062/22/A

Date of revision: 20. 4. 2022