



PHAN[®]

Test Strips for Urine Analysis



Making Automation Affordable
for Labs Everywhere

www.erba.com

PHAN[®] test strips for urine analysis

Quick Overview of Patient's Health

| Parameter | Abbreviation | Units | Evaluation time | Colour scale | Principle of the test |
|-------------------------|--------------|----------------------|-----------------|--------------|--|
| Haemoglobin | BLD | Ery/ μ l | ca 60 s | | oxidation of chromogene by organic hydroperoxide in the presence of the haemoglobine |
| Erythrocytes | | | | | |
| Ketones | KET | mmol/l mg/dl | ca 60 s | | sodium nitropruside in alkaline buffer (Legal's test) |
| Bilirubin | BIL | arb.u. | ca 60 s | | reaction of diazonium salt in acidic surroundings |
| Urobilinogen | UBG | μ mol/l mg/dl | ca 60 s | | reaction of diazonium salt in acidic surroundings |
| Glucose | GLU | mmol/l | ca 60 s | | enzymatic reaction – glucoseoxidase, peroxidase, chromogene |
| Protein | PRO | g/l mg/dl | ca 60 s | | protein error of pH indicator – mixed acido-basic indicator changes colour in the presence of proteins |
| pH | pH | | ca 60 s | | mixed acido-basic indicator |
| Nitrites | NIT | | ca 60 s | | modified Griess' reaction |
| Ascorbic Acid | AA | mmol/l mg/dl | ca 60 s | | reduction of molybdophosphoric acid into molybdenum blue |
| Specific Gravity | SG | | ca 60 s | | colour change of acidobasic indicator dependant on ion exchange |
| Leucocytes | LEU | Leu/ μ l | ca 120 s | | enzymatic reaction – esterase splits substrate into free indoxyl, which reacts with diazonium salt |
| Microalbumin | MA | g/l mg/l | ca 60 s | | acido-basic indicator changes colour in the presence of albumine |
| Creatinine | CRE | mmol/l g/l | ca 60 s | | Benedict-Behres' reaction |

| Sensitivity | | Specificity | Specificity Interference | |
|---------------------|-----------------|---|--|---|
| SI | Conv. | | Ascorbic Acid | Other |
| 5 Ery/ μ l | | specific for haemoglobin and myoglobin | All pads are protected against normal concentrations of ascorbic acid. | extremely high SG |
| 0.1–0.2 mmol/l | 1.0–2.0 mg/dl | high for acetoacetic acid, low for acetone, none for butyric acid | | drugs and diagnostics based on phenolphthalein or sulphophthalein |
| 4.3–5.2 μ mol/l | 0.25–0.30 mg/dl | specific for conjugated bilirubin | | high concentration of UBG and light |
| 6.0 μ mol/l | 0,35 mg/dl | urobilinogen and sterokobilinogen | | phenazopyridine, bilirubin and light |
| 0.9 mmol/l | 16 mg/dl | specific for D-glucose | | traces of detergents in the bases of peroxides and oxidizing agents |
| 0.15 g/l | 15 mg/dl | specific for albumin | | drugs based on quinine and quinoline, alkaline urine with pH > 8, traces of detergents and disinfectants based on quarternaryammonium salt and urine with high buffer capacity |
| | | | | foreign alkaline and/or acidic substances, old urine with pH about 9 |
| 11 mmol/l | 0.05 mg/dl | specific for nitrite (70 % of bacteriuria) | | diuresis and phenazopyridine |
| 0.2–0.3 mmol/l | 3.0–5.0 mg/dl | non specific oxidation – reduction reaction | | reducing agents present in the urine |
| | | | | pH > 6,5 |
| 10 Leu/ μ l | | granulocytes and histiocytes | | alkaline pH, higher SG and high concentration of bilirubin increase the intensity of colour reaction |
| 0.03 g/l | 30 mg/l | specific for albumine | | drugs based on quinine and quinoline, alkaline urine with pH > 8, traces of detergents and disinfectants based on quarternaryammonium salt and urine with a high buffer capacity, high concentration of creatinine (>26.5 mmol/l) |
| 0.4 mmol/l | 0.04 g/l | specific for creatinine | | urine with high buffer capacity decreases intensity of colour, high concentration of acetoacetic acid (>50 mmol/l) |

Visual strips PHAN®

| Product name | Cat. Nr. | Qty | Exp. | SG | NIT | pH | ASCO | PRO | GLU | KET | UBG | BIL | LEU | BLD | MA | CRE | CP |
|-----------------|----------|-----|------|----|-----|----|------|-----|-----|-----|-----|-----|-----|-----|----|-----|----|
| AlbuPHAN | URPH0001 | 50 | 24 | | | | | ● | | | | | | | | | |
| GlukoPHAN | URPH0002 | 50 | 27 | | | | | | ● | | | | | | | | |
| HemoPHAN | URPH0003 | 50 | 24 | | | | | | | | | | | ● | | | |
| KetoPHAN | URPH0004 | 50 | 30 | | | | | | | ● | | | | | | | |
| DiaPHAN | URPH0005 | 50 | 27 | | | | | | ● | ● | | | | | | | |
| IktoPHAN | URPH0006 | 50 | 24 | | | | | | | | ● | ● | | | | | |
| TriPHAN | URPH0007 | 50 | 27 | | | ● | | ● | ● | | | | | | | | |
| TriPHAN | URPH0008 | 100 | 27 | | | ● | | ● | ● | | | | | | | | |
| TetraPHAN dia | URPH0009 | 50 | 24 | | | ● | | ● | ● | ● | | | | | | | |
| PentaPHAN | URPH0010 | 50 | 24 | | | ● | | ● | ● | ● | | | | ● | | | |
| HexaPHAN | URPH0011 | 50 | 24 | | | ● | | ● | ● | ● | ● | | | ● | | | |
| HexaPHAN | URPH0012 | 100 | 24 | | | ● | | ● | ● | ● | ● | | | ● | | | |
| HeptaPHAN | URPH0013 | 50 | 24 | | | ● | | ● | ● | ● | ● | ● | | ● | | | |
| HeptaPHAN | URPH0014 | 100 | 24 | | | ● | | ● | ● | ● | ● | ● | | ● | | | |
| NonaPHAN SG | URPH0015 | 100 | 24 | ● | ● | ● | | ● | ● | ● | ● | ● | | ● | | | |
| NefroPHAN leuco | URPH0016 | 50 | 15 | | ● | ● | | ● | | | | | ● | ● | | | |
| DekaPHAN leuco | URPH0017 | 50 | 15 | ● | ● | ● | | ● | ● | ● | ● | ● | ● | ● | | | |
| DekaPHAN leuco | URPH0018 | 100 | 15 | ● | ● | ● | | ● | ● | ● | ● | ● | ● | ● | | | |
| UndekaPHAN | URPH0019 | 50 | 15 | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | | | |
| MicroalbuPHAN | URPH0020 | 50 | 21 | | | | | | | | | | | | ● | ● | |

Objective strips PHAN® LAURA

| | | | | | | | | | | | | | | | | | |
|---------------------|----------|-----|----|---|---|---|--|---|---|---|---|---|---|---|---|---|---|
| DiaPHAN LAURA | URPH0024 | 100 | 21 | | | | | | ● | ● | | | | | | | ● |
| TetraPHAN SG Laura | URPH0025 | 100 | 21 | ● | | ● | | ● | ● | | | | | | | | ● |
| PentaPHAN LAURA | URPH0026 | 100 | 21 | | | ● | | ● | ● | ● | | | | ● | | | ● |
| HeptaPHAN LAURA | URPH0027 | 100 | 21 | | | ● | | ● | ● | ● | ● | ● | | ● | | | ● |
| DekaPHAN LAURA | URPH0028 | 100 | 15 | ● | ● | ● | | ● | ● | ● | ● | ● | ● | ● | | | ● |
| MicroalbuPHAN LAURA | URPH0029 | 50 | 15 | | | | | | | | | | | | ● | ● | ● |

Specifications are subject to change without prior notice.
 Pictures used only for representation purposes - subject to change without prior notice.
 The devices are in compliance with the IVDR requirements for CE marking.

Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno, Czech Republic
 Phone: +420 517 077 111

✉ sales@erba.com

www.erba.com

