## EC cartridge XS plus iCa

Cat. No.	Pack Name	Number of Tests
REG00066	EC cartridge XS plus iCa	250

### INTENDED USE

IVD kit intended for calibration and quantitative determination of sodium (Na<sup>+</sup>), potassium (K<sup>+</sup>), ionized calcium (iCa<sup>2+</sup>) and chlorides (Cl<sup>-</sup>) in human serum, plasma, whole blood or diluted urine on the EC 90 analyser. Intended for screening, monitoring and diagnosis of diseases associated with ionic imbalance as kidney disease, disorders of acid-base balance, nutritional disorders and digestive tract function. For professional use in clinical laboratory only.

### PRINCIPLE

The electrolyte sensors are potentiometric sensors which measure the potential difference between an ion-selective electrode (ISE) and a reference electrode immersed in a sample liquid. The electrochemical potential between the sample side and the solid contact of the ion-selective membrane is proportional to the concentration of the corresponding ion in the sample. These potentials are recorded, electronically amplified and calculated using the Nernst Equation.<sup>1</sup>

### DESCRIPTION AND COMPOSITION

#### 1) EC cartridge

##### Calibration solution 1:

Na<sup>+</sup> 150 mmol/L, K<sup>+</sup> 5 mmol/L, iCa<sup>2+</sup> 1.25 mmol/L, Cl<sup>-</sup> 103 mmol/L, buffer, preservative, wetting agent

##### Calibration solution 2:

Na<sup>+</sup> 105 mmol/L, K<sup>+</sup> 2 mmol/L, iCa<sup>2+</sup> 2.34 mmol/L, Cl<sup>-</sup> 71.1 mmol/L, buffer, preservative, wetting agent

##### Reference solution:

K<sup>+</sup> 1.2 mol/L, buffer, preservative, wetting agent

2) EC Biosensor: Cat. No. REG00064 (for ordering extra pieces)

### PREPARATION

Insert EC Biosensor into the cartridge according to the EC 90 user manual.

### MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

EC 90 analyser, Cat. No. INS00079

EC urine diluent, Cat. No. REG00059

Control material: ERBA NORM 4 × 5, Cat. No. BLT00080

ERBA NORM 10 × 5, Cat. No. XSYS0123

ERBA PATH 4 × 5, Cat. No. BLT00081

ERBA PATH 10 × 5, Cat. No. XSYS0124

### STABILITY AND STORAGE

Store at 2–30 °C. Note the expiration date on the labels.

On board stability: 60 days after opening.

### SPECIMEN COLLECTION AND HANDLING

Whole blood samples may be directly measured or further processed to create plasma and serum samples. Take care to prevent haemolysis in blood samples during handling and processing of specimens. Urine must be diluted 1:3 (1 part urine and 2 parts diluent v/v) to permit accurate measurements. Use only the recommended EC urine diluent. Use only sample containers spray-coated with sodium heparin. Do not use EDTA or citrate as anticoagulants. This may result in incorrect measurement results.

### QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. It is recommended to perform QC measurements according to your local laboratory guidelines. More information about quality control can be found in the EC 90 User manual.

### METROLOGICAL TRACEABILITY

Traceability is ensured by testing laboratory standards according to standard EN ISO 17511:2021 against primary and secondary certified reference materials.

### ASSAY PROCEDURE AND CALCULATION

- In the analyser software, in the maintenance screen, click on CHANGE EC CARTRIDGE button. Follow instructions on the screen.
- Unpack new cartridge, remove the cap and seal cap, insert EC Biosensor and place it in the analyser.
- In the analyser software click the OK button. The software will detect new cartridge and perform priming automatically.
- Calibration will be performed automatically.

The EC 90 analyser performs the evaluation automatically according to Nernst equation<sup>1</sup>:

$$E_i = E_i^0 + \frac{R \cdot T}{z_i \cdot F} \ln a_i$$

#### Where:

- E<sub>i</sub> Ion activity dependent potential  
 E<sup>0</sup> Standard potential  
 R universal gas constant  
 T absolute Temperature  
 z<sub>i</sub> valency of ion i  
 F Faraday constant  
 a<sub>i</sub> activity of ion i

### EXPECTED VALUES

The products are intended for any general population without any limitation.<sup>2</sup>

Blood/Plasma/Serum [mmol/L], 24h	Parameter			
	Na <sup>+</sup>	K <sup>+</sup>	iCa <sup>2+</sup>	Cl <sup>-</sup>
Child	135–145	3.7–5.30	1.06–2.70	98–107
Adult	135–145	3.50–5.10	1.15–1.30	98–107

Urine [mmol/L], 24h	Parameter		
	Na <sup>+</sup>	K <sup>+</sup>	Cl <sup>-</sup>
Child	20–115	35–100	100–250
Adult	40–250	35–100	100–250

### ANALYTICAL PERFORMANCE

#### Precision:

Precision was determined using human samples and controls in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP05-A3 requirements with repeatability (n=20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days).

The following results were obtained: (S1 = sample 1, S2 = sample 2)

Repeatability	Plasma/Serum/Whole blood		Urine	
	CV (%) *			
	S1	S2	S1	S2
Na <sup>+</sup>	<1 %	<1 %	<5 %	<5 %
K <sup>+</sup>	<1 %	<1 %	<5 %	<5 %
Ca <sup>2+</sup>	<0.02 SD mmol/L	<0.02 SD mmol/L	-	-
Cl <sup>-</sup>	<1 %	<1 %	<5 %	<5 %

Intermediate precision	Plasma/Serum/Whole blood		Urine	
	CV (%) *			
	S1	S2	S1	S2
Na <sup>+</sup>	<2 %	<2 %	<5 %	<5 %
K <sup>+</sup>	<2.5 %	<2.5 %	<5 %	<5 %
Ca <sup>2+</sup>	<0.03 SD mmol/L	<0.03 SD mmol/L	-	-
Cl <sup>-</sup>	<2.5 %	<2.5 %	<5 %	<5 %

\* acceptance criteria for iCa were based on SD

### COMPARISON

A comparison between EC 90 and competitor's system with 1300 patient samples in clinical setting.

#### Analysed EC 90

##### Pearson correlation

Parameter	Na <sup>+</sup> (mmol/L)	K <sup>+</sup> (mmol/L)	iCa <sup>2+</sup> (mmol/L)	Cl <sup>-</sup> (mmol/L)
R	0.968	0.994	0.954	0.920

#### Competitor's analyser

##### Pearson correlation

Parameter	Na <sup>+</sup> (mmol/L)	K <sup>+</sup> (mmol/L)	iCa <sup>2+</sup> (mmol/L)	Cl <sup>-</sup> (mmol/L)
R	0.963	0.994	0.971	0.974

### Interferences

The following substances were found to significantly interfere with respective measurement parameters of the instrument.

Substance	Test concentration	Significant interference with				Test matrix
		c Na <sup>+</sup>	c K <sup>+</sup>	c iCa <sup>2+</sup>	c Cl <sup>-</sup>	
Magnesium (Mg <sup>2+</sup> )	>2.2 mmol/L	–	–	✓	–	Plasma
Bromide (Br <sup>-</sup> )	>5.3 mmol/L	–	–	–	✓	Plasma
	>12.0 mmol/L	–	–	–	✓	Urine
Iodide (I <sup>-</sup> )	>3.0 mmol/L	–	–	–	✓	Plasma
Thiocyanate (SCN <sup>-</sup> )	>3.9 mmol/L	–	–	–	✓	Plasma
	>3.6 mmol/L	–	–	–	✓	Urine
Hydrogen-phosphate (HPO <sub>4</sub> <sup>2-</sup> )	>34.0 mmol/L	✓	–	–	–	Urine
	–	–	–	✓	–	Plasma
pH	6.9–7.9	–	–	✓	–	Whole Blood
		–	–	✓	–	Plasma
Citrate	>0.4 mmol/L	–	–	✓	–	Plasma
EDTA	>0.5 mmol/L	–	–	✓	–	Plasma
Ibuprofen	>5.7 mmol/L	–	✓	–	–	Urine

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

#### Additional information:

##### Calibration solution 1 and 2:

EUH 208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

Reference solution: The reagent is not classified as dangerous.

### WASTE MANAGEMENT

EC cartridge needs to be considered as a potentially infectious material. Please refer to local legal requirements.

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