C-REACTIVE PROTEIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0047	CRP	R1: 2 x 40 mL, R2: 2 x 10 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* quantitative determination of C-reactive protein in human serum on automatic systems ERBA XL by turbidimetric immunoassay. Intended for monitoring of cardiovascular disease and screening and monitoring of occurrence, progression, and thrombo-occlusive complications of atherosclerosis of all population. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

C-Reactive Protein (CRP) is a non-specific acute phase-reactive protein which appears in the blood during an inflammatory process. CRP is the most sensitive of the acute phase reactants and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the classical complement pathway. The CRP response frequently precedes clinical symptoms, including fever. In normal healthy individuals CRP is a trace protein with a range up to 5 mg/L. After onset of an acute phase response the serum CRP concentration rises rapidly and extensively. The increase begins within 6 to 12 hours and the peak value is reached within 24 to 48 hours. Levels above 100 mg/L are associated with severe stimuli such as major trauma and severe infection (sepsis). CRP response may be less pronounced in patients suffering from liver disease. CRP assays are used to detect systemic inflammatory processes; to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amnio rhexis; to differentiate between active and inactive forms of disease with concurrent infection; to therapeutically monitor rheumatic disease and assess anti-inflammatory therapy; to determine the presence of post-operative complications at an early stage, such as infected wounds, thrombosis and pneumonia, and to distinguish between infection and bone marrow rejection. Postoperative monitoring of CRP levels of patients can aid in the recognition of unexpected complications (persisting high or increasing levels). Persistence of a high serum CRP concentration is usually a grave prognostic sign which generally indicates the presence of an uncontrolled infection^{1,2,3,4,5,6,7,8}.

PRINCIPLE

The human serum sample reacts upon a specific antibody for human CRP and the turbidity induced by the formation of immune complexes is recorded at 340 nm. The turbidity measured is directly proportional to the CRP sample concentration⁴.

REAGENT DESCRIPTION AND COMPOSITION

R1	
Phosphate buffered saline (pH 7.43)	
Polyethylene glycol	40 g/L
Sodium azide	0.09 %

R2	
Phosphate buffered saline (pH 7.43)	
Polyclonal goat anti-human CRP	variable
Sodium azide	0.09 %

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

CRP CAL SH, Cat. No. XSYS0053
CRP CON L, Cat. No. BLT20014
CRP CON H, Cat. No. BLT20013
MULTICON L 1, Cat. No. BLT 20034
MULTICON L 2, Cat. No. BLT20035
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 60 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum⁹:	11 days at	20–25 °C
	2 months	4–8 °C
	3 years	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with CRP CAL SH is recommended.

6-point calibration: generate a 6-point calibration curve by diluting the calibrator 1:32, 1:16, 1:8, 1:4, 1:2 and undiluted in saline; saline is used as blank

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control CRP CON L, CRP CON H, MULTICONTROL LEVEL 1 and MULTICON LEVEL 2 are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the reference material BCR470/CRM470⁹.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com

ASSAY PARAMETERS FOR ERBA XL AUTOMATIC SYSTEMS

Assay type	2-Point
Curve type	Cubic Spline
Wavelength (prim. / sec.)	340/700 nm
Reading time 1	just before adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (mg/L)

Reagent volumes	
R1	180 µL
R2	45 µL
Sample volumes	12 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION

mg/dL × 10 = mg/L

EXPECTED VALUES^{6,7}

In Serum:	
Newborn up to 3 weeks	<0.41 mg/dL
Infants and children	<0.28 mg/gL
Adults	<0.5 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com

Limit of quantification: 0.11 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30)

Linearity: approx. 22 mg/dL

Linearity is given by concentration of the calibrator used for generating 6-point calibration curve.

Hook Effect: safety zone >110 mg/dL

Precision

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	2.05	0.032	1.59
Sample 2	4.21	0.015	0.37

Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.95	0.090	4.62
Sample 2	4.23	0.109	2.57

Accuracy

Two different validated control materials were used. Determined bias is and -14.3 % at the target value 1.84 mg/dL and -6.3 % at the target value 3.55 mg/dL.

Comparison

A comparison between XL-640 automatic system CRP (y) and a commercially available test (x) using 146 samples gave following results:

Linear regression:
y = 0.921x + 0.474 mg/dL r = 0.984

Passing-Bablok⁸:
y = 0.992x + 0.174 mg/dL r = 0.985

Interferences

Criterion: Recovery within ±10 % of initial value of albumin concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 9.3 g/L, bilirubin up to 30 mg/dL, triglycerides up to 200 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁹.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of CRP. See paragraph Interferences.

Contraindications

Genetic polymorphisms within the CRP gene, particularly single-nucleotide polymorphisms (SNPs) in the promoter region, can influence baseline CRP levels. These genetic variations may account for individual differences in CRP concentrations and should be considered when interpreting results¹⁰.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. The Summary of Safety and Performance is available at Eudamed see <https://ec.europa.eu/tools/eudamed>.

Health protection

Used animal-source materials are accompanied by an official state veterinarian health certificate. The animal components were collected from healthy animals under the care of a registered veterinarian and were apparently free from infectious and contagious diseases and injurious parasites. However, this material should be handled as thought capable of transmitting infectious diseases.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Johnson AM, Rohlfis EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER. editors. Tietz textbook of clinical chemistry. 3rd ed. Philadelphia: W. B. Saunders Company; 477-540, 1999.
2. Young B, Gleeson M, Cripps AW. C-reactive protein: A critical review. Pathology 23: 118-124, 1991.
3. Wasunna A, Whitelaw A, Gallimore R, et al. C-reactive protein and bacterial infection in preterm infants. Eur J Pediatr, Mar;149(6): 424-427, 1990.
4. Hafner G, Larkin T, Peetz D, Erbes H, Eyres H, Wrynn K, Lackner KJ. An improved assay for the determination of C-reactive protein. Clin Lab. 48(7-8): 369-376, 2002.
5. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
6. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 34: 517-520, 1996.
7. Schlebusch H, Liappis N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood. Poster presented at AACC/CSCC; July/August 2001, Chicago, Illinois.
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.
10. Lopez-Roblero A, Serrano-Guzmán E, Guerrero-Báez RS, Delgado-Enciso I, Gómez-Manzo S, Aguilar-Fuentes J, Ovando-Garay V, Hernández-Ochoa B, Quezada-Cruz IC, Lopez-Lopez N, Canseco-Ávila LM, Single-nucleotide polymorphisms in the promoter of the gene encoding for C-reactive protein associated with acute coronary syndrome. Biomed Rep. 21(5):150, 2024, doi: 10.3892/br.2024.1838.

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