

ADVISORY NOTICE

34_2026_Change in the Instructions of use - ANAEROTest 23

Dear customers,

Before using the suspension medium, carefully read the current Instructions for Use for **ANAEROTest 23**. The section "**Inoculum Preparation**" has been revised.

If you are currently using ANAEROTest 23 with a previous version of the Instructions for Use, please follow the updated inoculum preparation procedure described below.

Updated Inoculum Preparation:

Before preparing the bacterial suspension, **measure the optical density of the suspension medium alone (initial OD)**.

The turbidity of the prepared bacterial suspension must correspond to **McFarland 3**. Therefore, when measuring the bacterial suspension, **add the initial OD of the suspension medium to the OD value corresponding to McFarland 3** and prepare the bacterial suspension to the resulting final OD.

A bacterial suspension that is too dilute or too concentrated may lead to false reactions.

Please apply this updated procedure when preparing the bacterial suspension for use with **ANAEROTest 23**.

The revised Instructions for Use should be used for all future testing. If you are using a previous version of the Instructions for Use, please ensure that this Advisory Notice is kept together with the product documentation and is available to all users of the ANAEROTest 23 kit.

In case you have any additional questions, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o



Brno, 1st October 2026

Ref. No.AN_34_2026