MULTICONTROL LEVEL 2

Cat. No.	Pack Name	Packaging (Content)
BLT20035	MULTICON L2	1 × 1 mL, instruction for use

EN

CE 2797 IVD

INTENDED USE
MULTICONTROL LEVEL 2 is control material for *in vitro* quantitative determination of proteins in human serum on automatic systems ERBA XL by turbidimetric immunoassay. Intended for monitoring accuracy and precision for the quantitative methods as specified in the value sheet. For professional use in clinical laboratory only.

PRINCIPLE
This human plasma control is to be used in the same way as an unknown specimen that would be used in accordance with the direction of the assay procedure. The results obtained for the control are to be compared with the assigned values given in the "ASSAYED VALUES" section of the insert and an evaluation made by standard statistical techniques to determine if the procedure is within the control limits.

DESCRIPTION AND COMPOSITION
The control is diluted human plasma containing high levels of parameters with phosphate buffered saline, stabilized with 0.095 % sodium azide.

ASSAYED VALUES

Control Lot	Expiry	Parameter	Value	Interval	1 SD	Unit
xxxxxxxxx	yyyy-mm	ASO	00.0	00.0 – 00.0	0.00	IU/mL (kU/L)
		CRP	0.00	0.00 – 0.00	0.00	mg/dL
			00.0	00.0 – 00.0	0.00	mg/L
		RF	00.0	00.0 – 00.0	0.00	IU/mL (kU/L)

UNIT CONVERSION
IU/mL = kU/L
mg/dL × 10 = mg/L

CONTROL PREPARATION
Control is liquid, ready to use.

STORAGE AND HANDLING
The unopened MULTICON L2 stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
After first opening the vial, control can be used for 6 weeks if stored tightly closed at 2–8 °C after use.
Do not freeze.

MATERIALS REQUIRED, BUT NOT PROVIDED
ASO, Cat. No. XSYS0046
CRP, Cat. No. XSYS0047
CRP-HS, Cat. No. XSYS0084
RF, Cat. No. XSYS0048
ASO CAL SH, Cat. No. XSYS0051
CRP CAL SH, Cat. No. XSYS0053
CRP-HS CAL, Cat. No. XSYS0088
RF CAL SH, Cat. No. XSYS0052
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

ASSAY PROCEDURE
Pour the required volume into a sample cup and analyze in the same way as patient samples.
The controls should be run daily in parallel with patient samples and after every calibration. Control intervals must be adapted to individual laboratory's requirements.
Follow the applicable government regulations and local guidelines for quality control.

EVALUATION
The XL-analysers perform the evaluation automatically.

EXPECTED RESULTS
The target value stated in the "ASSAYED VALUES" section. Determination was performed under strictly standardized conditions on Erba analyzers using Erba system reagents and appropriate calibrators.
The target value specified is the mean of all values obtained. The corresponding control range is calculated as the target value ±3 standard deviations. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material. Results should be within the defined ranges.
For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range.

TRACEABILITY
Traceability information of individual analytes is given in the respective instructions for use of the appropriate kits.

PERFORMANCE CHARACTERISTICS
Performance characteristics of individual analytes on the Erba XL automatic systems are available on websites www.erba.com.

LIMITATIONS
Erroneous results can occur from technical errors associated with the assay procedure.
Improper storage or handling of the control can also affect the results. Do not use the control if there is visible evidence of microbial growth in the vial.

WARNING AND PRECAUTIONS
For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Health protection
Human material used in the manufacture of this product has been tested and found non-reactive for hepatitis B surface antigen (HbsAg), antibody to Hepatitis C (HCV), HIV-1 p24 antigen, and antibody to HIV-1/2. The test procedures do not guarantee that all infectious agents will be detected. Because no test method can offer complete assurance that Hepatitis B virus Hepatitis C virus and HIV-1/2 or other infectious agents are absent, the material should be handled as potentially infectious. During the work all principles of general safety as well as the principles of work with biological materials are to be maintained. Above all observe the principles of personal hygiene, use personal protection aids (clothes, shoes, etc.), do not eat, drink and smoke. After finishing the work and before a meal wash your hands with soap and warm water and apply a suitable regenerating cream. When the bottle is broken, the contaminated surface must be disinfected with Persteril or Chloramine.

Hazards identification in accordance with Regulation (EC) No 1272/2008
The reagent is not classified as dangerous.

WASTE MANAGEMENT
Please refer to local legal requirements.

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