C-REACTIVE PROTEIN HIGH SENSITIVE

Cat. No.	Pack Name	Packaging (Content)
XSYS0084	CRP-HS	R1: 2 × 40 mL, R2: 2 × 11 mL, RFID tag, instruction for use

EN

CE 2797 IVD

INTENDED USE
The kit is intended for *in vitro* quantitative determination of low concentration of C-reactive protein in human serum on automatic systems ERBA XL by turbidimetric immunoassay. Intended for monitoring of cardiovascular disease and screening and monitoring of occurrence, progression, and thrombo-occlusive complications of atherosclerosis. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE
C-Reactive Protein (CRP) is an acute marker of inflammatory processes. In case of an acute inflammation the concentration of CRP increases and decreases more quickly than the red cell sedimentation rate. The increase of CRP occurs in a non-specific way in different kinds of tissular aggression, as for example in infectious states, rheumatoid arthritis, myocardial infarct, malignant tumour, etc.
Although not diagnostic it is very useful for following-up and monitoring such illnesses, as well as for differential diagnosis in certain cases.
Routinely available immunochemical assay methods for CRP have limited sensitivity, and until recently, CRP concentrations below 1.0 mg/dL could not be measured precisely, leading to the widespread adoption of this value as the upper limit of the health-associated reference range. This is satisfactory for most purposes in general medicine. However, in neonatal paediatric practice, a highly sensitive CRP immunoassay shows that health-associated reference values are below 0.1–0.2 mg/dL and that any rise above such values is associated with serious disease, usually bacterial infection.
More recently, application of sensitive CRP assays to studies of adult cardiovascular disease has revealed important prognostic relationships between modest increase of CRP and the occurrence, progression, and thrombo-occlusive complications of atherosclerosis[1,2,3,4,5,6,7,8].

PRINCIPLE
Determination of the concentration of CRP by photometric measurement of antigen-antibody reaction of antibodies to human CRP bound to latex particles with CRP present in the sample. The turbidity of immune complex, measured at 600 nm, is directly proportional to the CRP sample concentration.

REAGENT DESCRIPTION AND COMPOSITION
R1

Sodium chloride	9 g/L
Detergent	0.1%
Sodium azide	0.09%

R2

Glycine buffer, pH 8.42	
Rabbit anti-human CRP sensitized latex	0.20%
Sodium azide	0.09%

REAGENT PREPARATION
Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE
CRP-HS CAL, Cat. No. XSYS0088
CRP CON L, Cat. No. BLT20014,
CRP CON H, Cat. No. BLT20013,
MULTICON L 1, Cat. No. BLT 20034
MULTICON L 2, Cat. No. BLT20035
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

STABILITY AND STORAGE
The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
On board stability: min. 60 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING
It is recommended to follow ISO 15189 and laboratory instruction.
For specimen collection and preparation only use suitable tubes or collection containers.
Only the specimens listed below were tested and found acceptable.
Serum.
The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.
Centrifuge samples containing precipitates before performing the assay.
See the limitations and interferences section for details about possible sample interferences.

Stability in serum[9]:

11 days at	20–25 °C
2 months at	4–8 °C
3 years at	-20 °C

Discard contaminated specimens.

CALIBRATION
Calibration with CRP-HS CAL is recommended.
6-point calibration: generate a 6-point calibration curve by diluting the calibrator 1:32, 1:16, 1:8, 1:4, 1:2 and undiluted in saline; saline is used as blank
Calibration frequency: 30 days
Calibration is needed:
• after reagent lot change
• as required by internal quality control procedures
• calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL
For quality control CRP CON L, CRP CON H, MULTICON L1 and MULTICON L2 are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY
This method, calibrator CRP-HS CAL and controls CRP CON L, CRP CON H, MULTICON L1 and MULTICON L2 have been standardized against the reference material ERM-DA474/IFCC[10].

ASSAY PROCEDURE AND CALCULATION
ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

ASSAY PARAMETERS FOR ERBA XL AUTOMATIC SYSTEMS

Assay type	2-Point
Curve type	5P Calibration Logit-Log
Wavelength (prim. / sec.)	600 / - nm
Reading time 1	just before adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (mg/L)
Reagent volumes	
R1	250 µL
R2	60 µL
Sample volumes	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION
mg/dL × 10 = mg/L

EXPECTED VALUES[11,12]
Cardiovascular Risk Classification:

Low risk	<0.10 mg/dL
Average risk	0.10–0.30 mg/dL
High risk	>0.30 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE
Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 0.058 mg/dL
Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20% (n = 30)

Linearity: approx. 15 mg/dL
Linearity is given by concentration of the calibrator used for generating 6-point calibration curve.

Hook Effect: safety zone >225 mg/dL

Precision
Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	2.13	0.022	1.02
Sample 2	3.90	0.029	0.75

Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	2.14	0.056	2.64
Sample 2	3.90	0.083	2.12

Accuracy
Two different validated control materials were used. Determined bias is and 11.5% at the target value 1.84 mg/dL and 7.2% at the target value 3.55 mg/dL.

Comparison
A comparison between XL-640 automatic system CRP-HS (y) and a commercially available test (x) using 150 samples gave following results:
Linear regression:
y = 1.013x - 0.075 mg/dL r = 0.996
Passing-Bablok[13]:
y = 1.009x - 0.050 mg/dL r = 0.994

Interferences
Criterion: Recovery within ±10% of initial value of CRP concentration in the sample (serum) without interfering substance.
Following substances do not interfere: haemoglobin up to 13.3 g/L, bilirubin up to 5 mg/dL, triglycerides up to 1910 mg/dL.
Drugs: No interference was found at therapeutic concentrations using common drug panels[14].

Limitations:
- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of CRP. See paragraph Interferences.

Contraindications
Genetic polymorphisms within the CRP gene, particularly single-nucleotide polymorphisms (SNPs) in the promoter region, can influence baseline CRP levels. These genetic variations may account for individual differences in CRP concentrations and should be considered when interpreting results[15].

WARNING AND PRECAUTIONS
For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Health protection
Used animal-source materials are accompanied by an official state veterinarian health certificate. The animal components were collected from healthy animals under the care of a registered veterinarian and were apparently free from infectious and contagious diseases and injurious parasites. However, this material should be handled as thought capable of transmitting infectious diseases.

Hazards identification in accordance with Regulation (EC) No 1272/2008
R1, R2
The reagents are not classified as hazardous.

WASTE MANAGEMENT
Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com
IC/IFU/004/26/A *Datum revize: 8. 10. 2026*

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1. Claus DR, Osmand AP, Gewurz H, Radioimmunoassay of human C-reactive protein and levels in normal sera, J Lab Clin Med 87, 120-128, 1976.
2. Wasunna A, Whitelaw A, Gallimore R, Hawkins PN, Pepys MB, Creactive protein and bacterial infection in preterm infants, Eur J Pediatr 149, 424-427, 1990
3. Heinrich J, Schulte H, Schönfeld R, Köhler E, Assmann G, Association of variables of coagulation, fibrinolysis and acute-phase with atherosclerosis in coronary and peripheral arteries and those arteries supplying the brain, Thromb Haemostas 73, 374-379, 1995.
4. Grau AJ, Buggle F, Becher H, Werle E, Hacke W, The association of leukocyte count, fibrinogen and C-reactive protein with vascular risk factors and ischemic vascular diseases, Thromb Res 82, 245-255, 1996.
5. Dixon JS, Bird HA, Sitton NG et al, C-Reactive protein in the serial assessment of disease activity in rheumatoid arthritis, Scand J Rheum 13, 39-44, 1984.
6. Hind CRH, and Pepys MB, The role of serum C-Reactive protein (CRP) measurement in clinical practice, Int Med 5, 112-151, 1984.
7. Hanson LÅ, Wadsworth Ch, Das C-reacktive Protein und sein diagnosticher Wert, insbesondere bei infektionen. Laboratoriumblätter 29, 58-68, 1979.
8. Macy EM, Hayes TE, Tracy RP, Variability in the measurement of C-reactive protein in healthy subjects: implications reference intervalsand epidemiological applications. Clin Chem 43, 52-58, 1997.
9. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
10. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 34: 517-520, 1996.
11. Ridker PM, et al, C-reactive protein adds to the predictive value of total and HDL cholesterol in determining risk of first myocardial infarction, Circulation 97, 2007-2011, 1998.
12. Pearson T A, Mensah GA,.Alexander RW, Anderson JL, Canon RO 3rd, Criqui M, Fadl YY, Fortmann SP, Hong Y, Myers GL, Rifai N, Smith SC Jr, Taubert K, Tracy RP, Vinicor F, Markers of Inflammation and Cardiovascular Disease: Application to Clinical and Public Health Practice: A Statement for Healthcare Professionals From the Centers for Disease Control and Prevention and the American Heart Association, Circulation 107, 499-511, 2003.
13. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
14. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.
15. Lopez-Roblero A, Serrano-Guzmán E, Guerrero-Báez RS, Delgado-Enciso I, Gómez-Manzo S, Aguilar-Fuentes J, Ovando-Garay V, Hernández-Ochoa B, Quezada-Cruz IC, Lopez-Lopez N, Canseco-Ávila LM, Single-nucleotide polymorphisms in the promoter of the gene encoding for C-reactive protein associated with acute coronary syndrome. Biomed Rep. 21(5):150, 2024, doi: 10.3892/br.2024.1838.

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