BILIRUBIN TOTAL DCA

Cat. No.	Pack Name	Packaging (Content)
XSYS0087	BIL T DCA 330	R1: 6 × 44 mL, R2: 6 × 11 mL, RFID tag, instruction for use

EN

CE 2797 IVD

INTENDED USE
The kit is intended for *in vitro* photometric quantitative determination of bilirubin total in human serum and plasma on automatic systems ERBA XL. In combination with other clinical and laboratory findings, the assay is intended as an aid in the diagnosis and monitoring of liver and biliary disorders. The assay is also intended for monitoring bilirubin levels in neonates as an aid in the clinical management of neonatal hyperbilirubinemia. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE
Red blood cells at the end of their circulating lives are broken down in the reticuloendothelial system, mainly the spleen. The resulting heme is converted to bilirubin upon removal of iron. This process accounts for about 80% of the 500 µmol (292 mg) of bilirubin formed daily. Other sources of bilirubin include the breakdown of myoglobin and cytochromes and the catabolism of immature red blood cells in the bone marrow.
Once formed, bilirubin is transported to the liver bound to albumin. This fraction of bilirubin is referred to as unconjugated (indirect) bilirubin. In the liver, bilirubin is conjugated to glucuronic acid (mono- and diglucuronides) by the enzyme uridyl diphosphate glucuronyl transferase to form conjugated bilirubin. Conjugated (direct) bilirubin is excreted via the biliary system into the intestine, where it is metabolized by bacteria to a group of products known collectively as stercobilinogen. Elimination is almost complete and serum levels are normally negligible.
Total bilirubin is the sum of the unconjugated and conjugated fractions. Total bilirubin is elevated in hepatitis, cirrhosis, hemolytic disorders, several inherited enzyme deficiencies, and conditions causing hepatic obstruction.

TEST PRINCIPLE
Total bilirubin is coupled with diazotized 2,4-dichloroaniline in strongly acidic medium to produce an intensely coloured diazo compound, Azobilirubin. Free bilirubin is not soluble in aqueous media, but this reagent contains an association of surfactant and accelerators in order to provide an accurate measurement of total bilirubin[1,2,3,4].

Bilirubin + Diazotized 2,4-dichloroaniline —(Acidic medium)→ Azobilirubin

Absorbance of the Azobilirubin measured at 546 nm is proportional to the concentration of the total bilirubin in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
HCl	55 mmol/L
EDTA	0.1 mmol/L
Surfactant	18 mmol/L
R2	
HCl	176 mmol/L
2,4-dichloroaniline	1.0 mmol/L
EDTA	0.1 mmol/L
Sodium nitrite	0.8 mmol/L
COMPOSITION OF REACTION MIXTURE	
HCl	75.4 mmol/L
EDTA	0.095 mmol/L
Surfactant	13.7 mmol/L
2,4-dichloroaniline	0.19 mmol/L
Sodium nitrite	0.15 mmol/L

REAGENT PREPARATION
Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE
XL MULTICAL 4×3, Cat. No. XSYS0034
XL MULTICAL 10×3, Cat. No. XSYS0122
ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

STABILITY AND STORAGE
The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING
It is recommended to follow ISO 15189 and laboratory instruction.
For specimen collection and preparation only use suitable tubes or collection containers.
Only the specimens listed below were tested and found acceptable.
Serum.
Plasma: K_2-EDTA.
The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.
Centrifuge samples containing precipitates before performing the assay.
See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma[5]: 1 day at 15–25 °C
7 days at 4–8 °C
6 months at -20 °C
Protect sample from light.
Discard contaminated specimens.

CALIBRATION
Calibration with calibrator XL MULTICAL is recommended.
2 point calibration (blank and calibrator); distilled water is recommended as blank
Calibration frequency: 30 days
Calibration is needed:
• after reagent lot change
• as required by internal quality control procedures
• calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL
For quality control ERBA NORM and ERBA PATH are recommended.
The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY
This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the reference material NIST SRM 916.

ASSAY PROCEDURE AND CALCULATION
ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type		2-Point
Curve type		Linear
Wavelength (prim./sec.)		546/660 nm
Reading time 1		just before adding of R2
Reading time 2		10 min after adding of R1
Reaction direction		Increase
Unit		mg/dL (µmol/L)
Reagent volumes		
R1		200 µL
R2		50 µL
Sample volumes:	adults	12.5 µL
	neonates	12.5 µL diluted 5 times

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION
mg/dL × 17.1 = µmol/L

EXPECTED VALUES[9]

In serum:	Premature	Full term	
Cord	<2.0	<2.0	mg/dL
0–1 d	1.0–8.0	2.0–6.0	mg/dL
1–2 d	6.0–12.0	6.0–10.0	mg/dL
3–5 d	10.0–14.0	4.0–8.0	mg/dL
Adult		0–2.0	mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE
Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com

Limit of quantification: Adults: 0.01 mg/dL
Neonates: 0.10 mg/dL
Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: Adults, neonates: 30 mg/dL
Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:
Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (Adults)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.47	0.013	0.87
Sample 2	5.21	0.023	0.45

Intermediate precision (Adults)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.47	0.041	2.79
Sample 2	5.23	0.072	1.38

Repeatability (Neonates)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.41	0.038	2.68
Sample 2	4.27	0.119	2.79

Intermediate precision (Neonates)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.51	0.037	2.45
Sample 2	4.39	0.128	2.91

Accuracy
Two different validated control materials were used. Adults: determined bias is 1.3 % at the target value 0.65 mg/dL and 3.1 % at the target value 1.91 mg/dL. Neonates: determined bias is 5.2 % at the target value 4.76 mg/dL and 7.0 % at the target value 7.31 mg/dL.

Comparison
Adults: A comparison between XL-640 automatic system BILIRUBIN TOTAL DCA (y) and a commercially available test (x) using 97 samples gave following results:
Linear regression:
y = 0.934x + 0.0656 mg/dL r = 0.999
Passing-Bablok[10]:
y = 0.953x + 0.0234 mg/dL r = 0.992

Neonates: A comparison between XL-1000 automatic system BILIRUBIN TOTAL DCA (y) and a commercially available test (x) using 100 samples gave following results:
Linear regression:
y = 0.954x + 0.930 mg/dL r = 0.988
Passing-Bablok[10]:
y = 0.982x + 0.472 mg/dL r = 0.978

Interferences
Criterion: Recovery within ±10 % of initial value of total bilirubin concentration in the sample without interfering substance.
Following substances do not interfere: Adults: haemoglobin up to 2 g/L, triglycerides up to min. 1700 mg/dL. Neonates: haemoglobin up to 1 g/L, triglycerides up to min. 1700 mg/dL.
Paraproteins can affect assay results[11,12,13].

Limitations:
- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin and triglycerides in sample can interfere with determination of total bilirubin. See paragraph Interferences.

WARNING AND PRECAUTIONS
For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety and Performance is available at Eudamed see:
https://ec.europa.eu/tools/eudamed

Hazards identification in accordance with Regulation (EC) No 1272/2008
R1, R2
Reagents are not classified as dangerous.

WASTE MANAGEMENT
Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com
CC/IFU/056/26/B *Date of revision: 11. 11. 2026*

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1. Wahlefeld AW, Herz G, Bernt E. Modification of the Malloy-Evelyn method for a simple, reliable determination of total bilirubin in serum. Scand J Clin Lab Invest 29: Supplement 126: Abstract 11.12, 1972.
2. Rand RN, di Pasqua A, A new diazo method for the determination of bilirubin. Clin Chem 6: 570-578, 1962.
3. Mc Phaual L, Kershaw M, Tilque D, J.H. Eckfeldt JH, A 2,4-dichlorophenyl diazonium-based method for total bilirubin without interference from indican in uremic sera, Clin. Chem. 31: 1229-1231, 1985.
4. Dapper GM, Sundberg MW, Wu TW, Babb BE, J.R. Schaeffer JR, A diazo-based dry film for determination of total bilirubin in serum, Clin. Chem. 29: 37-41, 1983.
5. Quality of Diagnostic Samples, Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory medicine, 3rd completely revised ed. 2010.
6. Balistreri WF, Shaw LM. Liver function. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders: 729-761, 1987.
7. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 192-202, 1998.
8. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1125-77, 1999.
9. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
10. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988
11. Sarkar R, Paraprotein Interferences: Insights from a Short Study Involving Multiple Platforms and Multiple Measurands, EJIFCC 35 (2), 120-131, 2024.
12. Yilmaz NS, Sen B, Gulbahar O, Contribution of the laboratory to a diagnosis process by sequential reflective testing: Paraprotein interference on a direct bilirubin assay, Biochem Med (Zagreb 31:020801, 2021, doi: 10.11613/BM.2021.020801.
13. Chauhan G, Dutta B, Chandra A, Paraproteins: A Potential Source of Fallacy in Bilirubin Assay, J Med Sci 11(1 -4):00385, 2025.

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