APOLIPOPROTEIN A1

Cat. No.	Pack Name	Packaging (Content)
XSYS0102	APO A1	R1: 2 × 27.5 mL, R2: 2 × 6.7 mL, RFID tag, instruction for use

EN

CE 2797 IVD

INTENDED USE
The kit is intended for *in vitro* quantitative determination of apolipoprotein A1 in human serum on automatic systems ERBA XL by turbidimetric immunoassay. Intended for screening, monitoring and diagnosis of atherosclerotic vascular changes, dyslipoproteinemias of all population. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE
Apolipoprotein A1 (Apo A1) is the main protein component of HDL (High Density Lipoprotein). Apo A1 activates lecithin cholesterol acyltransferase which catalyses the esterification of cholesterol. The resulting esterified cholesterol can then be transported to the liver, metabolised, and excreted.
Apo A1 levels increase in liver disease, pregnancy and as a result of estrogen administration (e.g. oral contraceptives). Apo A1 levels decrease in inherited hypo-α-lipoproteinemia (e.g. Tangier disease), cholestasis, sepsis and atherosclerosis. The liver also synthesizes very low-density lipoproteins (VLDL) which mainly contain triglycerides and cholesterol. In the presence of lipoprotein lipase, the triglycerides are hydrolyzed and LDL-particles with a high proportion of cholesterol are formed. Apolipoprotein B (Apo B) is the main constituent of LDL.
The combined determination of Apo A1/Apo B and the calculation of the Apo B: Apo A1 ratio can reflect a lipid metabolism disorder and the risk of developing atherosclerosis or coronary heart disease particularly well, thus providing an excellent addition to the classical HDL/LDL-cholesterol determination. A high level of Apo A1 (HDL) and a low level of Apo B (LDL) correlate best with a low risk for these diseases[1,2].

PRINCIPLE
The Apo A1 assay is an immunoturbidimetric procedures that measure increasing sample turbidity caused by the formation of insoluble immune complexes. The human serum (plasma) sample reacts upon a specific antibody for human Apo A1 and the turbidity induced by the formation of immune complexes is recorded at 600 nm. The turbidity measured is directly proportional to the Apo A1 sample concentration.

REAGENT DESCRIPTION AND COMPOSITION
R1

TRIS Buffer (pH 8)	200 mmol/L
PEG	56 g/L
Sodium chloride	50 mmol/L
Sodium azide	0.9 g/L

R2
Sterile delipidated goat serum against human Apo A1

HEPES buffer (pH 7.4)	50 mmol/L
EDTA	9 mmol/L
Sodium azide	0.9 g/L

REAGENT PREPARATION
Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE
APO CAL SET, Cat. No. XSYS0105
ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

STABILITY AND STORAGE
The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING
It is recommended to follow ISO 15189 and laboratory instruction.
For specimen collection and preparation only use suitable tubes or collection containers.
Only the specimens listed below were tested and found acceptable[3].
Serum, heparinised, EDTA or citrated plasma.
The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.
Centrifuge samples containing precipitates before performing the assay.
See the limitations and interferences section for details about possible sample interferences.

Stability in serum[3]: 1 day at 20–25 °C
8 days at 4–8 °C
3 months at -20 °C

Discard contaminated specimens.

CALIBRATION
Calibration with APO CAL SET is recommended.
5-point calibration: generate a 5-point calibration curve by using of set of calibrators APO CAL SET; saline is used as blank
Calibration frequency: 30 days
Calibration is needed:
- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL
For quality control ERBA NORM and EDRBA PATH are recommended.
The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY
This method, calibrator APO CAL SET and controls ERBA NORM and ERBA PATH have been standardized against the reference material IFCC SP1-01 (WHO-IRP)[4,5,6,7].

ASSAY PROCEDURE AND CALCULATION
ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	2-Point
Curve type	Cubic Spline
Wavelength (prim. / sec.)	600 / - nm
Reading time 1	just before adding of R2
Reading time 2	10 min after adding of R1
Reaction direction	Increase
Unit	mg/dL (g/L)
Reagent volumes	
R1	250 µL
R2	50 µL
Sample volumes	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION
mg/dL × 0.01 = g/L

EXPECTED VALUES[8]
Serum:

4–5 y	Male	109–172 mg/dL	40–49 y	Male	103–178 mg/dL
	Female	104–163 mg/dL		Female	115–195 mg/dL
6–11 y	Male	111–177 mg/dL	50–59 y	Male	107–173 mg/dL
	Female	110–166 mg/dL		Female	117–211 mg/dL
12–19 y	Male	99–165 mg/dL	60–69 y	Male	111–184 mg/dL
	Female	105–180 mg/dL		Female	120–205 mg/dL
20–29 y	Male	105–173 mg/dL	>69 y	Male	109–180 mg/dL
	Female	111–209 mg/dL		Female	118–199 mg/dL
30–39 y	Male	105–173 mg/dL			
	Female	110–189 mg/dL			

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE
Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com

Limit of quantification: 1.70 mg/dL
Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: approx. 210 mg/dL
Linearity is given by concentration of the calibrator used for generating 5-point calibration curve.

Hook Effect: safety zone >2300 mg/dL

Precision:
Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	173.5	3.13	1.80
Sample 2	219.7	2.81	1.28

Intermediate precision	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	176.2	6.67	3.78
Sample 2	222.7	4.09	1.84

Accuracy
Two different validated control materials were used. Determined bias is and 8.3 % at the target value 104 mg/dL and 3.1 % at the target value 214 mg/dL.

Comparison
A comparison between XL-640 automatic system APO A1 (y) and a commercially available test (x) using 130 samples gave following results:
Linear regression:
y = 0.917x + 7.076 mg/dL r = 0.981
Passing-Bablok[9]:
y = 0.937x + 5.062 mg/dL r = 0.979

Interferences
Criterion: Recovery within ±10 % of initial value of APO A1 concentration in the sample (serum) without interfering substance.
Following substances do not interfere: haemoglobin up to 4.0 g/L, bilirubin up to 50 mg/dL, triglycerides up to 2600 mg/dL.
Drugs: No interference was found at therapeutic concentrations using common drug panels[10].

Limitations:
- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of APO A1. See paragraph Interferences.

WARNING AND PRECAUTIONS
For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Health protection
Used animal-source materials are accompanied by an official state veterinarian health certificate. The animal components were collected from healthy animals under the care of a registered veterinarian and were apparently free from infectious and contagious diseases and injurious parasites. However, this material should be handled as thought capable of transmitting infectious diseases.

Hazards identification in accordance with Regulation (EC) No 1272/2008
R1, R2
Reagents are not classified as dangerous.

WASTE MANAGEMENT
Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com
IC/IFU/025/26/A *Date of revision: 16. 10. 2026*

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1. Riesen WF, Apolipoproteine. In: Thomas L, ed. Labor und Diagnose, 5th ed. Frankfurt 171-190, 1998.
2. Brewer HB Jr, Gregg RE, Hoeg JM, et al, Apolipoproteins and lipoproteins in human plasma: an overview. Clin Chem. 34, B4-B8, 1988.
3. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
4. Marcovina SM, Albers JJ, Dati F, et al, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B, Clin Chem 37, 1676-1682, 1991.
5. Albers JJ, Marcovina SM, Kennedy H, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. II. Evaluation and Selection of Candidate Reference Materials, Clin Chem 38, 658-662, 1992.
6. Marcovina SM, Albers JJ, Henderson LO, et al, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. III. Comparability of Apolipoprotein A-I Values by Use of International Reference Material, Clin Chem 39, 773-781, 1993.
7. Marcovina SM, Albers JJ, Kennedy H, et al, International Federation of Clinical Chemistry Standardization Project for Measurements of Apolipoproteins A-I and B. IV. Comparability of Apolipoprotein B Values by Use of International Reference Material, Clin Chem 40, 586-592, 1994.
8. Burtis CA, Ashwood ER, Bruns DE, Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th edition: W.B. Saunders Company Ltd., 2012.
9. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
10. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38: 376-385, 2001.

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