

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

SECTION 1: Identification of the substance/mixture and of the company/undertaking

- 1.1. Product identifier** URIC ACID (SINGLE REAGENT)_R1
Substance / mixture mixture
Number BLT00062; XSYS0042
Other mixture names
UA SINGLE 200; UA 440
- 1.2. Relevant identified uses of the substance or mixture and uses advised against**
Mixture's intended use
Diagnostic reagent for quantitative in vitro determination of Uric Acid in human serum, plasma and urine.
Main intended use
PC-MED-OTH Other medical devices
Secondary uses
PC-TEC-19 Reagents and laboratory chemicals
The use descriptors
PC 21 Laboratory chemicals
Mixture uses advised against
not available
- 1.3. Details of the supplier of the safety data sheet**
Manufacturer
Name or trade name Erba Lachema s.r.o.
Address Karásek 2219/1d , Brno, 62100
Czech Republic
Identification number (CRN) 26918846
VAT number CZ26918846
Phone +420 517 077 111
Email msds@erba.com
Web address www.erba.com
- Competent person responsible for the safety data sheet**
Name Erba Lachema s.r.o.
Email msds@erba.com
- 1.4. Emergency telephone number**
European emergency number: 112 112

SECTION 2: Hazards identification

- 2.1. Classification of the substance or mixture**
Classification of the mixture in accordance with Regulation (EC) No 1272/2008
The mixture is not classified as dangerous according to Regulation (EC) No 1272/2008.
- 2.2. Label elements**
Signal word
none
- 2.3. Other hazards**
The mixture does not contain substances with endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605. Mixture does not contain any substance meet the criteria for PBT or vPvB in accordance with Annex XIII of Regulation (EC) No. 1907/2006 (REACH) as amended. Does not contain any PMT or vPvM components.

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

SECTION 3: Composition/information on ingredients**3.2. Mixtures****Chemical characterization**

Mixture of substances and additives specified below.

Mixture contains these hazardous substances and substances with the highest permissible concentration in the working environment

Identification numbers	Substance name	Content in % weight	Classification according to Regulation (EC) No 1272/2008	Note
Index: 011-002-00-6 CAS: 1310-73-2 EC: 215-185-5 Registration number: 01-2119457892-27-0000	sodium hydroxide	<0.4	Met. Corr. 1, H290 Skin Corr. 1A, H314 Specific concentration limit: Skin Corr. 1B, H314: 2 % ≤ C < 5 % Skin Corr. 1A, H314: C ≥ 5 % Eye Irrit. 2, H319: 0.5 % ≤ C < 2 % Skin Irrit. 2, H315: 0.5 % ≤ C < 2 %	1
Index: 011-004-00-7 CAS: 26628-22-8 EC: 247-852-1	sodium azide	<0.1	Acute Tox. 2, H300 Aquatic Acute 1, H400 Aquatic Chronic 1, H410 EUH032	1
CAS: 58-15-1 EC: 200-365-8	4-Dimethylaminoantipyrine	<0.1	Acute Tox. 4, H302+H312 Skin Irrit. 2, H315 Eye Irrit. 2, H319 STOT SE 3, H335	
CAS: 6381-92-6 EC: 205-358-3	Disodium EDTA	<0.04	Acute Tox. 4, H332 STOT RE 2, H373	
CAS: 14348-40-4	2,4,6-Tribromo-3-hydroxybenzoic Acid	<0.02	Acute Tox. 4, H302+H312+H332 Skin Irrit. 2, H315 Eye Irrit. 2, H319 STOT SE 3, H335	
CAS: 83-07-8 EC: 201-452-3	4-Aminoantipyrine	<0.01	Acute Tox. 4, H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319 STOT SE 3, H335	
Index: 604-014-00-3 CAS: 59-50-7 EC: 200-431-6	chlorocresol	<0.01	Acute Tox. 4, H302 Skin Corr. 1C, H314 Skin Sens. 1B, H317 Eye Dam. 1, H318 STOT SE 3, H335 Aquatic Acute 1, H400 (M=1) Aquatic Chronic 3, H412	

Notes

1 A substance for which exposure limits are set.

Full text of all classifications and hazard statements is given in the section 16.

SECTION 4: First aid measures**4.1. Description of first aid measures**

Take care of your own safety. If any health problems are manifested or if in doubt, inform a doctor and show him information from this safety data sheet.

If inhaled

Terminate the exposure immediately; move the affected person to fresh air.

If on skin

Remove contaminated clothes. After contact with skin, wash with soap and water.

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

If in eyes

Rinse eyes immediately with a flow of running water, open the eyelids (also using force if needed); remove contact lenses immediately if worn by the affected person.

If swallowed

Rinse out the mouth with clean water. In the event of issues, find medical help.

4.2. Most important symptoms and effects, both acute and delayed**If inhaled**

Not expected.

If on skin

Not expected.

If in eyes

Not expected.

If swallowed

Irritation, nausea.

4.3. Indication of any immediate medical attention and special treatment needed

Symptomatic treatment.

SECTION 5: Firefighting measures**5.1. Extinguishing media****Suitable extinguishing media**

Alcohol-resistant foam, carbon dioxide, powder, water spray jet, water mist.

Unsuitable extinguishing media

Water - full jet.

5.2. Special hazards arising from the substance or mixture

Hazardous fumes may be released during combustion. Inhalation of hazardous degradation (pyrolysis) products may cause serious health damage.

5.3. Advice for firefighters

Self-Contained Breathing Apparatus (SCBA) with a chemical protection suit only where personal (close) contact is likely. Use a self-contained breathing apparatus and full-body protective clothing. Do not allow run-off of contaminated fire extinguishing material to enter drains or surface and ground water.

SECTION 6: Accidental release measures**6.1. Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment for work. Follow the instructions in the Sections 7 and 8. Prevent contact with skin and eyes. Observe the principles of work safety in chemical laboratories.

6.2. Environmental precautions

Prevent contamination of the soil and entering surface or ground water.

6.3. Methods and material for containment and cleaning up

After removal of the product, wash the contaminated site with plenty of water.

6.4. Reference to other sections

See the Section 7, 8 and 13.

SECTION 7: Handling and storage**7.1. Precautions for safe handling**

Prevent formation of gases and vapours in concentrations exceeding the occupational exposure limits. Prevent contact with skin and eyes. Use personal protective equipment as per Section 8. Observe valid legal regulations on safety and health protection.

7.2. Conditions for safe storage, including any incompatibilities

Store in tightly closed containers in cold, dry and well ventilated areas designated for this purpose.

Storage temperature

min 2 °C, max 8 °C

7.3. Specific end use(s)

For in vitro diagnostic devices.

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

SECTION 8: Exposure controls/personal protection**8.1. Control parameters**

The mixture contains substances for which occupational exposure limits are set.

Czech Republic**Government Regulation 20/2025 Coll.**

Substance name (component)	Type	Value
sodium hydroxide (CAS: 1310-73-2)	PEL	1 mg/m ³
	NPK-P	2 mg/m ³

Notes

Irritating to mucous membranes (eyes, respiratory system) and skin.

Czech Republic**Government Regulation 20/2025 Coll.**

Substance name (component)	Type	Value
sodium azide (CAS: 26628-22-8)	PEL	0,1 mg/m ³
	NPK-P	0,3 mg/m ³

Notes

Skin penetration is significantly involved during exposure.

Irritating to mucous membranes (eyes, respiratory system) and skin.

European Union**Commission Directive 2000/39/EC**

Substance name (component)	Type	Value
sodium azide (CAS: 26628-22-8)	OEL 8 hours	0,1 mg/m ³
	OEL 15 minutes	0,3 mg/m ³

Notes

Skin.

8.2. Exposure controls

Do not eat, drink and smoke during work. Wash your hands thoroughly with water and soap after work and before breaks for a meal and rest.

Eye/face protection

It is not needed.

Skin protection

When handling in long-term or repeatedly, use protective gloves.

Respiratory protection

Not required with adequate ventilation.

Thermal hazard

Not available.

Environmental exposure controls

Observe usual measures for protection of the environment, see Section 6.2.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties**

Physical state	liquid
Colour	colourless
Odour	without fragrance
Melting point/freezing point	data not available
Boiling point or initial boiling point and boiling range	data not available
Flammability	The product is non-flammable.
Lower and upper explosion limit	data not available
Flash point	data not available
Auto-ignition temperature	data not available
Decomposition temperature	data not available
pH	7.1-7.5 (undiluted at 20 °C)
Kinematic viscosity	data not available

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

Solubility in water	data not available
Solubility in fats	data not available
Partition coefficient n-octanol/water (log value)	data not available
Vapour pressure	data not available
Density and/or relative density	data not available
Relative vapour density	data not available
Particle characteristics	data not available
Form	Clear colourless liquid

9.2. Other information

Evaporation rate	data not available
Oxidising properties	The product has no oxidizing properties.
Explosive properties	The product does not have explosive properties.

SECTION 10: Stability and reactivity**10.1. Reactivity**

not available

10.2. Chemical stability

The product is stable under normal conditions.

10.3. Possibility of hazardous reactions

Unknown.

10.4. Conditions to avoid

The product is stable and no degradation occurs under normal use. Protect against flames, sparks, overheating and against frost.

10.5. Incompatible materials

Protect against strong acids, bases and oxidizing agents.

10.6. Hazardous decomposition products

Not developed under normal uses.

SECTION 11: Toxicological information**11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008**

No toxicological data is available for the mixture.

Acute toxicity

Based on available data the classification criteria are not met.

URIC ACID (SINGLE REAGENT)_R1

Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	ATE	51870 mg/kg				Calculation of value
Dermal	ATE	3143000 mg/kg				Calculation of value
Inhalation (vapor)	ATE	22270 mg/l				Calculation of value

4-Aminoantipyrine

Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	1700 mg/kg		Rat		

4-Dimethylaminoantipyrine

Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	285 mg/kg		Rat		

URIC ACID (SINGLE REAGENT)_R1

Creation date 06th October 2015 Version 5.0
 Revision date 08th January 2026

Disodium EDTA						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	2800 mg/kg		Rat		
Inhalation	LC ₅₀	1000-5000 mg/m ³	4 hours	Rat		

chlorocresol						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	1830 mg/kg		Rat		
Dermal	LD ₅₀	>2000 mg/kg		Rat		
Inhalation	LC ₅₀	2871 mg/l	48 hours	Rat		

sodium azide						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	27 mg/kg bw				
Inhalation	LC ₅₀	54 mg/m ³	4 hours	Rat		

sodium hydroxide						
Route of exposure	Parameter	Value	Exposure time	Species	Sex	Value determination
Oral	LD ₅₀	325 mg/kg		Rat		
Dermal	LD ₅₀	1350 mg/kg		Rabbit		
Oral	LD ₅₀	500 mg/kg		Rabbit		
Dermal	LD ₅₀	40 mg/kg		Mouse		

Skin corrosion/irritation

Based on available data the classification criteria are not met.

Serious eye damage/irritation

Based on available data the classification criteria are not met.

Respiratory or skin sensitisation

Based on available data the classification criteria are not met.

Germ cell mutagenicity

Based on available data the classification criteria are not met.

Carcinogenicity

Based on available data the classification criteria are not met.

Reproductive toxicity

Based on available data the classification criteria are not met.

Toxicity for specific target organ - single exposure

Based on available data the classification criteria are not met.

Toxicity for specific target organ - repeated exposure

Based on available data the classification criteria are not met.

Aspiration hazard

Based on available data the classification criteria are not met.

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

11.2. Information on other hazards**Endocrine disrupting properties**

Based on the available data, the criteria for classification of the mixture are not met. Does not contain any components that may cause endocrine disruption for humans.

Other information

not available

SECTION 12: Ecological information**12.1. Toxicity**

not available

Acute toxicity

Disodium EDTA				
Parameter	Value	Exposure time	Species	Environment
EC ₅₀	140 mg/l	48 hours	Daphnia (Daphnia magna)	
LC ₅₀	320 mg/l	96 hours	Fish (Poecilia reticulata)	
EC ₅₀	56 mg/l	8 hours	Bacteria (Pseudomonas putida)	
NOEC	25.7 mg/l	35 days	Fish (Branchydanio rerio)	
NOEC	25 mg/l	21 days	Daphnia (Daphnia magna)	

chlorocresol				
Parameter	Value	Exposure time	Species	Environment
LC ₅₀	3.11-5.27 mg/l	96 hours	Fish (Oncorhynchus mykiss)	Fresh water
EC ₅₀	2 mg/l	48 hours	Crustaceans (Buchanka)	
EC ₅₀	4.2 mg/l	72 hours	Algae	Fresh water
EC ₅₀	10 mg/l	96 hours	Algae	Fresh water

sodium azide				
Parameter	Value	Exposure time	Species	Environment
LC ₅₀	680 µg/l		Fish	Fresh water
EC ₅₀ /LC ₅₀	400 µg/l		Invertebrates	Fresh water
EC ₅₀ /LC ₅₀	150 µg/l		Invertebrates	Salt water
EC ₅₀ /LC ₅₀	348 µg/l		Algae	Fresh water
EC ₅₀ /LC ₅₀	5.6 mg/l		Microorganisms	
NOEC	30 µg/l		Microorganisms	

sodium hydroxide				
Parameter	Value	Exposure time	Species	Environment
EC ₅₀	76 mg/l	24 hours	Daphnia (Daphnia magna)	
EC ₅₀	145 mg/l	24 hours	Fish (Poecilia reticulata)	
EC ₅₀	40.4 mg/l	48 hours	Daphnia (Daphnia magna)	
LC ₅₀	160 mg/l	24 hours	Fish (Carassius auratus)	

12.2. Persistence and degradability

Data not available.

12.3. Bioaccumulative potential

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

Not available.

12.4. Mobility in soil

Not available.

12.5. Results of PBT and vPvB assessment

Based on the available data, the criteria for classification of the mixture are not met. Does not contain any PBT or vPvB components. Product does not contain any substance meeting the criteria for PBT or vPvB in accordance with the Annex XIII of Regulation (EC) No 1907/2006 (REACH) as amended.

12.6. Endocrine disrupting properties

Based on the available data, the criteria for classification of the mixture are not met. Does not contain any components that may cause endocrine disruption in the environment.

12.7. Other adverse effects

Not available.

SECTION 13: Disposal considerations**13.1. Waste treatment methods**

Hazard of environmental contamination; dispose of the waste in accordance with the local and/or national regulations. Any unused product and contaminated packaging should be put in labelled containers for waste collection and submitted for disposal to a person authorised for waste removal (a specialized company) that is entitled for such activity. Do not empty unused product in drainage systems. The product must not be disposed of with municipal waste. Empty containers may be used at waste incinerators to produce energy or deposited in a dump with appropriate classification. Perfectly cleaned containers can be submitted for recycling.

Waste management legislation

Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste, as amended. Decision 2000/532/EC establishing a list of wastes, as amended.

Waste type code

16 05 06* laboratory chemicals, consisting of or containing hazardous substances, including mixtures of laboratory chemicals

Packaging waste type code

15 01 10* packaging containing residues of or contaminated by hazardous substances

(*) - Hazardous waste according to Directive 2008/98/EC on hazardous waste

SECTION 14: Transport information**14.1. UN number or ID number**

not subject to transport regulations

14.2. UN proper shipping name

not relevant

14.3. Transport hazard class(es)

not relevant

14.4. Packing group

not relevant

14.5. Environmental hazards

not relevant

14.6. Special precautions for user

Reference in the Sections 4 to 8.

14.7. Maritime transport in bulk according to IMO instruments

not relevant

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

Regulation (EC) No. 1907/2006 of the European Parliament and of the Council of 18th December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing the European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No. 793/93 and Commission Regulation (EC) No. 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC, as amended. REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16th December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No. 1907/2006, as amended. Commission Regulation (EU) 2020/878 of 18 June 2020 amending Annex II to Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).

15.2. Chemical safety assessment

not available

SECTION 16: Other information**A list of standard risk phrases used in the safety data sheet**

EUH032	Contact with acids liberates very toxic gas.
H290	May be corrosive to metals.
H300	Fatal if swallowed.
H302	Harmful if swallowed.
H302+H312	Harmful if swallowed or in contact with skin.
H302+H312+H332	Harmful if swallowed, in contact with skin or if inhaled.
H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.
H335	May cause respiratory irritation.
H373	May cause damage to organs through prolonged or repeated exposure.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.

Other important information about human health protection

The product must not be - unless specifically approved by the manufacturer/importer - used for purposes other than as per the Section 1. The user is responsible for adherence to all related health protection regulations.

Key to abbreviations and acronyms used in the safety data sheet

Acute Tox.	Acute toxicity
ADR	Agreement concerning the international carriage of dangerous goods by road
Aquatic Acute	Hazardous to the aquatic environment
Aquatic Chronic	Hazardous to the aquatic environment (chronic)
ATE	Acute toxicity estimate
BCF	Bioconcentration Factor
CAS	Chemical Abstracts Service
CLP	Regulation (EC) No 1272/2008 on classification, labelling and packaging of substance and mixtures
EC	Identification code for each substance listed in EINECS
EC ₅₀	Concentration of a substance when it is affected 50 % of the population
EINECS	European Inventory of Existing Commercial Chemical Substances
EmS	Emergency Response Procedures for Ships Carrying Dangerous Goods
EU	European Union
EuPCS	European Product Categorisation System
Eye Dam.	Serious eye damage
Eye Irrit.	Eye irritation

URIC ACID (SINGLE REAGENT)_R1

Creation date	06th October 2015	Version	5.0
Revision date	08th January 2026		

IATA	International Air Transport Association
IBC	International Code For The Construction And Equipment of Ships Carrying Dangerous Chemicals
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
INCI	International Nomenclature of Cosmetic Ingredients
ISO	International Organization for Standardization
IUPAC	International Union of Pure and Applied Chemistry
LC ₅₀	Lethal concentration of a substance in which it can be expected death of 50% of the population
LD ₅₀	Lethal dose of a substance in which it can be expected death of 50% of the population
log K _{ow}	Octanol-water partition coefficient
Met. Corr.	Corrosive to metals
NOEC	No observed effect concentration
NPK	Maximum admissible concentration
OEL	Occupational Exposure Limits
PBT	Persistent, bioaccumulative and toxic
PEL	Permissible Exposure Limit
PMT	Persistent, mobile and toxic
ppm	Parts per million
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
RID	Regulation concerning the International Carriage of Dangerous Goods by Rail
Skin Corr.	Skin corrosion
Skin Irrit.	Skin irritation
Skin Sens.	Skin sensitization
STOT RE	Specific target organ toxicity - repeated exposure
STOT SE	Specific target organ toxicity - single exposure
UN number	Four-figure identification number of the substance or article taken from the UN Model Regulations
UVCB	Substances of unknown or variable composition, complex reaction products or biological materials
VOC	Volatile organic compounds
vPvB	Very persistent and very bioaccumulative
vPvM	Very persistent and very mobile

Training guidelines

Inform the personnel about the recommended ways of use, mandatory protective equipment, first aid and prohibited ways of handling the product.

Recommended restrictions of use

not available

Information about data sources used to compile the Safety Data Sheet

REGULATION (EC) No. 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL (REACH) as amended.
REGULATION (EC) No. 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL as amended. Data from the manufacturer of the substance / mixture, if available - information from registration dossiers.

The changes (which information has been added, deleted or modified)

The version 5.0 replaces the SDS version from Wednesday, 16 August 2023. Changes were made in sections 2, 11, 12, 13 and 16.

Statement

The safety data sheet provides information aimed at ensuring safety and health protection at work and environmental protection. The provided information corresponds to the current status of knowledge and experience and complies with valid legal regulations. The information should not be understood as guaranteeing the suitability and usability of the product for a particular application.