

ADVISORY NOTICE

21_2025_Maintenance kit for EC 90 change

Dear customer,

We hereby inform you that the current **EC 90 Maintenance Kit** (Cat. No. 10020754) needed for the EC 90 Periodic maintenances will be from **7th May 2025** replaced by **2 new separate kits**:

Cat. No.10020959 – EC 90 REPAIR KIT / 3 YEARS MAINTENANCE KIT

Content:

- Pump Tube
- Tube Set
- Barb to Female adapter
- Four-port Cross packed
- Sample Needle packed
- Needle block packed
- Needle guide packed
- Set Fluidic connector
- Set Reagent fluidic connector
- Rubber fan packed
- Other (part list, tube lengths incl. fluidic scheme, etc.)

Cat. No.10020960 – EC 90 MAINTENANCE KIT / 1 YEAR









Content:

- Pump Tube
- Soft EC-90 Tube set and Needle to Valve 2 FEP
- Set Reagent fluidic Connector
- Set Fluidic Connector
- Needle block packed
- Other (part list, tube lengths incl. fluidic scheme, etc.)







Note: The detailed contents of both kits are listed in the Annex below.

The User Manual will be updated in accordance with the new kit details, however there is no change in the periodic maintenance procedure itself.

Annex I - EC 90 REPAIR KIT / 3 YEARS MAINTENANCE KIT

No.	Part Name	Pcs./Set	Picture
1	Pump tube	1	
2	Tube set	1	
3	Barb to Female adapter	2	
4	Four-port Cross packed	1	
5	Sample Needle packed	1	
6	Needle block packed	1	
7	Needle guide packed	1	
8	Set Fluidic connector	1	
9	Set Reagent fluidic connector	1	
10	Rubber fan packed	5	
11	Part list + tube lengths incl. fluidic scheme	1	

Annex II – EC 90 MAINTENANCE KIT / 1 YEAR

No.	Part Name	Pcs./Set	Picture
1	Pump tube	1	
2	Soft EC 90 Tube set and Needle to Valve 2 FEP	1	
3	Set Reagent fluidic connector	1	
4	Set Fluidic connector	1	
5	Needle block packed	1	
6	Part list + tube lengths incl. fluidic scheme	1	

In case you have any additional question, please do not hesitate to contact us.

Best regards,

Erba Lachema s.r.o.

Brno, 23rd May 2025

Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno

DIČ: CZ26918846 -7-



Ref. No.AN_21_2025